


ELMED RIGID & ROTATABLE 1-PIECE LANDRENEAU STYLE THORACIC INSTRUMENTS



MAINTENANCE AND CARE

CLEANING AND CARE: EVERY INSTRUMENT MUST BE CLEANED AND STERILIZED BEFORE BEING USED FOR THE FIRST TIME AND AFTER EVERY SUBSEQUENT USE. APPROPRIATE CLEANING, INSPECTION AND MAINTENANCE; HELP TO ENSURE THE SERVICEABILITY OF SURGICAL INSTRUMENTS. CLEAN, INSPECT AND TEST ALL INSTRUMENTS THOROUGHLY, AND STERILIZE BEFORE USE. EFFECTIVE CLEANING AND MAINTENANCE WILL PROLONG SERVICE LIFE.

CLEANING AND RINSING SHOULD BE DONE PROMPTLY AFTER EVERY USE. OTHERWISE TISSUE PARTICLES OR DRIED SECRETIONS MAY ADHERE TO IT, WHICH MAY MAKE SUBSEQUENT CLEANING AND STERILIZATION DIFFICULT, IF NOT IMPOSSIBLE. INSTRUMENTS MUST BE ENTIRELY FREE OF ANY FOREIGN BODIES.

HANDLING, STORAGE AND STERILIZATION: EXTREME CARE MUST BE EXERCISED WHILE USING, TRANSPORTING, CLEANING, STERILIZING, AND STORING. DAMAGE TO THE FUNCTION AND SAFETY OF THE INSTRUMENT MAY OCCUR IF THE INSTRUMENT IS HANDLED ROUGHLY OR IMPROPERLY, OR USED FOR SOMETHING OTHER THAN ITS INTENDED PURPOSE.

CLEANING

MANUAL CLEANING: THE INSTRUMENT IS CLEANED THROUGH THE FLUSHPORT. USE THE APPROPRIATE CONNECTIONS; WE SUGGEST THE ELMED CLEANING GUN. PREPARE A CLEANING BATH ACCORDING TO THE MANUFACTURER'S INSTRUCTIONS.

1. RINSE ALL PARTS WITH DEIONIZED COLD TAP WATER TO REMOVE AND DEBRIS.
2. HAND WASH ALL PARTS USING NEUTRAL OR ENZYMATIC CLEANER*¹ AND A SOFT BRUSH FOR CLEANING EXTERNAL PARTS AND ALL LUMENS AND CHANNELS. THOROUGHLY CLEAN THE INSIDE OF THE INSULATED OR INSULATED STAINLESS STEEL SHAFT USING A SOFT ROUND BRUSH*² AND FLUSHING SOLUTION THROUGH IT. DO NOT USE STEEL WOOL, WIRE BRUSHES AND/OR ABRASIVE DETERGENTS. CLEAN*³ ENTIRE SURFACE OF INSERT WITH SOFT BRUSH PAYING SPECIAL ATTENTION TO THE TIP AREA TO BE SURE ALL BLOOD AND DEBRIS IS REMOVED FROM THE JAW LINKAGE, THE INSIDE AND OUTSIDE OF THE JAW HOUSING, THE ENTIRE SHAFT AND RINSE WITH FLUSHING SOLUTION.
3. THOROUGHLY RINSE INSTRUMENT WITH DEIONIZED WATER AFTER CLEANING TO REMOVE ANY RESIDUAL DEBRIS OR CLEANING SOLUTION.

*¹ RECOMMENDED MERITZ PLUS DISINFECTANT/DECONTAMINANT: ITEM NO. **52003-90**.

*² RECOMMENDED ELMED CLEANING BRUSHES, ITEM NO. **52003-31** & **52003-42**.

*³ RECOMMENDED ELMED WATER JET CLEANING KIT, ITEM NO. **52003-80** & **52003-81**.

MACHINE REPROCESSING:

PLACE THE INSTRUMENTS IN A BASKET ON THE INSERT MODULE OR ON THE INSERTS OF THE MIS MODULE AND START THE CLEANING PROCESS.

1. PRE-RINSE FOR 1 MINUTE WITH COLD WATER
2. DISCHARGING
3. PRE-RINSE FOR 3 MINUTES WITH COLD WATER
4. DISCHARGING
5. WASH FOR 5 MINUTES AT 55°C WITH A 0.5% ALKALINE OR AT 45°C WITH AN ENZYMATIC CLEANING AGENT
6. DISCHARGING
7. NEUTRALIZE FOR 2 MINUTES WITH WARM TAP WATER (>40°C) AND A NEUTRALIZING AGENT
8. DISCHARGING
9. RINSE FOR 2 MINUTES WITH WARM TAP WATER (>40°C)
10. DISCHARGE

DISINFECTION: MACHINE OPERATED THERMAL DISINFECTION HAS TO BE CARRIED OUT IN CONSIDERATION OF THE NATIONAL REQUIREMENTS WITH REGARD TO THE A0 VALUE (SEE ISO15883).

DRYING: DRY THE OUTSIDE OF THE INSTRUMENTS BY CARRYING OUT A DRYING CYCLE OF THE CLEANING / DISINFECTION MACHINE. IF NECESSARY, MANUAL DRYING MAY ADDITIONALLY BE CARRIED OUT USING A LINT FREE CLOTH. DRY CAVITIES BY BLOWING WITH STERILE COMPRESSED AIR.

STERILIZATION: THE INSTRUMENT MUST BE STERILIZED FULLY ASSEMBLED AFTER IT HAS BEEN CLEANED.

IMPORTANT: THE FLUSHPORT OR RINSING CHANNEL CAP MUST BE REMOVED WHEN THE INSTRUMENT IS BEING STERILIZED TO KEEP THE INNER INSTRUMENT CAVITY OPEN AT BOTH ENDS.

1. **STEAM AUTOCLAVING:** WHEN USING THE WRAP METHOD, MAKE CERTAIN THAT ALL COMPLETE INSTRUMENTS AND ALL PARTS ARE INDIVIDUALLY WRAPPED OR SEALED IN A STERILE PACK. METAL OBJECTS SHOULD NEVER COME IN CONTACT WITH THE INSULATING MATERIAL OF FORCEPS AND HANDLES, OR WITH RF-CONNECTION CABLES. SUCH POINTS OF CONTACT WILL CAUSE SEVERE MELTING.

WE RECOMMEND THE FOLLOWING VALUES/PARAMETERS, BUT ALSO SUGGEST FOLLOWING THE MANUFACTURER'S STERILIZATION UNIT INSTRUCTIONS FOR STEAM STERILIZATION:

CYCLE	STERILIZING TEMP.	STERILIZING TIME	DRYING TIME* ³
PRE VACUUM/WRAPPED	270° F (132° C)	4 MINUTES	30 MINUTES
GRAVITY/WRAPPED	250° F (121° C)	30 MINUTES	45 MINUTES
GRAVITY/WRAPPED	270° F (132° C)	30 MINUTES	45 MINUTES

2. **ETO STERILIZATION:** INSTRUMENTS CAN BE STERILIZED BY ETHYLENE OXIDE IN ANY STANDARD CYCLE. PRESSURE READING SHOULD NOT EXCEED 12 PSI. TEMPERATURE SHOULD NOT EXCEED 68.3Â°C (155Â°F). IT IS RECOMMENDED TO FOLLOW THE MANUFACTURER'S INSTRUCTIONS FOR THE ETO STERILIZATION UNIT CONCERNING HUMIDITY, VACUUM, CYCLE TIME, GAS CONCENTRATION AND TEMPERATURE.
3. **STERRAD STERILIZATION PROCESS INCLUDING STERRAD NX:** THE STERILIZATION PROCESS IS A MULTIPLE STERILIZATION PROCESS THAT UTILIZES A COMBINATION OF EXPOSURE TO HYDROGEN PEROXIDE VAPOR AND PLASMA TO AFFECT STERILIZATION. THE STERRAD NX STERILIZER CAN STERILIZE INSTRUMENTS WHICH HAVE DIFFUSION RESTRICTED SPACES, SUCH AS HINGED PORTIONS OF FORCEPS AND SCISSORS. ADHERE TO THE STERILIZATION INSTRUCTIONS PROVIDED BY THE MANUFACTURER. (ADVANCED STERILIZATION PRODUCTS A JOHNSON & JOHNSON COMPANY)
4. **FLASH AUTOCLAVING** (FAST HEATING/COOLING CYCLE): FLASH AUTOCLAVING WILL REDUCE THE USEFUL LIFE OF THE INSTRUMENT, PARTICULARLY WHEN IT IS CONSTRUCTED OF VARIOUS MATERIALS, ENCOMPASSING DIFFERENT EXPANSION RATES.

CAUTION: CONTINUED USE OF THIS PROCESS WILL CAUSE PREMATURE FAILURE OF INSTRUMENTS, ESPECIALLY THE INSULATION.

5. **CHEMICLAVING—SOAKING:** NOT RECOMMENDED: THIS IS DESTRUCTIVE TO THE INSULATING AND SILICONE MATERIALS OF ELECTROSURGICAL ACCESSORIES AND CAN CAUSE RAPID DETERIORATION AND FAILURE.

*³ **IMPORTANT:** ADHERE TO PROPER DRYING CYCLE TO MAKE SURE THAT INSTRUMENTS ARE COMPLETELY DRY ON THE INSIDE (SAME APPLIES FOR CABLES). MOISTURE WILL PREVENT PROPER ELECTROSURGICAL CONDUCTIVITY.

WARRANTY:

- ALL ELMED INSTRUMENTS ARE WARRANTED AGAINST DEFECTS IN MATERIAL & WORKMANSHIP FOR 1 YEAR. ROUTINE RESHARPENING, REINSULATING, AND ADJUSTMENTS AS WELL AS REPAIR OF DAMAGED INSTRUMENTS IS AVAILABLE FOR A MODEST CHARGE.
- ANY SERVICING OF ELMED INSTRUMENTS PERFORMED BY ANY SOURCE OTHER THAN ELMED, INC. **VOIDS** THIS WARRANTY.
- INSTRUMENTS ARE ONLY TO BE USED AS DESIGNATED IN THEIR SPECIALIZED MEDICAL FIELDS BY THE RESPECTIVE TRAINED AND QUALIFIED PERSONNEL. THEY ARE NOT INTENDED TO BE USED ON THE CENTRAL CIRCULATORY OR CENTRAL NERVOUS SYSTEM.
- ELMED, INC., AS A MANUFACTURER AND SELLER OF THE PRODUCT, ACCEPTS NO LIABILITY FOR INDIRECT DAMAGE OR SUBSEQUENT DAMAGES THAT RESULT BY IMPROPER USAGE, IMPROPER HANDLING OR THROUGH IMPROPER PREPARATION, STERILIZATION OR MAINTENANCE.

WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS, IL
60139
USA

WE SUBSCRIBE TO COST
CONTAINMENT AND
PROTECTION OF THE
ENVIRONMENT



MADE IN THE USA

THEREFORE, WE
MANUFACTURE
REUSABLE PRODUCTS
FOR A CLEANER
WORLD

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM
WWW.ELMED.COM