




## LAPAROSCOPIC TROCAR PORT-SITE CLOSURE SYSTEM

In conjunction with leading laparoscopic surgeons, we are featuring another cost-effective, reusable, American-made instrument to help you control your spiraling operating room expenses.

***Completely Reusable – Quality American made by Elmed***



**IMPORTANT**  
**Procedure to be performed**  
**under direct vision.**

Precision Micro-Pointed and Polished Tip with Platform

### **Features:**

- Micro polished, precision-pointed, front loading, spring-action tip with platform for precise puncture and suture grasping.
- Durable stainless steel construction for long lasting durability and minimum maintenance.
- Flush port for simple and thorough cleaning.
- 3mm diameter and longer shaft makes it ideal for all laparoscopic surgeries including bariatric.
- Spring-activated handle for self closing jaw for automatic suture grasping.
- Multiple uses- The suture passer facilitates the closure of large and small trocar incision wounds, secures hernia mesh in laparoscopic vertical hernia repairs, ligates vessels. It is the reusable instrument that is needed to ensure proper placement of sutures during abdominal laparoscopy, gynecological laparoscopy and pelviscopy.

**SEE BACK PAGE FOR...**

- Complete listing of reusable Laparoscopic Trocar Port Closure Guide Cones

<b>List No. 52203-41</b>	ELMED LAPAROSCOPIC TROCAR INCISION WOUND CLOSURE & HERNIA SUTURE RETRIEVER 3.0 MM DIAMETER, 17 CM WORKING LENGTH, SPRING HANDLE, CLEANING FLUSH PORT, PRECISION MICRO-POINTED & POLISHED TIP W/ PLATFORM
<b>List No. 52203-41C</b>	ELMED LAPAROSCOPIC TROCAR INCISION WOUND CLOSURE & HERNIA SUTURE RETRIEVER 3.0 MM DIAMETER, 17 CM WORKING LENGTH, SPRING HANDLE, MODULAR 3 PIECE TAKE A PART, PRECISION MICRO-POINTED & POLISHED TIP W/ PLATFORM
<b>List No. 52213-41</b>	ELMED BARIATRIC LAPAROSCOPIC TROCAR INCISION WOUND CLOSURE & HERNIA SUTURE RETRIEVER 3.0 MM DIAMETER, 20 CM WORKING LENGTH, SPRING HANDLE, CLEANING FLUSH PORT, PRECISION MICRO-POINTED & POLISHED TIP W/ PLATFORM

Made in the USA



# LAPAROSCOPIC TROCAR PORT CLOSURE GUIDE CONES

*Completely Reusable – Quality American made by Elmed*

## Stainless Steel Trocar Guide Cones



Not Pictured – 5mm Long-52204-05 and 15mm Long-52204-15

6CM LENGTH		10CM LENGTH	
52203-05	5MM DIAMETER	52204-05	5MM DIAMETER
52203-10/11	10/11MM DIAMETER	52204-10/11	10/11MM DIAMETER
52203-12	12MM DIAMETER	52204-12	12MM DIAMETER
52203-15	15MM DIAMETER	52204-15	15MM DIAMETER

5mm – 12mm  
2-Port precision  
guide ports



15mm  
4-Port precision  
guide ports



**As A Set (Buy the Suture Retriever and Any 2 Cones, Get the Sterilization Container Free)**



### List No. 52203-41-03M

**ELMED LAPAROSCOPIC TROCAR INCISION WOUND CLOSURE & HERNIA SUTURE RETRIEVER**

3.0 MM DIAMETER, 17 CM WORKING LENGTH, SPRING HANDLE, CLEANING FLUSH PORT, PRECISION MICRO-POINTED & POLISHED TIP W/ PLATFORM PLUS ANY (2) STAINLESS STEEL TROCAR GUIDE CONES

***WE DESIGN, MANUFACTURE & SELL THE TOOLS THAT SURGEONS USE***

Made in the USA



ELMED INCORPORATED  
35 N Brandon Dr  
Glendale Heights, IL 60139  
USA

We subscribe to Cost Containment and Protection of the Environment.



Therefore, we manufacture reusable products for a cleaner world.

PH : (224) 353-6446  
FX : (224) 653-8178  
E: [medical@elmed.com](mailto:medical@elmed.com)  
[www.elmed.com](http://www.elmed.com)