



ARTHROSCOPY



WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS,
IL 60139
USA

WE SUBSCRIBE TO
COST CONTAINMENT
AND PROTECTION
OF THE
ENVIRONMENT



THEREFORE, WE
MANUFACTURE
REUSABLE PRODUCTS
FOR A CLEANER
WORLD

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM
WWW.ELMED.COM

TABLE OF CONTENTS

ARTHROSCOPY HANDLES	
NON-RATCHET HANDLES	1
RATCHET HANDLES	1
ARTHROSCOPIC HAND INSTRUMENTS	PAGE
SPECIALTY PUNCHES	2-4
SPECIALTY GRASPERS	5
SPECIALTY ROTARY PUNCHES	6
SCISSORS, PUNCHES & WIDE JAW PUNCHES	7
WIDE JAW PUNCHES, 150 UP-MILLED	8
SCISSOR PUNCHES	8
45° HOOK SCISSORS & PUNCHES	9
90° HOOK PUNCHES & WIDE JAW PUNCHES	9
GRASPING FORCEPS	9
BIRD BEAK SUTURE FORCEPS	10
SUTURE GRASPING FORCEPS	11
SUTURE HOOKS	12
ARTHROSCOPIC PUNCHES	PAGE
CUSHING	13
LOVE - GRUNWALD	13
SPURLING	13
MENISCUS, TENDON GRASPING FORCEPS	PAGE
MENISCUS GRASPING FORCEPS	13
MENISCUS RETRACTING FORCEPS	13
TENDON GRASPING FORCEPS	14
ARTHROSCOPIC INSTRUMENTS	PAGE
2-PIECE TAKE APART	15-16
1-PIECE W/ TRIANGLE HANDLE	15-16
1-PIECE ECONOMY LINE	19
ARTHROSCOPIC SCOPES & CANNULAS	PAGE
ELMED / STORZ™ STYLE	18
STRYKER™ STYLE	19
DYONICS™ STYLE	20
ARTHROSCOPIC IRRIGATION CANNULAS	PAGE
ELMED STYLE	21
DYONICSTM STYLE	21
MONOPOLAR ELECTRODES & PENCIL FOR ARTHROSCOPY	PAGE
MONOPOLAR ELECTRODES	22
HAND-SWITCHING PENCIL	22
ARTHROSCOPIC SHAVER SYSTEM	PAGE
ELMED TM 930	23
PLASTIC SHAFT SHAVER BLADES	24
STAINLESS STEEL SHAVER BLADES	25
SERVICE AND REPAIR	PAGE
SERVICE AND REPAIR	26



ARTHROSCOPY HANDLES



MADE IN THE USA

ELMED OFFERS THE MOST EXTENSIVE LINE OF HANDLE CHOICES WITH OUR ARTHROSCOPY LINE. BELOW YOU WILL FIND DIFFERENT HANDLES TO CHOOSE FROM THAT SUIT THE NEEDS OF YOUR SURGEON'S.

NON-RATCHET HANDLES



RING HANDLE (RH)



LOOP HANDLE (LH)



PRO HANDLE (PH)



WISHBONE HANDLE (WH)



CIGAR HANDLE (CH)

RATCHET HANDLES



SLIDE LOCK HANDLE (SLH)



WISHBONE W/ RATCHET HANDLE
(WHR)



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

PUNCHES



RING
HANDLE
(RH)



LOOP
HANDLE
(LH)



PRO HANDLE
(PH)



WISHBONE HANDLE
(WH)



BASKET PUNCH -
SCOOP

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
STRAIGHT	1.8 MM	23-2011RH	23-2011LH	23-2011PH	23-2011WH
STRAIGHT	2.0 MM	23-2012RH	23-2012LH	23-2012PH	23-2012WH
UPBITER	1.6 MM	23-2020RH	23-2020LH	23-2020PH	23-2020WH
STRAIGHT	1.6 MM	23-2029RH	23-2029LH	23-2029PH	23-2029WH
UPBITER	1.6 MM	23-2030RH	23-2030LH	23-2030PH	23-2030WH
UPSWEPT	1.6 MM	23-2031RH	23-2031LH	23-2031PH	23-2031WH
UPSWEPT	1.6 MM	23-2032RH	23-2032LH	23-2032PH	23-2032WH
DOWNBITER	1.6 MM	23-2033RH	23-2033LH	23-2033PH	23-2033WH
STRAIGHT	1.6 MM	23-2034RH	23-2034LH	23-2034PH	23-2034WH
CURVED LEFT	1.6 MM	23-2040RH	23-2040LH	23-2040PH	23-2040WH
CURVED RIGHT	1.6 MM	23-2041RH	23-2041LH	23-2041PH	23-2041WH
UPBITER -CURVED LEFT	1.6 MM	23-2042RH	23-2042LH	23-2042PH	23-2042WH
UPBITER -CURVED RIGHT	1.6 MM	23-2043RH	23-2043LH	23-2043PH	23-2043WH



BASKET PUNCH -
DUCKBILL

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
UPSWEPT	3.2 MM	23-2000RH	23-2000LH	23-2000PH	23-2000WH
STRAIGHT	3.2 MM	23-2013RH	23-2013LH	23-2013PH	23-2013WH
UPBITER	3.2 MM	23-2014RH	23-2014LH	23-2014PH	23-2014WH
UPBITER -CURVED LEFT	3.2 MM	23-2015RH	23-2015LH	23-2015PH	23-2015WH
UPBITER -CURVED RIGHT	3.2 MM	23-2016RH	23-2016LH	23-2016PH	23-2016WH



TO ORDER CALL: (224) 353-6446




ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

PUNCHES



BASKET PUNCH -
DUCKLING

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
STRAIGHT	3.2 MM	23-2044RH	23-2044LH	23-2044PH	23-2044WH
UPBITER	3.2 MM	23-2045RH	23-2045LH	23-2045PH	23-2045WH
CURVED LEFT	3.2 MM	23-2046RH	23-2046LH	23-2046PH	23-2046WH
CURVED RIGHT	3.2 MM	23-2047RH	23-2047LH	23-2047PH	23-2047WH
UPSWEPT	3.2 MM	23-2048RH	23-2048LH	23-2048PH	23-2048WH




BASKET PUNCH -
BLUNT NOSE

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
STRAIGHT	1.7 MM	23-2035RH	23-2035LH	23-2035PH	23-2035WH



OVAL PUNCH

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
UPBITER	4.4 MM	23-7200RH	23-7200LH	23-7200PH	23-7200PH
STRAIGHT	4.4 MM	23-7201RH	23-7201LH	23-7201PH	23-7201PH
CURVED LEFT	4.4 MM	23-7294RH	23-7294LH	23-7294PH	23-7294PH
CURVED RIGHT	4.4 MM	23-7295RH	23-7295LH	23-7295PH	23-7295PH
UPSWEPT	4.4 MM	23-7296RH	23-7296LH	23-7296PH	23-7296PH



BASKET PUNCH -
NARROWLINE

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
UPBITER	1.7 MM	23-2051RH	23-2051LH	23-2051PH	23-2051WH
STRAIGHT	1.7 MM	23-2052RH	23-2052LH	23-2052PH	23-2052WH
UPBITER-CURVED LEFT	1.7 MM	23-2053RH	23-2053LH	23-2053PH	23-2053WH
UPBITER-CURVED RIGHT	1.7 MM	23-2054RH	23-2054LH	23-2054PH	23-2054WH
UPSWEPT	1.7 MM	23-2055RH	23-2055LH	23-2055PH	23-2055WH



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

PUNCHES



straight

upbiter

upswept

curved left / right



RING
HANDLE
(RH)



LOOP
HANDLE
(LH)



PRO HANDLE
(PH)



WISHBONE HANDLE
(WH)



SCISSOR PUNCH

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
STRAIGHT	0.7 MM	23-2036RH	23-2036LH	23-2036PH	23-2036WH
UPBITER	0.7 MM	23-2037RH	23-2037LH	23-2037PH	23-2037WH
UPBITER-CURVED LEFT	0.7 MM	23-2038RH	23-2037LH	23-2037PH	23-2037WH
UPBITER- CURVED RIGHT	0.7 MM	23-2039RH	23-2038LH	23-2038PH	23-2038WH



BASKET PUNCH -
BACKBITER

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
RIGHT	2.4 MM	23-2058RH	23-2058LH	23-2058PH	23-2058WH
LEFT	2.4 MM	23-2059RH	23-2059LH	23-2059PH	23-2059WH



POSTERIOR PUNCH

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
UPBITER	2.2 MM	23-2201RH	23-2201LH	23-2201PH	23-2201WH
UPSWEPT	2.2 MM	23-2202RH	23-2202LH	23-2202PH	23-2202WH
STRAIGHT	2.2 MM	23-2203RH	23-2203LH	23-2203PH	23-2203WH
CURVED LEFT	2.2 MM	23-2204RH	23-2204LH	23-2204PH	23-2204WH
CURVED RIGHT	2.2 MM	23-2205RH	23-2205H	23-2205PH	23-2205WH



MENISCAL PUNCH -
OVERBITER

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE HANDLE
OVERBITE	—	23-9878RH	23-9878LH	23-9878PH	23-9878WH



TO ORDER CALL: (224) 353-6446




ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

GRASPING FORCEPS



STRAIGHT
UPBITER
UPSWEPT
CURVED LEFT / RIGHT



RING
HANDLE
(RH)



LOOP
HANDLE
(LH)



PRO HANDLE
(PH)



WISHBONE HANDLE
W/ RATCHET
(WH)




SLIDE LOCK
HANDLE (SLH)



GRASPING FORCEPS -
PITBULL

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE W/ RATCHET HANDLE	SLIDE LOCK HANDLE
STRAIGHT	3.4MM	23-1022RH	23-1022LH	23-1022PH	23-1022WHR	23-1022SLH



GRASPING FORCEPS
-CUPPED

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE W/ RATCHET HANDLE	SLIDE LOCK HANDLE
STRAIGHT	2.7MM	23-1012RH	23-1012LH	23-1012PH	23-1012WHR	23-1012SLH
STRAIGHT	3.4MM	23-1014RH	23-1014LH	23-1014PH	23-1014WHR	23-1014SLH



GRASPING FORCEPS -
RONGEUR

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE W/ RATCHET HANDLE	SLIDE LOCK HANDLE
RAPTOR	—	23-4720RH	23-4720LH	23-4720PH	23-4720WHR	23-4720SLH
RAPTOR JR.	—	23-4843RH	23-4843LH	23-4843PH	23-4843WHR	23-4843SLH



GRASPING FORCEPS
- ALLIGATOR

SHAFT DESIGN	BITE WIDTH	RING HANDLE	LOOP HANDLE	PRO HANDLE	WISHBONE W/ RATCHET HANDLE	SLIDE LOCK HANDLE
STRAIGHT	3.0MM	23-1020RH	23-1020LH	23-1020PH	23-1020WHR	23-1020SLH
UPCURVED	3.0MM	23-1021RH	23-1021LH	23-1021PH	23-1021WHR	23-1021SLH
STRAIGHT	3.5MM	23-1024RH	23-1024LH	23-1024PH	23-1024WHR	23-1024SLH



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ROTARY PUNCHES



STRAIGHT

UPBITER

UPSWEPT


CURVED LEFT / RIGHT



CIGAR HANDLE (CH)



LOOP HANDLE



BASKET PUNCH - 90° ROTARY

SHAFT DESIGN	TIP HEIGHT	BITE WIDTH	CIGAR HANDLE	LOOP HANDLE
90° - LEFT	3.9MM	3.4MM	23-0900CH	23-0900LH
90° - RIGHT	3.9MM	3.4MM	23-0901CH	23-0901LH
90° - LEFT	3.4MM	2.2MM	23-0906CH	23-0906LH
90° - RIGHT	3.4MM	2.2MM	23-0907CH	23-0907LH



HOOK ROTARY SCISSORS

SHAFT DESIGN	TIP HEIGHT	BITE WIDTH	CIGAR HANDLE	LOOP HANDLE
20° - LEFT	---	---	23-0800CH	23-0800LH
20° - RIGHT	---	---	23-0801CH	23-0801LH
60° - LEFT	---	---	23-0802CH	23-0802LH
60° - RIGHT	---	---	23-0803CH	23-0803LH



ROTARY SCISSORS

SHAFT DESIGN	TIP HEIGHT	BITE WIDTH	CIGAR HANDLE	LOOP HANDLE
20° SERRATED - LEFT	---	---	23-0806CH	23-0806LH
20° SERRATED - RIGHT	---	---	23-0807CH	23-0807LH



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

HOOK SCISSORS, HOOK PUNCHES & WIDE JAW PUNCHES

AVAILABLE HANDLES FOR THE
INSTRUMENTS BELOW...

PLEASE SEE EXAMPLE FOR ORDERING

23-8903XX + RING HANDLE = 23-8903RH



SHAFT DESIGN	SHAFT DIAMETER	HOOK SCISSORS	HOOK PUNCH W/ SCOOP	HOOK PUNCH, SET DOWN	HOOK PUNCH, ANGULAR
	3.4 MM	23-2503XX	23-2903XX	23-3503XX	23-3703XX
	2.7 MM	23-2502XX	23-2902XX	23-3502XX	---
	3.4 MM	23-2513XX	23-2913XX	23-3513XX	23-3713XX
	2.7 MM	23-2512XX	23-2912XX	23-3512XX	---
	3.4 MM	23-2533XX	23-2933XX	23-3533XX	23-3733XX
	2.7 MM	23-2532XX	23-2932XX	23-3532XX	---
	3.4 MM	23-2553XX	23-2953XX	23-3553XX	23-3753XX
	2.7 MM	23-2552XX	23-2952XX	23-3552XX	---
	3.4 MM	23-2573XX	23-2973XX	23-3573XX	23-3773XX
	2.7 MM	23-2572XX	23-2972XX	23-3572XX	---
	3.4 MM	23-2593XX	23-2993XX	23-3593XX	23-3793XX
	2.7 MM	23-2592XX	23-2992XX	23-3592XX	---
	3.4 MM	23-2613XX	23-3013XX	23-3613XX	23-3813XX
	2.7 MM	23-2612XX	23-3012XX	23-3612XX	---
	3.4 MM	23-2633XX	23-3033XX	23-3633XX	23-3833XX
	2.7 MM	23-2632XX	23-3032XX	23-3632XX	---
	3.4 MM	23-2653XX	23-3053XX	23-3653XX	23-3853XX
	2.7 MM	23-2652XX	23-3052XX	23-3652XX	---



SHAFT DESIGN	SHAFT DIAMETER	HOOK PUNCH, 150	WIDE JAW PUNCH, 4.5MM BITE	WIDE JAW PUNCH, 3.5MM BITE	WIDE JAW PUNCH, 2.5 MM BITE
	3.4 MM	23-3903XX	23-4903XX	23-5603XX	23-6303XX
	3.4 MM	23-3913XX	23-4913XX	23-5613XX	23-6313XX
	3.4 MM	23-3933XX	23-4933XX	23-5633XX	23-6333XX
	3.4 MM	23-3953XX	23-4953XX	23-5653XX	23-6353XX
	3.4 MM	23-3973XX	23-4973XX	23-5673XX	23-6373XX
	3.4 MM	23-3993XX	23-4993XX	23-5693XX	23-6393XX
	3.4 MM	23-4013XX	23-5013XX	23-5713XX	23-6413XX
	3.4 MM	23-4033XX	23-5033XX	23-5733XX	23-6433XX
	3.4 MM	23-4053XX	23-5053XX	23-5753XX	23-6453XX



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

AVAILABLE HANDLES FOR THE INSTRUMENTS
BELOW...

PLEASE SEE EXAMPLE FOR ORDERING
23-8903XX + RING HANDLE = 23-8903RH



WIDE JAW PUNCH - 15° UP-MILLED



	SHAFT DESIGN	SHAFT DIAMETER	WIDE JAW PUNCH, 4.5 MM	WIDE JAW PUNCH, 3.5 MM	WIDE JAW PUNCH, 2.5 MM	WIDE JAW PUNCH 45° UP, 2.0 MM
	STRAIGHT	3.4MM	23-5303XX	23-6003XX	23-6803XX	23-6983XX
	30°- RIGHT	3.4MM	23-5313XX	23-6013XX	23-6813XX	---
	30°- LEFT	3.4MM	23-5333XX	23-6033XX	23-6833XX	---
	45°- RIGHT	3.4MM	23-5353XX	23-6053XX	23-6853XX	---
	45°- LEFT	3.4MM	23-5373XX	23-6073XX	23-6873XX	---
	7°- UP	3.4MM	23-5393XX	23-6093XX	23-6893XX	---
	7°- DOWN	3.4MM	23-5413XX	23-6113XX	23-6913XX	---
	15°- UP	3.4MM	23-5433XX	23-6133XX	23-6933XX	---
	15°- DOWN	3.4MM	23-5453XX	23-6153XX	23-6953XX	---

SCISSOR PUNCHES



	SHAFT DESIGN	SHAFT DIAMETER	SCISSOR PUNCH	SCISSOR PUNCH, 15° UP-MILLED
	STRAIGHT	3.4MM	23-5303XX	23-6003XX
	30°- RIGHT	3.4MM	23-5313XX	23-6013XX
	30°- LEFT	3.4MM	23-5333XX	23-6033XX
	45°- RIGHT	3.4MM	23-5353XX	23-6053XX
	45°- LEFT	3.4MM	23-5373XX	23-6073XX
	7°- UP	3.4MM	23-5393XX	23-6093XX
	7°- DOWN	3.4MM	23-5413XX	23-6113XX
	15°- UP	3.4MM	23-5433XX	23-6133XX
	15°- DOWN	3.4MM	23-5453XX	23-6153XX



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

AVAILABLE HANDLES FOR THE INSTRUMENTS
BELOW...

PLEASE SEE EXAMPLE FOR ORDERING
23-8903XX + RING HANDLE = 23-8903RH



RING HANDLE
(RH)



LOOP HANDLE
(LH)



PRO HANDLE
(PH)



WISHBONE HANDLE
(WH)

45° HOOK SCISSORS & PUNCHES



TIP DESIGN	SHAFT DIA.	HOOK SCISSORS	HOOK PUNCH W/ SCOOP	HOOK PUNCH, BLUNT
45° RIGHT	3.4 MM	23-2703XX	23-4103XX	23-4203XX
45° LEFT	3.4 MM	23-2803XX	23-4303XX	23-4403XX

90° HOOK PUNCHES & WIDE JAW PUNCHES



TIP DESIGN	SHAFT DIA.	HOOK PUNCH W/ SCOOP	HOOK PUNCH, BLUNT	WIDE JAW PUNCH, 3.0 MM	WIDE JAW PUNCH, 1.8 MM
90° RIGHT	3.4 MM	23-4503XX	23-46503XX	23-7403XX	23-7603XX
90° LEFT	3.4 MM	23-47033XX	23-4803XX	23-7503XX	23-7703XX

GRASPING FORCEPS



SHAFT DESIGN	SHAFT DIA.	GRASPING FORCEPS	GRASPING FORCEPS 1 X 2 TEETH, FINE	GRASPING FORCEPS 1 X 2 TEETH, CRASS
STRAIGHT	2.7 MM	23-8002XX	---	23-3202XX
STRAIGHT	3.4 MM	23-8003XX	23-8203XX	23-3203XX



SHAFT DESIGN	SHAFT DIA.	GRASPING & BIOPSY FORCEPS	BULLDOG GRASPING FORCEPS	BIOPSY FORCEPS	BIOPSY FORCEPS W/ SPIKE
STRAIGHT	2.7 MM			23-3602XX	23-3602-01XX
STRAIGHT	3.4 MM	23-8603XX	---	23-3603XX	23-3603-01XX
STRAIGHT	4.5 MM	---	23-8454XX	---	---



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

AVAILABLE HANDLES FOR THE INSTRUMENTS
BELOW...

PLEASE SEE EXAMPLE FOR ORDERING
23-8903XX + **R**ING **H**ANDLE = 23-8903RH



RING HANDLE
(RH)



LOOP HANDLE
(LH)




PRO HANDLE
(PH)



WISHBONE HANDLE
(WH)

BIRD BEAK SUTURE FORCEPS



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
STRAIGHT	3.4 MM	23-9003XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
10° - UP MILLED	3.4 MM	23-8953XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
20° - UP MILLED	3.4 MM	23-3122XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
40° - UP MILLED	3.4 MM	23-3125XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
45° - UP MILLED	3.4 MM	23-8933XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
45° - LEFT	3.4 MM	23-3116XX



TIP DESIGN	SHAFT DIA.	BIRD BEAK SUTURE HOOK FORCEPS
45° - RIGHT	3.4 MM	23-3119XX



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

AVAILABLE HANDLES FOR THE INSTRUMENTS
BELOW...

PLEASE SEE EXAMPLE FOR ORDERING
23-8903XX + **R**ING **H**ANDLE = 23-8903RH




SUTURE GRASPING FORCEPS




TIP DESIGN	SHAFT DIA.	SUTURE GRASPING FORCEPS
STRAIGHT	3.4 MM	23-8933XX




TIP DESIGN	SHAFT DIA.	SERRATED SUTURE GRASPING FORCEPS
STRAIGHT	3.4 MM	23-3136XX




TIP DESIGN	SHAFT DIA.	SERRATED SUTURE GRASPING FORCEPS
45°- LEFT	3.4 MM	23-3138XX



TIP DESIGN	SHAFT DIA.	SERRATED SUTURE GRASPING FORCEPS
45°- RIGHT	3.4 MM	23-3140XX



TIP DESIGN	SHAFT DIA.	SERRATED SUTURE GRASPING FORCEPS
15°- UP MILLED	3.4 MM	23-3142XX



TIP DESIGN	SHAFT DIA.	KING SUTURE GRASPING FORCEPS
STRAIGHT	3.4 MM	23-3142XX



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

AVAILABLE HANDLES FOR THE INSTRUMENTS
BELOW...

PLEASE SEE EXAMPLE FOR ORDERING
23-8903XX + **R**ING **H**ANDLE = 23-8903RH



SUTURE HOOKS



TIP DESIGN	SHAFT DIA.	SUTURE HOOK
POINTED LEFT	3.4 MM	23-3128XX



TIP DESIGN	SHAFT DIA.	SUTURE HOOK
POINTED RIGHT	3.4 MM	23-3130XX



TIP DESIGN	SHAFT DIA.	SUTURE HOOK
TURNED LEFT	3.4 MM	23-3132XX



TIP DESIGN	SHAFT DIA.	SUTURE HOOK
TURNED RIGHT	3.4 MM	23-3134XX



TO ORDER CALL: (224) 353-6446




ARTHROSCOPIC HAND INSTRUMENTS




MADE IN THE USA

ARTHROSCOPIC PUNCHES



	WORKING LENGTH	BITE DIMENSION	STRAIGHT	45° UP CURVED	45° DOWN CURVED
CUSHING 2 X 10MM	X=180 MM	2 X 10 MM	23-0118	23-1118	23-2118
	X=150 MM	2 X 10 MM	23-0115	23-1115	23-2115
	X=120MM	2 X 10 MM	23-0112	23-1112	23-2112
LOVE - GRUWALD 3 X 10MM	X=180 MM	3 X 10 MM	23-0218	23-1218	23-2218
	X=150 MM	3 X 10 MM	23-0215	23-1215	23-2215
	X=120MM	3 X 10 MM	23-0212	23-1212	23-2212
SPURLING 4 X 10MM	X=180 MM	4 X 10 MM	23-0318	23-1318	23-2318
	X=150 MM	4 X 10 MM	23-0315	23-1315	23-2315
	X=120MM	4 X 10 MM	23-0312	23-1312	23-2312

MENISCUS GRASPING FORCEPS




WORKING LENGTH	BITE DIMENSIONS	STRAIGHT
X=120 MM	4.5 X 15 MM	23-1104

MENISCUS RETRACTING FORCEPS



STANDARD

1x2 TEETH



WORKING LENGTH	OVERALL	BITE DIMENSION	CURVED UP
85 CM	15 CM	3 X 10 MM STD	23-1105
120 MM	19 CM	4 X 15 MM STD	23-1106
85 MM	15 CM	3 X 10 MM 1x2 TEETH	23-1107
120 MM	19 CM	4 X 15 MM 1x2 TEETH	23-1108



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

TENDON GRASPING FORCEPS - STRAIGHT, RIGHT & LEFT

		(X)	
WORKING LENGTH	STRAIGHT		
X=190 MM	23-0045		

		(X)	
WORKING LENGTH	RIGHT		
X=190 MM	23-1045		

		(X)	
WORKING LENGTH	LEFT		
X=190 MM	23-2045		

CAROLL TENDON GRASPING FORCEPS

		
WORKING LENGTH	STRAIGHT	
X=190 MM	23-0101	



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC INSTRUMENTS




2 PIECE INSTRUMENT W/ ROUND HANDLE

240MM OVERALL
LENGTH



1 PIECE INSTRUMENT W/ TRIANGLE HANDLE

23-0012C



HOOK PROBE &
RETRACTOR, POINTED

23-0012A

23-0022C



KNIFE, 45° ANGLED

23-0022A


23-0032C



BAYONET KNIFE

23-0032A

23-0042C



BANANA KNIFE,
SERRATED

23-0042A


23-0052C



BANANA KNIFE, SMOOTH

23-0052A

23-0062C



SPOON, CURVED 30°,
3MM

23-0062A

23-0072C



RING CURETTE, 15°
ANGLED, 3MM

23-0072A

23-0082C



RING CURETTE, 15°
ANGLED, 7MM

23-0082A

23-0092C



HOOK PROBE,
GRADUATED, 3MM

23-0092A



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC INSTRUMENTS



2 PIECE INSTRUMENT W/ ROUND HANDLE



1 PIECE INSTRUMENT W/ TRIANGLE HANDLE

23-0102C



HOOK PROBE,
GRADUATED, 5MM

23-0102A

23-0112C



SMILLIE MENISCOTOME,
3MM

23-0112A

23-0122C



SMILLIE MENISCOTOME,
5MM

23-0122A

23-0132C



SMILLIE MENISCOTOME,
7MM

23-0132A

23-0142C



SCALPEL KNIFE

23-0142A

23-0152C



ROSETTE, STRAIGHT,
SERRATED

23-0152A

23-0162C



ROSETTE, CURVED LEFT
30°, SERRATED

23-0162A

23-0172C



ROSETTE, CURVED RIGHT
30°, SERRATED

23-0172A

23-0182C



ROSETTE, CURVED RIGHT
30°, SERRATED

23-0182A

23-0192C



HOOK KNIFE, 90°,
LATERAL CUTTING

23-0192A



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC INSTRUMENTS-ECONOMY LINE



190MM OVERALL
LENGTH



HOOK PROBE &
RETRACTOR, POINTED

23-0012E




KNIFE, 45° ANGLED

23-0022E



BAYONET KNIFE

23-0032E




BANANA KNIFE,
SERRATED

23-0042E



BANANA KNIFE,
SMOOTH

23-0052E



SPOON, CURVED 30°,
3MM

23-0062E



RING CURETTE, 15°
ANGLED, 3MM

23-0072E



RING CURETTE, 15°
ANGLED, 7MM

23-0082E



HOOK PROBE,
GRADUATED, 3MM

23-0092E



HOOK PROBE,
GRADUATED, 5MM

23-0102E



SMILLIE MENISCOTOME,
3MM

23-0112E



SMILLIE MENISCOTOME,
5MM

23-0122E



SMILLIE MENISCOTOME,
7MM

23-0132E



SCALPEL KNIFE

23-0142E



ROSETTE, STRAIGHT,
SERRATED

23-0152E



ROSETTE, CURVED LEFT
30°, SERRATED

23-0162E



ROSETTE, CURVED RIGHT
30°, SERRATED

23-0172E



SICKLE KNIFE 30°

23-0182E



HOOK KNIFE, 90°,
LATERAL CUTTING

23-0192E



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




ARTHROSCOPIC HAND INSTRUMENTS




MADE IN THE USA

ARTHROSCOPIC SCOPES & CANNULAS - ELMED / STORZ™ STYLE



ITEM NO.	DESCRIPTION
37270-00	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 0° VIEWING
37270-30	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 30° VIEWING
37270-45	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 45° VIEWING
37270-70	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 70° VIEWING

ITEM NO.	DESCRIPTION
37400-00	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 0° VIEWING
37400-30	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 30° VIEWING
37400-45	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 45° VIEWING
37400-70	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 70° VIEWING



23-1755	2.7MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STORZ™ STYLE
23-1756	4.0MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STORZ™ STYLE
23-1757	5.8MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STORZ™ STYLE
23-1758	6.5MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STORZ™ STYLE
23-1759	2.7MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STORZ™ STYLE
23-1760	4.0MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STORZ™ STYLE
23-1761	5.8MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STORZ™ STYLE
23-1762	6.5MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STORZ™ STYLE
23-1751PT	2.7MM, 175MM, PYRAMIDAL TROCAR, STORZ™ STYLE
23-1752PT	4.0MM, 175MM, PYRAMIDAL TROCAR, STORZ™ STYLE
23-1753PT	5.8MM, 175MM, PYRAMIDAL TROCAR, STORZ™ STYLE
23-1554PT	6.5MM, 175MM, PYRAMIDAL TROCAR, STORZ™ STYLE
23-1751BT	2.7MM, 175MM, BLUNT TROCAR, STORZ™ STYLE
23-1752BT	4.0MM, 175MM, BLUNT TROCAR, STORZ™ STYLE
23-1753BT	5.8MM, 175MM, BLUNT TROCAR, STORZ™ STYLE
23-1754BT	6.5MM, 175MM, BLUNT TROCAR, STORZ™ STYLE



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC SCOPES & CANNULAS - STRYKER™ STYLE



ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
37270-00-STRYKER	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 0° VIEWING, STRYKER™ STYLE	37400-00-STRYKER	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 0° VIEWING, STRYKER™ STYLE
37270-30-STRYKER	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 30° VIEWING, STRYKER™ STYLE	37400-30-STRYKER	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 30° VIEWING, STRYKER™ STYLE
37270-45-STRYKER	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 45° VIEWING, STRYKER™ STYLE	37400-45-STRYKER	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 45° VIEWING, STRYKER™ STYLE
37270-70-STRYKER	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 70° VIEWING, STRYKER™ STYLE	37400-70-STRYKER	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 70° VIEWING, STRYKER™ STYLE

23-1755-STRYKER

2.7MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STRYKER™ STYLE

23-1756-STRYKER

4.0MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STRYKER™ STYLE

23-1757-STRYKER

5.8MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STRYKER™ STYLE

23-1758-STRYKER

6.5MM, 175MM, ARTHROSCOPY SHEATH, SINGLE ROTATING STOPCOCK, STRYKER™ STYLE

23-1759-STRYKER

2.7MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STRYKER™ STYLE

23-1760-STRYKER

4.0MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STRYKER™ STYLE

23-1761-STRYKER

5.8MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STRYKER™ STYLE

23-1762-STRYKER

6.5MM, 175MM, ARTHROSCOPY SHEATH, DUAL ROTATING STOPCOCK, STRYKER™ STYLE

23-1751PT-STRYKER

2.7MM, 175MM, PYRAMIDAL TROCAR, STRYKER™ STYLE

23-1752PT-STRYKER

4.0MM, 175MM, PYRAMIDAL TROCAR, STRYKER™ STYLE

23-1753PT-STRYKER

5.8MM, 175MM, PYRAMIDAL TROCAR, STRYKER™ STYLE

23-1754PT-STRYKER

6.5MM, 175MM, PYRAMIDAL TROCAR, STRYKER™ STYLE

23-1751BT-STRYKER

2.7MM, 175MM, BLUNT TROCAR STRYKER™ STYLE

23-1752BT-STRYKER

4.0MM, 175MM, BLUNT TROCAR STRYKER™ STYLE

23-1753BT-STRYKER

5.8MM, 175MM, BLUNT TROCAR STRYKER™ STYLE

23-1754BT-STRYKER

6.5MM, 175MM, BLUNT TROCAR STRYKER™ STYLE



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC SCOPES & CANNULAS - DYONICS™ STYLE



ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
37270-00-DYONICS	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 0° VIEWING, DYONICS SCOPE LOCK	37400-00-DYONICS	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 0° VIEWING, DYONICS SCOPE LOCK
37270-30-DYONICS	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 30° VIEWING, DYONICS SCOPE LOCK	37400-30-DYONICS	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 30° VIEWING, DYONICS SCOPE LOCK
37270-45-DYONICS	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 45° VIEWING, DYONICS SCOPE LOCK	37400-45-DYONICS	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 45° VIEWING, DYONICS SCOPE LOCK
37270-70-DYONICS	ELMED EOS™ HD SINUSCOPE: 2.7MM Ø, 70° VIEWING, DYONICS SCOPE LOCK	37400-70-DYONICS	ELMED EOS™ HD SINUSCOPE: 4.0MM Ø, 70° VIEWING, DYONICS SCOPE LOCK
	23-1755-DYONICS	2.7MM, 175MM, ARTHROSCOPIC SHEATH, SINGLE ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1756-DYONICS	4.0MM, 175MM, ARTHROSCOPIC SHEATH, SINGLE ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1757-DYONICS	5.8MM, 175MM, ARTHROSCOPIC SHEATH, SINGLE ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1758-DYONICS	6.5MM, 175MM, ARTHROSCOPIC SHEATH, SINGLE ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1759-DYONICS	2.7MM, 175MM, ARTHROSCOPIC SHEATH, DUAL ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1760-DYONICS	4.0MM, 175MM, ARTHROSCOPIC SHEATH, DUAL ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1761-DYONICS	5.8MM, 175MM, ARTHROSCOPIC SHEATH, DUAL ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1762-DYONICS	6.5MM, 175MM, ARTHROSCOPIC SHEATH, DUAL ROTATING STOPCOCK, DYONICS™ SCOPE LOCK	
	23-1751PT-DYONICS	2.7MM, 175MM, CONICAL TROCAR, DYONICS™ STYLE	
	23-1752PT-DYONICS	4.0MM, 175MM, CONICAL TROCAR, DYONICS™ STYLE	
	23-1753PT-DYONICS	5.8MM, 175MM, CONICAL TROCAR, DYONICS™ STYLE	
	23-1754PT-DYONICS	6.5MM, 175MM, CONICAL TROCAR, DYONICS™ STYLE	
	23-1751BT-DYONICS	2.7MM, 175MM, SHARP TROCAR, DYONICS™ STYLE	
	23-1752BT-DYONICS	4.0MM, 175MM, SHARP TROCAR, DYONICS™ STYLE	
	23-1753BT-DYONICS	5.8MM, 175MM, SHARP TROCAR, DYONICS™ STYLE	
	23-1754BT-DYONICS	6.5MM, 175MM, SHARP TROCAR, DYONICS™ STYLE	



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




ARTHROSCOPIC HAND INSTRUMENTS



MADE IN THE USA

ARTHROSCOPIC IRRIGATION CANNULAS - ELMED STYLE



23-1032

3.2MM, 7.0CM, IRRIGATION CANNULA W/ 4 HOLES, STOPCOCK, LUER-LOCK CONNECTION



23-1045

4.5MM, 7.0CM, IRRIGATION CANNULA W/ 4 HOLES, STOPCOCK, LUER-LOCK CONNECTION

23-3032

3.2MM, 7.0CM, TROCAR SHARP

23-3045

4.5MM, 7.0CM, TROCAR SHARP

23-4032

3.2MM, 7.0CM, OBTURATOR

23-4045

3.2 MM, 7.0 CM, OBTURATOR

ARTHROSCOPIC IRRIGATION CANNULAS - DYONICS™ STYLE



23-3000

4.5MM, 8.3CM, IRRIGATION CANNULA W/ 2 HOLES, LUER-LOCK CONNECTION



23-3001

4.5MM, 8.3CM, SHARP OBTURATOR



23-3001

4.5MM, 8.3CM, BLUNT OBTURATOR



23-3003

IRRIGATION ADAPTER W/ TUBE CONNECTION



TO ORDER CALL: (224) 353-6446



MONOPOLAR ELECTRODES & PENCIL FOR ARTHROSCOPY



BLADE ELECTRODE, 45° ANGLED
No. **23-1110**



BLADE ELECTRODE, 90° ANGLED
No. **23-2110**



BALL ELECTRODE, 4MM
No. **23-3110**



SCALPEL BLADE, "CLASSIC"
No. **23-4110**



SCALPEL BLADE, "ROUND"
No. **23-5110**



5750

HAND SWITCHING HANDLE, DUAL MICRO-TOUCH SWITCH FOR CUTTING COAGULATING, ATTACHED W/10FT. CABLE W/THREE PRONG PLUG FOR VALLEYLAB AND OTHER HAND SWITCHING UNITS. FOR USE W/STANDARD ELECTRODES 3MM (3/32") DIA. SHAFT





ARTHROSCOPIC SHAVER SYSTEM



MADE IN THE USA

ELMED TM 930



TECHNICAL DATA:

DIMENSIONS (W x D x H) 323 x 267 x 118 MM

WEIGHT : 5.5 KG

VOLTAGE 90-230 VAC (50Hz/60Hz)

MAX. POWER CONSUMPTION: 100 VA

PRIMARY FUSE: 2 x T6.3A

-PROTECTION: I

-FIELD OD APPLICATION: TYPE B

-OPERATING MODE: CONTINUOUS OPERATION

-PROTECTION CLASS: IP XO

-EN-CLASS: LLA

-EMV: EN 60601-1-2

- MODEL TM 930 IS DESIGNED FOR SURGICAL INTERVENTIONS ON SMALL AND LARGE JOINTS. ERGONOMIC HAND PIECES GUARANTEE EASY WORKING WITH LESS FATIGUE. QUIET RUNNING AND LARGE SUCTION CAPACITY ARE FURTHER ATTRIBUTES OF THE UNITS
- THE MODEL TM 930 IS CLEARLY STRUCTURED FOR USER FRIENDLY APPLICATION. THE USER MAY CHANGE SETTINGS FOR ROTATION SPEED AND DIRECTION EITHER ON THE FRONT PANEL, OR VIA THE FOOTSWITCH. THE FOOTSWITCH IS NONSKID AND WATERPROOF.
- THE MODEL TM 930 OFFERS THREE CONNECTORS FOR DIFFERENT HAND PIECES, E.G. SHAVER, DRILL OR SAW. ALL HAND PIECES ARE AUTOCLAVABLE AT 1340 C.



ITEM NO.	DESCRIPTION
23-81001	MODEL TM 930 SHAVER SYSTEM, CONTROL UNIT W/ 3 CONNECTIONS FOR HAND PIECES



ITEM NO.	DESCRIPTION
23-80007	HAND PIECE FOR MODEL TM 930 SHAVER SYSTEM (1000 / 2000 / 3000 / 4000 / 5000 / 6000 / 7000 / 8000 U/MIN +/-10%)



ITEM NO.	DESCRIPTION
23-81012	FOOTSWITCH FOR MODEL 930 SHAVER SYSTEM (FUNCTIONS: RIGHT, LEFT, OSCILLATION, SPEED-CHANGE, HAND PIECE SELECTION)

MODEL TM 930 ACCESSORIES



ITEM NO.	DESCRIPTION
23-81004	DRILL, CANNULATED, AUTOCLAVABLE, FOR USE WITH MODEL TM 930 (SPEED: 150 / 300 / 600 / 900 / 1100 R/MIN +/-10%)



ITEM NO.	DESCRIPTION
23-81005	OSCILLATING SAW, AUTOCLAVABLE, W/OUT SAW BLADE FOR USE WITH MODEL TM 930 (SPEED: 20,000 OSZ/MIN +/-10%)



ITEM NO.	DESCRIPTION
23-81006	SAW BLADE FOR OSCILLATING SAW, 7MM, 50MM LENGTH



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ARTHROSCOPIC SHAVER SYSTEM



MADE IN THE USA

COMPATIBLE SHAVER BLADES IN ALL STYLES



ELMED /
LINVATEC /
ARTHREX



Dyonis Shaver PS3500 &
PS 3500EP
(ADD "A" AFTER ITEM #)



DYONIS SHAVER EP1
(ADD "B" AFTER ITEM #)



STRYKER
(ADD "C" AFTER ITEM #)

PLASTIC SHAFT SHAVER BLADES (PKG. OF 5)

AUTOCLAVABLE 134° C / 273° F



***FOR ELMED TM 930, LINVATEC/CONCEPT, & ARTHREX



SHAFT DIA.	END CUTTER	FULL RADIUS RESECTOR	AGGRESSIVE FULL RADIUS RESECTOR	MENISCUS CUTTER
3.5 MM	—	23-20330	23-20310	23-20320
4.2 MM	23-20406	23-20430	23-20410	23-20420
5.5 MM	23-20506	23-20530	23-20510	23-20520



SHAFT DIA.	AGGRESSIVE MENISCUS CUTTER	ROUND BURR	TAPERED BURR	OVAL BURR	DOUBLE AGGRESSIVE FULL RADIUS
3.5 MM	23-20304	—	—	—	—
4.2 MM	23-20404	23-20460	—	23-20480	23-20750
5.5 MM	23-20504	23-20560	23-20570	23-20580	23-20850



TO ORDER CALL: (224) 353-6446



ARTHROSCOPIC SHAVER SYSTEM



MADE IN THE USA

COMPATIBLE SHAVER BLADES IN ALL STYLES



ELMED /
LINVATEC /
ARTHREX



DYONIS SHAVER PS3500 &
PS 3500EP
(ADD "A" AFTER ITEM #)



DYONIS SHAVER EP1
(ADD "B" AFTER ITEM #)



Stryker
(ADD "C" AFTER ITEM #)

STAINLESS STEEL SHAVER BLADES (PKG. OF 5)

AUTOCLAVABLE 1340 C / 2730 F



***FOR ELMED, LINVATEC/CONCEPT, & ARTHREX



SHAFT DIA.	TRIMMER	FULL RADIUS RESECTOR	AGGRESSIVE FULL RADIUS RESECTOR	MENISCUS CUTTER
3.5 MM	23-20603	23-20630	23-20610	23-20620
4.2 MM	23-20703	23-20730	23-20710	23-20720
5.5 MM	23-20803	23-20830	23-20810	23-20820



SHAFT DIA.	AGGRESSIVE MENISCUS CUTTER	ROUND BURR	TAPERED BURR	OVAL BURR
3.5 MM	23-20605	—	—	—
4.2 MM	23-20705	23-20760	—	—
5.5 MM	23-20805	23-20860	23-20870	23-20880



ITEM NO.	DESCRIPTION
23-20200	CLEANING ADAPTER FOR SHAVING BLADES



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



SERVICE AND REPAIR



MADE IN THE USA

IN NEED OF INSTRUMENT REPAIRS? SEND THEM TO THE *ER!*

ELMED REPAIR

60 YEARS OF EXPERIENCE AND KNOWLEDGE IN DESIGNING, ENGINEERING, MANUFACTURING AND REPAIRING FOR THE MEDICAL COMMUNITY HAS GIVEN ELMED, INC. ONE OF THE BEST REPUTATIONS IN THE INDUSTRY. WE REPAIR WHAT OTHERS ARE UNABLE OR UNWILLING TO DO AND WE REPAIR ALL COMPETITIVE BRANDS. OUR WORK IS PERFORMED BY SKILLED TECHNICIANS, ALL OF WHICH HAVE BEEN WITH ELMED FOR THE PAST 10-30 YEARS.

COMPLETE IN-HOUSE MACHINE SHOP, POWDER COATING, LASER WELDING, & MARKING FINISHING, ELECTRONICS AND OPTICAL DEPARTMENTS

AREAS OF SPECIALTY

- LAPAROSCOPIC AND ENDOSCOPIC INSTRUMENTS
- ELECTROSURGICAL UNITS AND ACCESSORIES
- FIBER OPTIC LIGHT CABLES AND RETRACTORS
- RIGID ENDOSCOPES & FLEXIBLE SINUSCOPES
- STAINLESS STEEL HAND INSTRUMENTS
- LIGHT SOURCES AND HEADLIGHTS
- DEBRIDER UNITS, HAND PIECES AND BLADES
- SELF RETAINING RETRACTOR SYSTEMS

HOW TO PACKAGE AND SHIP INSTRUMENTS

- IMPORTANT:** ELMED CANNOT ACCEPT INSTRUMENTS THAT HAVE NOT BEEN PROPERLY CLEANED, DISINFECTED OR STERILIZED.
- INCLUDE CUSTOMER NAME AND ALL CONTACT INFORMATION (NAME, ADDRESS, PHONE NO., FAX NO., EMAIL ADDRESS, ETC.)
- INCLUDE CUSTOMER P.O. (PURCHASE ORDER), IF USED
- INCLUDE REASON FOR RETURN - BRIEF DESCRIPTION OF REPAIR PROBLEM IF AVAILABLE
- WRAP INSTRUMENT FOR PROTECTION DURING TRANSPORT
- USE A COMMON CARRIER SUCH AS UPS OR FEDEX
- INSURE THE PACKAGE AND RECORD THE TRACKING NUMBER

PLEASE FILL OUT, DETACH, AND INCLUDE THE FORM BELOW WITH YOUR REPAIR:

FROM:	
FACILITY:	CONTACT PERSON:
ADDRESS:	
CITY, STATE, ZIP	
EMAIL:	PHONE NO.
PO. NO.	QTY.
INSTRUMENT DESCRIPTION:	

(NOTE: PLEASE INCLUDE MODEL NO. AND SERIAL NO. IN INSTRUMENT DESCRIPTION FOR SCOPES)

FAST SERVICE: INCOMING REPAIRS COMING IN BY NEXT DAY AIR WILL BE COMPLETED WITHIN 24 HRS (WHEN PARTS ARE AVAILABLE)

SHIP TO: **ELMED** 35 N BRANDON DR., GLENDALE HEIGHTS, IL 60139