



ELMED



Clear Flush®

PATENTED & VALIDATED CLEAN TECHNOLOGY

Flushable Kerrison Rongeur



The Patented Elmed Clear Flush® Kerrison Rongeur Provides:

- *Validated* to Remove 99.9998% of Bioburden with a Simple Flush*
- *Eliminates the need to disassemble and reassemble for cleaning*
- *Designed to be Clean, Sterile and Moisture-Free for Every Patient, Every Time*
- *Increased Steam Penetration to Allow 100% Bacteria Kill*
- *Patented Internal Lubrication Channel Provides Smoother Performance*
- *IFU Independently Validated* Using FDA and AAMI Testing Protocols*

**Nelson Laboratories and Midbrook Biomedical Cleaning Validation Reports Available at www.elmed.com*

We subscribe to Cost
Containment and Protection
of the Environment.



Therefore, we manufacture
reusable products for a cleaner
world.

Elmed Incorporated
35 N Brandon Drive
Glendale Heights, IL 60139
USA

Phone - (224) 353-6446
Fax - (224) 653-8178
Email - medical@elmed.com
Web - www.elmed.com



Elmed's Validated* Clear Flush® Kerrisons Rongeurs

Elmed's patented, validated* Clear Flush® Kerrisons provide surgeons with a sharp, clean cut every time. That is because the internal bioburden that builds up over time and starts to block the cutting surface from closing completely against the footplate of the Kerrison is flushed out* on every reprocessing cycle

- The unique, exclusive design provides a smooth, tactile response for surgeons due to the removal of friction-causing bioburden combined with the complete internal lubrication of the instrument through the exclusive flush port. Elmed's Clear Flush® Kerrisons are superior in every way to the conventional Kerrisons that most surgeons trained on and still prefer over modular Kerrisons.
- The Kerrisons also provide surgeons with the distinctive ability to be flushed during a surgical case to maintain optimal performance throughout the length of the case. With the patented Elmed Clear Flush® technology, surgeons can now have the validated* confidence of a Kerrison that is not contaminated after reprocessing with bioburden, biofilm or waterborne pathogen, all of which can potentially be the source of a dangerous surgical infection..



Elmed's Clear Flush® Kerrison Being Flushed with an Enzymatic Detergent in an Automatic Washer.

- O.R. Supervisors and their personnel can now have the documented assurance* of clean, sterile, moisture-free Kerrisons for every patient, every time. Hospitals are under tremendous pressure to reduce their patients' risk of a dangerous and costly surgical infection. The exclusive Clear Flush® design represents the newest, validated technology to help fight infections caused by Kerrisons that are difficult, if not impossible, to thoroughly decontaminate, clean and sterilize.
- SPD/Sterile Reprocessing Managers can finally have the confidence* of reprocessing Kerrisons whose IFUs have been validated using FDA testing protocols. With a simple flush, the unique, patented Elmed Clear Flush® technology helps them reduce the risk of a surgical infection without the time consuming labor and effort required to manually disassemble, manually reprocess and then correctly reassemble modular Kerrisons.

*Nelson Laboratories and Midbrook Biomedical Validation Reports Available at www.aimclearflush.com

Clear Flush

40°

Aluminum Titanium Nitride (AlTiN) or Titanium Aluminum Nitride (TiAlN) is a super-hard, thin film coating developed from Titanium Nitride. This superior coating provides an extremely hard, smooth surface that reflects less light than stainless steel, which can be an inconvenience during a case. The AlTiN coating protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced surface friction helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
CC7-1MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	1mm	40 Degree
CC7-2MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	2mm	40 Degree
CC7-3MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	3mm	40 Degree
CC7-4MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	4mm	40 Degree
CC7-5MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	5mm	40 Degree
CC7-6MM40	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC8-1MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	1mm	40 Degree
CC8-2MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	2mm	40 Degree
CC8-3MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	3mm	40 Degree
CC8-4MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	4mm	40 Degree
CC8-5MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	5mm	40 Degree
CC8-6MM40	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC9-1MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	1mm	40 Degree
CC9-2MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	2mm	40 Degree
CC9-3MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	3mm	40 Degree
CC9-4MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	4mm	40 Degree
CC9-5MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	5mm	40 Degree
CC9-6MM40	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC10-1MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	1mm	40 Degree
CC10-2MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	2mm	40 Degree
CC10-3MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	3mm	40 Degree
CC10-4MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	4mm	40 Degree
CC10-5MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	5mm	40 Degree
CC10-6MM40	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC11-1MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	1mm	40 Degree
CC11-2MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	2mm	40 Degree
CC11-3MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	3mm	40 Degree
CC11-4MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	4mm	40 Degree
CC11-5MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	5mm	40 Degree
CC11-6MM40	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

40°

Aluminum Titanium Nitride (AlTiN) or Titanium Aluminum Nitride (TiAlN) is a super-hard, thin film coating developed from Titanium Nitride. This superior coating provides an extremely hard, smooth surface that reflects less light than stainless steel, which can be an inconvenience during a case. The AlTiN coating protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced surface friction helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
CC7-2MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	40 Degree
CC7-3MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	40 Degree
CC7-4MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	40 Degree
CC7-5MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	40 Degree
CC7-6MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC8-2MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	40 Degree
CC8-3MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	40 Degree
CC8-4MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	40 Degree
CC8-5MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	40 Degree
CC8-6MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC9-2MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	40 Degree
CC9-3MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	40 Degree
CC9-4MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	40 Degree
CC9-5MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	40 Degree
CC9-6MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC10-2MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	40 Degree
CC10-3MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	40 Degree
CC10-4MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	40 Degree
CC10-5MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	40 Degree
CC10-6MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
CC11-2MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	40 Degree
CC11-3MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	40 Degree
CC11-4MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	40 Degree
CC11-5MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	40 Degree
CC11-6MM40-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

40°

Titanium Nitride (TiN) is a very hard, inert, thin film coating that is applied on the highest quality surgical instruments. TiN has an ideal combination of hardness, toughness and superior adhesion. TiN's inert surface barrier protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced friction coefficient provided by the coating helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
SC7-1MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	1mm	40 Degree
SC7-2MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	2mm	40 Degree
SC7-3MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	3mm	40 Degree
SC7-4MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	4mm	40 Degree
SC7-5MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	5mm	40 Degree
SC7-6MM40	Super Coat Kerrison Rongeur	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC8-1MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	1mm	40 Degree
SC8-2MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	2mm	40 Degree
SC8-3MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	3mm	40 Degree
SC8-4MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	4mm	40 Degree
SC8-5MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	5mm	40 Degree
SC8-6MM40	Super Coat Kerrison Rongeur	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC9-1MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	1mm	40 Degree
SC9-2MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	2mm	40 Degree
SC9-3MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	3mm	40 Degree
SC9-4MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	4mm	40 Degree
SC9-5MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	5mm	40 Degree
SC9-6MM40	Super Coat Kerrison Rongeur	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC10-1MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	1mm	40 Degree
SC10-2MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	2mm	40 Degree
SC10-3MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	3mm	40 Degree
SC10-4MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	4mm	40 Degree
SC10-5MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	5mm	40 Degree
SC10-6MM40	Super Coat Kerrison Rongeur	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC11-1MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	1mm	40 Degree
SC11-2MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	2mm	40 Degree
SC11-3MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	3mm	40 Degree
SC11-4MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	4mm	40 Degree
SC11-5MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	5mm	40 Degree
SC11-6MM40	Super Coat Kerrison Rongeur	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

40°

Titanium Nitride (TiN) is a very hard, inert, thin film coating that is applied on the highest quality surgical instruments. TiN has an ideal combination of hardness, toughness and superior adhesion. TiN's inert surface barrier protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced friction coefficient provided by the coating helps to prevent tissue from adhering to the Kerrison.



w/ Pin Ejector



Item #	Description	Style	Length	Bite	Angle
SC7-2MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	40 Degree
SC7-3MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	40 Degree
SC7-4MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	40 Degree
SC7-5MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	40 Degree
SC7-6MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC8-2MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	40 Degree
SC8-3MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	40 Degree
SC8-4MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	40 Degree
SC8-5MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	40 Degree
SC8-6MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC9-2MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	40 Degree
SC9-3MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	40 Degree
SC9-4MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	40 Degree
SC9-5MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	40 Degree
SC9-6MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC10-2MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	40 Degree
SC10-3MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	40 Degree
SC10-4MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	40 Degree
SC10-5MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	40 Degree
SC10-6MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SC11-2MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	40 Degree
SC11-3MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	40 Degree
SC11-4MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	40 Degree
SC11-5MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	40 Degree
SC11-6MM40-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

40°

The original Clear Flush® Kerrisons are manufactured with the highest quality, surgical grade German Stainless Steel. Our traditional Clear Flush® Kerrisons remain extremely popular with surgeons due to their unique, patented ability to be flushed during a case.



Item #	Description	Style	Length	Bite	Angle
SS7-1MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	1mm	40 Degree
SS7-2MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	2mm	40 Degree
SS7-3MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	3mm	40 Degree
SS7-4MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	4mm	40 Degree
SS7-5MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	5mm	40 Degree
SS7-6MM40	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS8-1MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	1mm	40 Degree
SS8-2MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	2mm	40 Degree
SS8-3MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	3mm	40 Degree
SS8-4MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	4mm	40 Degree
SS8-5MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	5mm	40 Degree
SS8-6MM40	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS9-1MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	1mm	40 Degree
SS9-2MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	2mm	40 Degree
SS9-3MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	3mm	40 Degree
SS9-4MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	4mm	40 Degree
SS9-5MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	5mm	40 Degree
SS9-6MM40	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS10-1MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	1mm	40 Degree
SS10-2MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	2mm	40 Degree
SS10-3MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	3mm	40 Degree
SS10-4MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	4mm	40 Degree
SS10-5MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	5mm	40 Degree
SS10-6MM40	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS11-1MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	1mm	40 Degree
SS11-2MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	2mm	40 Degree
SS11-3MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	3mm	40 Degree
SS11-4MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	4mm	40 Degree
SS11-5MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	5mm	40 Degree
SS11-6MM40	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

40°

The original Clear Flush® Kerrisons are manufactured with the highest quality, surgical grade German Stainless Steel. Our traditional Clear Flush® Kerrisons remain extremely popular with surgeons due to their unique, patented ability to be flushed during a case.



Item #	Description	Style	Length	Bite	Angle
SS7-2MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	40 Degree
SS7-3MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	40 Degree
SS7-4MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	40 Degree
SS7-5MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	40 Degree
SS7-6MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS8-2MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	40 Degree
SS8-3MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	40 Degree
SS8-4MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	40 Degree
SS8-5MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	40 Degree
SS8-6MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS9-2MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	40 Degree
SS9-3MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	40 Degree
SS9-4MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	40 Degree
SS9-5MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	40 Degree
SS9-6MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS10-2MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	40 Degree
SS10-3MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	40 Degree
SS10-4MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	40 Degree
SS10-5MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	40 Degree
SS10-6MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	40 Degree

Item #	Description	Style	Length	Bite	Angle
SS11-2MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	40 Degree
SS11-3MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	40 Degree
SS11-4MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	40 Degree
SS11-5MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	40 Degree
SS11-6MM40-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	40 Degree

Clear Flush

90°

Aluminum Titanium Nitride (AlTiN) or Titanium Aluminum Nitride (TiAlN) is a super-hard, thin film coating developed from Titanium Nitride. This superior coating provides an extremely hard, smooth surface that reflects less light than stainless steel, which can be an inconvenience during a case. The AlTiN coating protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced surface friction helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
CC7-1MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	1mm	90 Degree
CC7-2MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	2mm	90 Degree
CC7-3MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	3mm	90 Degree
CC7-4MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	4mm	90 Degree
CC7-5MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	5mm	90 Degree
CC7-6MM90	Ceramic Coat Kerrison Rongeur	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC8-1MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	1mm	90 Degree
CC8-2MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	2mm	90 Degree
CC8-3MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	3mm	90 Degree
CC8-4MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	4mm	90 Degree
CC8-5MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	5mm	90 Degree
CC8-6MM90	Ceramic Coat Kerrison Rongeur	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC9-1MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	1mm	90 Degree
CC9-2MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	2mm	90 Degree
CC9-3MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	3mm	90 Degree
CC9-4MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	4mm	90 Degree
CC9-5MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	5mm	90 Degree
CC9-6MM90	Ceramic Coat Kerrison Rongeur	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC10-1MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	1mm	90 Degree
CC10-2MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	2mm	90 Degree
CC10-3MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	3mm	90 Degree
CC10-4MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	4mm	90 Degree
CC10-5MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	5mm	90 Degree
CC10-6MM90	Ceramic Coat Kerrison Rongeur	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC11-1MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	1mm	90 Degree
CC11-2MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	2mm	90 Degree
CC11-3MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	3mm	90 Degree
CC11-4MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	4mm	90 Degree
CC11-5MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	5mm	90 Degree
CC11-6MM90	Ceramic Coat Kerrison Rongeur	Micro	11" (300mm)	6mm	90 Degree

Clear Flush

90°

Aluminum Titanium Nitride (AlTiN) or Titanium Aluminum Nitride (TiAlN) is a super-hard, thin film coating developed from Titanium Nitride. This superior coating provides an extremely hard, smooth surface that reflects less light than stainless steel, which can be an inconvenience during a case. The AlTiN coating protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced surface friction helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
CC7-2MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	90 Degree
CC7-3MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	90 Degree
CC7-4MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	90 Degree
CC7-5MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	90 Degree
CC7-6MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC8-2MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	90 Degree
CC8-3MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	90 Degree
CC8-4MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	90 Degree
CC8-5MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	90 Degree
CC8-6MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC9-2MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	90 Degree
CC9-3MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	90 Degree
CC9-4MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	90 Degree
CC9-5MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	90 Degree
CC9-6MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC10-2MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	90 Degree
CC10-3MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	90 Degree
CC10-4MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	90 Degree
CC10-5MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	90 Degree
CC10-6MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
CC11-2MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	90 Degree
CC11-3MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	90 Degree
CC11-4MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	90 Degree
CC11-5MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	90 Degree
CC11-6MM90-PIN	Ceramic Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	90 Degree

Clear Flush

90°

Titanium Nitride (TiN) is a very hard, inert, thin film coating that is applied on the highest quality surgical instruments. TiN has an ideal combination of hardness, toughness and superior adhesion. TiN's inert surface barrier protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced friction coefficient provided by the coating helps to prevent tissue from adhering to the Kerrison.



Item #	Description	Style	Length	Bite	Angle
SC7-1MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	1mm	90 Degree
SC7-2MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	2mm	90 Degree
SC7-3MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	3mm	90 Degree
SC7-4MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	4mm	90 Degree
SC7-5MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	5mm	90 Degree
SC7-6MM90	Super Coat Kerrison Rongeur	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC8-1MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	1mm	90 Degree
SC8-2MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	2mm	90 Degree
SC8-3MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	3mm	90 Degree
SC8-4MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	4mm	90 Degree
SC8-5MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	5mm	90 Degree
SC8-6MM90	Super Coat Kerrison Rongeur	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC9-1MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	1mm	90 Degree
SC9-2MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	2mm	90 Degree
SC9-3MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	3mm	90 Degree
SC9-4MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	4mm	90 Degree
SC9-5MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	5mm	90 Degree
SC9-6MM90	Super Coat Kerrison Rongeur	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC10-1MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	1mm	90 Degree
SC10-2MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	2mm	90 Degree
SC10-3MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	3mm	90 Degree
SC10-4MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	4mm	90 Degree
SC10-5MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	5mm	90 Degree
SC10-6MM90	Super Coat Kerrison Rongeur	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC11-1MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	1mm	90 Degree
SC11-2MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	2mm	90 Degree
SC11-3MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	3mm	90 Degree
SC11-4MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	4mm	90 Degree
SC11-5MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	5mm	90 Degree
SC11-6MM90	Super Coat Kerrison Rongeur	Micro	11" (300mm)	6mm	90 Degree

Clear Flush

90°

Titanium Nitride (TiN) is a very hard, inert, thin film coating that is applied on the highest quality surgical instruments. TiN has an ideal combination of hardness, toughness and superior adhesion. TiN's inert surface barrier protects the Kerrison from corrosion, improves wear resistance and maintains the integrity of the cutting edge longer. The reduced friction coefficient provided by the coating helps to prevent tissue from adhering to the Kerrison.



w/ Pin Ejector

Item #	Description	Style	Length	Bite	Angle
SC7-2MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	90 Degree
SC7-3MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	90 Degree
SC7-4MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	90 Degree
SC7-5MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	90 Degree
SC7-6MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC8-2MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	90 Degree
SC8-3MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	90 Degree
SC8-4MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	90 Degree
SC8-5MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	90 Degree
SC8-6MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC9-2MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	90 Degree
SC9-3MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	90 Degree
SC9-4MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	90 Degree
SC9-5MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	90 Degree
SC9-6MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC10-2MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	90 Degree
SC10-3MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	90 Degree
SC10-4MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	90 Degree
SC10-5MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	90 Degree
SC10-6MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SC11-2MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	90 Degree
SC11-3MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	90 Degree
SC11-4MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	90 Degree
SC11-5MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	90 Degree
SC11-6MM90-PIN	Super Coat Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	90 Degree



The original Clear Flush® Kerrisons are manufactured with the highest quality, surgical grade German Stainless Steel. Our traditional Clear Flush® Kerrisons remain extremely popular with surgeons due to their unique, patented ability to be flushed during a case.

Item #	Description	Style	Length	Bite	Angle
SS7-1MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	1mm	90 Degree
SS7-2MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	2mm	90 Degree
SS7-3MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	3mm	90 Degree
SS7-4MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	4mm	90 Degree
SS7-5MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	5mm	90 Degree
SS7-6MM90	Stainless Steel Kerrison Rongeur	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS8-1MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	1mm	90 Degree
SS8-2MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	2mm	90 Degree
SS8-3MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	3mm	90 Degree
SS8-4MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	4mm	90 Degree
SS8-5MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	5mm	90 Degree
SS8-6MM90	Stainless Steel Kerrison Rongeur	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS9-1MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	1mm	90 Degree
SS9-2MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	2mm	90 Degree
SS9-3MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	3mm	90 Degree
SS9-4MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	4mm	90 Degree
SS9-5MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	5mm	90 Degree
SS9-6MM90	Stainless Steel Kerrison Rongeur	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS10-1MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	1mm	90 Degree
SS10-2MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	2mm	90 Degree
SS10-3MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	3mm	90 Degree
SS10-4MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	4mm	90 Degree
SS10-5MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	5mm	90 Degree
SS10-6MM90	Stainless Steel Kerrison Rongeur	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS11-1MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	1mm	90 Degree
SS11-2MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	2mm	90 Degree
SS11-3MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	3mm	90 Degree
SS11-4MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	4mm	90 Degree
SS11-5MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	5mm	90 Degree
SS11-6MM90	Stainless Steel Kerrison Rongeur	Micro	11" (300mm)	6mm	90 Degree



The original Clear Flush® Kerrisons are manufactured with the highest quality, surgical grade German Stainless Steel. Our traditional Clear Flush® Kerrisons remain extremely popular with surgeons due to their unique, patented ability to be flushed during a case.



Item #	Description	Style	Length	Bite	Angle
SS7-2MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	2mm	90 Degree
SS7-3MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	3mm	90 Degree
SS7-4MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	4mm	90 Degree
SS7-5MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	5mm	90 Degree
SS7-6MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	7" (175mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS8-2MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	2mm	90 Degree
SS8-3MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	3mm	90 Degree
SS8-4MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	4mm	90 Degree
SS8-5MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	5mm	90 Degree
SS8-6MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	8" (200mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS9-2MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	2mm	90 Degree
SS9-3MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	3mm	90 Degree
SS9-4MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	4mm	90 Degree
SS9-5MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	5mm	90 Degree
SS9-6MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	9" (225mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS10-2MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	2mm	90 Degree
SS10-3MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	3mm	90 Degree
SS10-4MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	4mm	90 Degree
SS10-5MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	5mm	90 Degree
SS10-6MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	10" (250mm)	6mm	90 Degree

Item #	Description	Style	Length	Bite	Angle
SS11-2MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	2mm	90 Degree
SS11-3MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	3mm	90 Degree
SS11-4MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	4mm	90 Degree
SS11-5MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	5mm	90 Degree
SS11-6MM90-PIN	Stainless Steel Kerrison Rongeur w/ Pin Ejector	Micro	11" (300mm)	6mm	90 Degree