

RECTAL & GASTROENTEROLOGY



WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS,
IL 60139
USA

WE SUBSCRIBE TO
COST CONTAINMENT
AND PROTECTION
OF THE
ENVIRONMENT



THEREFORE, WE
MANUFACTURE
REUSABLE
PRODUCTS FOR A
CLEANER WORLD

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM

TABLE OF CONTENTS

	PAGE
ENDOSCOPIC FLEXIBLE SPOT COAGULATION ELECTRODE REUSABLE “OLYMPUS SNARE”	1
REUSABLE POLYPECTOMY SNARES & FLEXIBLE HOT BIOPSY FORCEPS	2
ELECTROSURGICAL SNARES	3
SUCTION COAGULATORS	4
FLEXIBLE BIOPSY INSTRUMENTS	5-6
PRIMAX RECTAL INSTRUMENTS	7-8
RECTAL BIOPSY INSTRUMENTS	9
SERVICE AND REPAIR	10



ELMED ENDOSCOPIC FLEXIBLE SPOT COAGULATION ELECTRODE REUSABLE “OLYMPUS STYLE”



MADE IN THE USA

FEATURES OF THE FLEXIBLE SPOT COAGULATOR:

- REUSABLE AND AUTOCLAVABLE
- HIGHEST GRADE STAINLESS STEEL WIRE STRAND AND TIP
- TEFLON COATED OUTER JACKET
- CONNECTED WITH A MACHINED PLUG
- FLEXIBLE AND STRONG, PREVENTING THE WIRE FROM KINKING AFTER MULTIPLE USES



FLEXIBLE SPOT COAGULATOR

30005-230	REUSABLE, 2.3MM DIA., 230CM LONG
30005-160	REUSABLE, 2.3MM DIA., 160CM LONG
30005-110	REUSABLE, 2.3MM DIA., 110CM LONG



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ELMED REUSABLE POLYPECTOMY SNARES & FLEXIBLE HOT BIOPSY FORCEPS



MADE IN THE USA



30001-160
30001-230



30003-160
30003-230



30024-160
30024-230



30005-15

SNARES AND BIOPSY FORCEPS

30001-160	FLEXIBLE POLYPECTOMY SNARE: 15MM CRESCENT LOOP, 160CM LONG
30001-230	FLEXIBLE POLYPECTOMY SNARE: 15MM CRESCENT LOOP, 230CM LONG
30002-160	FLEXIBLE POLYPECTOMY SNARE: 25MM CRESCENT LOOP, 160CM LONG
30002-230	FLEXIBLE POLYPECTOMY SNARE: 25MM CRESCENT LOOP, 230CM LONG
30003-160	FLEXIBLE POLYPECTOMY SNARE: 25MM OVAL LOOP, 160CM LONG
30003-230	FLEXIBLE POLYPECTOMY SNARE: 25MM OVAL LOOP, 230CM LONG
30024-160	HOT BIPSY FORCEPS: ROUND CUP, 160CM LONG
30024-230	HOT BIPSY FORCEPS: ROUND CUP, 230CM LONG
30005-15	HANDLE: FOR ELMED/OLYMPUS SNARES & FLEXIBLE HOT BIOPSY FORCEPS

MONOPOLAR CABLES

- FULLY AUTOCLAVABLE
- 10' LONG



30006-300: MONOPOLAR RF-CONNECTION CABLE: FOR OLYMPUS HANDLE & UNITS W/ MALE BANANA PLUG



30007-300: MONOPOLAR RF-CONNECTION CABLE: FOR OLYMPUS HANDLE & UNITS W/ MALE BOVIE PLUG



TO ORDER CALL: (224) 353-6446



ELECTROSURGICAL SNARES



MADE IN THE USA

ELECTROSURGICAL SNARE REMOVAL OF POLYPS IS AN ACCEPTED AND SAFE TECHNIQUE, OFFERING THE FOLLOWING ADVANTAGES:

- STERILIZES THE WOUND AS IT CUTS, SEALING LYMPHATIC VESSELS AND CAPILLARIES, PREVENTING ABSORPTION OF TOXIC SUBSTANCES INTO THE BLOOD STREAM.
- SIMPLE HEMOSTASIS, REDUCED BLOOD LOSS, REDUCING THE RISK OF POSTOPERATIVE BLEEDING.
- REDUCES OPERATING TIME.
- CUTS EFFORTLESSLY, COAGULATES TISSUE SURFACE, PERMITTING FIRST INTENTION HEALING.



ITEM #	DESCRIPTION
5281	RECTAL SUCTION SNARE: 19" LONG, INSULATED SHAFT, COMPLETE WITH SELF OPENING LOOP.
5282	NASAL SUCTION SNARE: 4" LONG, INSULATED SHAFT, COMPLETE WITH SELF OPENING LOOP.
5283	TONSIL SUCTION SNARE: 9" LONG, INSULATED SHAFT, COMPLETE WITH SELF OPENING LOOP
5284	SELF OPENING LOOP FOR RECTAL SNARE, 19" LONG
5285	SELF OPENING LOOP FOR NASAL SNARE, 9" LONG
5286	SELF OPENING LOOP FOR NASAL SNARE, 4" LONG

ELECTROSURGICAL ACCESSORIES FOR PROCTOLOGY

FOR REMOVAL OF SMALL POLYPS AND BIOPSIES, LOOP ELECTRODES ARE AVAILABLE. FOR THE TREATMENT OF RECTAL BLEEDINGS, BALL ELECTRODES ARE OFFERED. THE FOLLOWING ELECTRODES FIT OUR STANDARD ELECTROSURGICAL HANDLE ITEM # 5014 AND ALL DISPOSABLE HANDLES.



ITEM #	DESCRIPTION
5014	HANDLE FOR FOOTSWITCH OPERATION, FOR ELECTRODES WITH 3/32" SHAFT
5191	WIRE LOOP ELECTRODE, INSULATED, 12" LONG, 10MM Ø, 3/32" SHAFT
5193	WIRE LOOP ELECTRODE, INSULATED, 12" LONG, 6MM Ø, 3/32" SHAFT
5194	BALL ELECTRODE, INSULATED, 12" LONG, 5MM Ø, 3/32" SHAFT
5197	BALL ELECTRODE, INSULATED, 12" LONG, 3MM Ø, 3/32" SHAFT



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




SUCTION COAGULATORS



MADE IN THE USA

FOR THE PURPOSE OF REMOVING BLOOD, MUCUS AND SMOKE FROM A BODY CAVITY AND COAGULATE IT AT THE SAME TIME, WE DEVELOPED A LINE OF SPECIALTY INSTRUMENTS. A WIDE RANGE OF SIZES AND STYLES ARE AVAILABLE FOR THE DIFFERENT SPECIALTIES.



PISTOL GRIP HANDLE W/ DETACHABLE EXTENSIONS. INSTRUMENTS ARE EASILY CLEANED AND STERILIZED BY CONVENTIONAL AUTOCLAVING

5239	HANDLE, PISTOL GRIP FOR SUCTION COAGULATOR FOR SINGLE BUTTON HANDSWITCH OPERATION, WITH RF-CABLE AND CONNECTOR FOR 3-PRONG OUTLET
5241	HANDLE, PISTOL GRIP FOR SUCTION COAGULATOR FOR FOOTSWITCH OPERATION
5245	HANDLE, PISTOL GRIP FOR SUCTION COAGULATOR WITH FINGERTIP SWITCH AND RF-CABLE FOR CONNECTION TO ELMED GENERATORS ONLY
5242	SUCTION COAGULATOR EXTENSION, NASAL 4" LONG
5244	SUCTION COAGULATOR EXTENSION,, TONSILLAR, 9" LONG
5246	SUCTION COAGULATOR EXTENSION, RECTAL, 13½" LONG

MONOPOLAR RF-CONNECTION CABLES

5021	RF-CABLE, SINGLE CONDUCTOR, MALE BANANA PLUG ON EACH END
5009	RF-CABLE, SINGLE CONDUCTOR, MALE BANANA PLUG ON ONE END AND MALE BOVIE PLUG ON THE OTHER

BIPOLAR SUCTION COAGULATORS

5980	BIPOLAR COAGULATORS, WITH IRRIGATION CHANNEL
5982	BIPOLAR COAGULATORS, WITH SUCTION AND IRRIGATION CHANNEL



TO ORDER CALL: (224) 353-6446



FLEXIBLE BIOPSY INSTRUMENTS



MADE IN THE USA



30020-160
30020-230



30021-160
30021-230



30022-160
30022-230



30023-160
30023-230



30024-160
30024-230

ITEM NO.	DESCRIPTION
30020-160	FORCEPS WITH ROUND CUP, 160CM LONG
30020-230	FORCEPS WITH ROUND CUP, 230CM LONG
30021-160	FORCEPS WITH ROUND CUP AND NEEDLE, 160CM LONG
30021-230	FORCEPS WITH ROUND CUP AND NEEDLE, 230CM LONG
30022-160	FORCEPS WITH OVAL CUP, 160CM LONG
30022-230	FORCEPS WITH OVAL CUP, 230CM LONG
30023-160	FORCEPS WITH OVAL CUP AND NEEDLE, 160CM LONG
30023-230	FORCEPS WITH OVAL CUP AND NEEDLE, 230CM LONG
30024-160	HOT BIOPSY FORCEPS, ROUND CUP, 160CM LONG
30024-230	HOT BIOPSY FORCEPS, ROUND CUP, 230CM LONG



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



FLEXIBLE BIOPSY INSTRUMENTS



MADE IN THE USA



30031-160
30031-230



30033-160
30033-230



30036-160
30036-230



30034-160
30034-230



30035-160
30035-230

LIST NO.	DESCRIPTION
30031-160	SCISSORS WITH ONE FINE SERRATED BLADE, 160CM LONG
30031-230	SCISSORS WITH ONE FINE SERRATED BLADE, 230CM LONG
30032-160	TENACULUM FORCEPS, 1X2 TEETH, 160 LONG
30032-230	TENACULUM FORCEPS, 1X2 TEETH, 230 LONG
30033-160	TENACULUM FORCEPS, 2X2 TEETH, 160 LONG
30033-230	TENACULUM FORCEPS, 2X2 TEETH, 230 LONG
30034-160	GRASPING FORCEPS, 3 PRONGS, 160CM LONG
30034-230	GRASPING FORCEPS, 3 PRONGS, 230CM LONG
30035-160	POLYPUS FORCEPS, 160CM LONG
30035-230	POLYPUS FORCEPS, 230CM LONG
30036-160	TENACULUM FORCEPS, RAT TOOTHED, 160CM LONG
30036-230	TENACULUM FORCEPS, RAT TOOTHED, 230CM LONG



TO ORDER CALL: (224) 353-6446



PRIMAX RECTAL INSTRUMENTS



MADE IN THE USA



CZERNY RECTAL SPECULUM

30-34-21	8¾", 22CM
-----------------	-----------




SIMS RECTAL SPECULUM

30-35-15	6", 15CM
30-36-15	FENESTRATED, 6", 15CM



PRATT RECTAL SPECULUM

30-37-16	8½" (21.6CM) LONG, 3½" BLADES
-----------------	-------------------------------------



SMITH-BUIE ANAL RETRACTOR

30-38-15	6", 15CM
-----------------	----------



BARR RECTAL SPECULUM

30-39-17	6¾", 18CM
-----------------	-----------



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



PRIMAX RECTAL INSTRUMENTS



MADE IN THE USA



IVES-FANSLER RECTAL SPECULUM

30-40-15	1 x 3 1/4"
-----------------	------------



HIRSCHMAN ANAOSCOPE

30-41-01	5/8" 2 3/8"
30-41-02	3/4" x 2 3/8"
30-41-03	7/8" x 2 3/8"



BRINKERHOFF RECTAL SPECULUM

30-42-01	7/8" x 3 1/2"
30-42-02	1 x 4"
30-42-03	1 1/8" x 4 3/8"



KELLY PROCTOSCOPE

30-43-01	7/8 x 2"
30-43-02	7/8 x 4"
30-43-03	7/8" x 5 1/2"
30-43-04	7/8 x 7 7/8"




TO ORDER CALL: (224) 353-6446



RECTAL BIOPSY INSTRUMENTS




MADE IN THE USA



YEOMAN, FIXED HANDLE BIOPSY FORCEPS

30-45-42 | 16½", 42CM



YEOMAN, ROTATING BIOPSY SHAFT AND HANDLE

30-46-46 | 18", 46CM





SERVICE AND REPAIR



MADE IN THE USA

IN NEED OF INSTRUMENT REPAIRS? SEND THEM TO THE *ER*!

ELMED REPAIR

60 YEARS OF EXPERIENCE AND KNOWLEDGE IN DESIGNING, ENGINEERING, MANUFACTURING AND REPAIRING FOR THE MEDICAL COMMUNITY HAS GIVEN ELMED, INC. ONE OF THE BEST REPUTATIONS IN THE INDUSTRY. WE REPAIR WHAT OTHERS ARE UNABLE OR UNWILLING TO DO AND WE REPAIR ALL COMPETITIVE BRANDS. OUR WORK IS PERFORMED BY SKILLED TECHNICIANS, ALL OF WHICH HAVE BEEN WITH ELMED FOR THE PAST 10-30 YEARS.

COMPLETE IN-HOUSE MACHINE SHOP, POWDER COATING, LASER WELDING, & MARKING FINISHING, ELECTRONICS AND OPTICAL DEPARTMENTS

AREAS OF SPECIALTY

- LAPAROSCOPIC AND ENDOSCOPIC INSTRUMENTS
- ELECTROSURGICAL UNITS AND ACCESSORIES
- FIBER OPTIC LIGHT CABLES AND RETRACTORS
- RIGID ENDOSCOPES & FLEXIBLE SINUSCOPES
- STAINLESS STEEL HAND INSTRUMENTS
- LIGHT SOURCES AND HEADLIGHTS
- DEBRIDER UNITS, HAND PIECES AND BLADES
- SELF RETAINING RETRACTOR SYSTEMS

HOW TO PACKAGE AND SHIP INSTRUMENTS

- IMPORTANT:** ELMED CANNOT ACCEPT INSTRUMENTS THAT HAVE NOT BEEN PROPERLY CLEANED, DISINFECTED OR STERILIZED.
- INCLUDE CUSTOMER NAME AND ALL CONTACT INFORMATION (NAME, ADDRESS, PHONE NO., FAX NO., EMAIL ADDRESS, ETC.)
- INCLUDE CUSTOMER P.O. (PURCHASE ORDER), IF USED
- INCLUDE REASON FOR RETURN - BRIEF DESCRIPTION OF REPAIR PROBLEM IF AVAILABLE
- WRAP INSTRUMENT FOR PROTECTION DURING TRANSPORT

PLEASE FILL OUT, DETACH, AND INCLUDE THE FORM BELOW WITH YOUR REPAIR:

FROM:	
FACILITY:	CONTACT PERSON:
ADDRESS:	
CITY, STATE, ZIP	
EMAIL:	PHONE NO.
PO. No.	QTY.
INSTRUMENT DESCRIPTION:	

(NOTE: PLEASE INCLUDE MODEL NO. AND SERIAL NO. IN INSTRUMENT DESCRIPTION FOR SCOPES)

FAST SERVICE: INCOMING REPAIRS COMING IN BY NEXT DAY AIR WILL BE COMPLETED WITHIN 24 HRS (WHEN PARTS ARE AVAILABLE)

SHIP TO: **ELMED** 35 N BRANDON DR., GLENDALE HEIGHTS, IL 60139