



ELECTROSURGERY

WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS, IL
60139

**WE SUBSCRIBE TO
COST CONTAINMENT
AND PROTECTION OF
THE ENVIRONMENT**



MADE IN THE USA

**THEREFORE, WE
MANUFACTURE
REUSABLE PRODUCTS
FOR A CLEANER
WORLD**

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM

TABLE OF CONTENTS

	PAGE
BIPOLAR FORCEPS FEATURES	1
BIPOLAR BAYONET FORCEPS	2-3
BIPOLAR STRAIGHT FORCEPS	4-5
BIPOLAR BAYONET FORCEPS WITH IRRIGATION	6-8
BIPOLAR STRAIGHT FORCEPS WITH IRRIGATION	8-9
SPECIAL MICRO & OPHTHALMIC BIPOLAR FORCEPS	9-10
BIPOLAR OUTPUT TESTER	10
BIPOLAR CABLES	10





ELMED REUSABLE BIPOLAR SURGICAL FORCEPS



MADE IN THE USA


Quality, value priced reusable American made surgical instruments by ELMED.

When precise atraumatic coagulation is essential, bipolar coagulation forceps play an important role. **Elmed's goal is to provide you with the most cost-effective, quality reusable, value priced instrument to help you reduce and control your spiraling surgical procedure costs.**

Elmed offers the largest selection of reusable, fully autoclavable bipolar surgical forceps for all surgical procedures, such as Neuro, ENT, Plastic, Reconstructive, Hand Microsurgery and Ophthalmic surgery.

ELMED BIPOLAR FORCEPS FEATURES:


- Exclusive Elmed Electro-Micro Polished tip for maximum non-stick coagulation at the proper power level setting
- Highest Grade medical Stainless Steel Construction
- Insulated forceps feature medical grade polymer coating applied electrostatically for maximum durability, no RF leakage, and precise design
- Non-insulated forceps feature a low sheen polished finish (For laser application – non-reflective black, chromium or silicone beaded finish available)
- Manufactured in-house by Elmed in Addison, Illinois, U.S.A.
- Most features can be interchanged to customize the forceps to suit the surgeon's needs.
- Special length and tip size on request.



2-Pin Connector Styles




Angled Tips



Straight Tips



Standard 2-Pin
Connector



2-Pin Connector W/
Irrigation



TO ORDER CALL: (224) 353-6446



BIPOLAR BAYONET FORCEPS—FOOT-SWITCHING, 2-PIN CONNECTION



MADE IN THE USA



CUSHING—STANDARD BAYONET	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST #	*TIP SIZE	LIST #	*TIP SIZE	LIST #	*TIP SIZE	LIST #	*TIP SIZE
BAYONET, 5 3/4" (14.5 CM)	5910-05		5911-05		5912-05		5913-05	
BAYONET, 6 3/4" (17 CM)	5910-06		5911-06		5912-06		5913-06	
BAYONET, 7 3/4" (19.5 CM)	5910-07		5911-07		5912-07		5913-07	
BAYONET, 8 3/4" (22 CM)	5910-08		5911-08		5912-08		5913-08	
BAYONET, 10 3/4" (27.3 CM)	5910-10		5911-10		5912-10		5913-10	
BAYONET, 12 3/4" (32.3 CM)	5910-12		5911-12		5912-12		5913-12	



SCOVILLE—GREENWOOD	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST #	*TIP SIZE	LIST #	*TIP SIZE
BAYONET, 6 3/4" (17 CM)	5940-06		5942-06	
BAYONET, 7 3/4" (19.5 CM)	5940-07		5942-07	
BAYONET, 8 3/4" (22 CM)	5940-08		5942-08	



YASARGIL	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST #	*TIP SIZE	LIST #	*TIP SIZE
BAYONET, 7" (17.8 CM)	5950-06		5952-06	
BAYONET, 8" (20.3 CM)	5950-07		5952-07	
BAYONET, 9" (22.8 CM)	5950-08		5952-08	
BAYONET, 10 3/4" (27.3 CM)	5950-10		5952-10	




GERALD	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST #	*TIP SIZE	LIST #	*TIP SIZE
BAYONET, 6 3/4" (17 CM)	5960-06	C	5962-06	C
BAYONET, 7 3/4" (19.5 CM)	5960-07	C	5962-07	C

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS


INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





MALIS	INSULATED STRAIGHT TIP	INSULATED ANG. TIP, CVD. UP		NON-INSULATED STRAIGHT TIP		NON-INSULATED, ANG. TIP, CVD UP		INSULATED, ANG. TIP, CVD DOWN		NON-INSULATED, ANG. TIP, CVD DOWN	
	LIST NO.	LIST NO.	TIP SIZE	LIST NO.	TIP SIZE	LIST NO.	TIP SIZE	LIST NO.	TIP SIZE	LIST NO.	TIP SIZE
BAYONET, 6 3/4" (17 CM)	5970-07	5971-07U	D	5972-07	D	5973-07U	D	5971-07D	D	5973-07D	D
BAYONET, 7 3/4" (19.5 CM)	5970-08	5971-08U	D	5972-08	D	5973-08U	D	5971-08D	D	5973-08D	D
BAYONET, 8 3/4" (22 CM)	5970-09	5971-09U	D	5972-09	D	5973-09U	D	5971-09D	D	5973-09D	D
BAYONET, 10 3/4" (27.3 CM)	5970-10	5971-10U	D	5972-10	D	5973-10U	D	5971-10D	D	5973-10D	D



RHOTON	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 6 3/4" (17 CM)	5980-06		5981-06		5982-06		5983-06	
BAYONET, 7 3/4" (19.5 CM)	5980-07		5981-07		5982-07		5983-07	
BAYONET, 8 3/4" (22 CM)	5980-08		5981-08		5982-08		5983-08	
BAYONET, 9 3/4" (24.7CM)	5980-09		5981-09		5982-09		5983-09	




HARDY	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 8 3/4", (20CM)	5990-08		5992-08	
BAYONET, 9 3/4", (24.7CM)	5990-09		5992-09	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS

INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





BIPOLAR STRAIGHT FORCEPS—FOOT-SWITCHING, 2-PIN CONNECTION



MADE IN THE USA



CUSHING	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 4¼" (12CM)	5814-04		5815-04		5816-04		5817-04	
STRAIGHT, 5¼" (14.5CM)	5814-05		5815-05		5816-05		5817-05	
STRAIGHT, 6¼" (17CM)	5814-06		5815-06		5816-06		5817-06	
STRAIGHT, 7¼" (19.5CM)	5814-07		5815-07		5816-07		5817-07	
STRAIGHT, 8¼" (22CM)	5814-08		5815-08		5816-08		5817-08	
STRAIGHT, 10¼" (25.3CM)	5814-10		5815-10		5816-10		5817-10	
STRAIGHT, 12¼" (32.3CM)	5814-12		5815-12		5816-12		5817-12	




GERALD	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 6¼" (17CM)	5860-06	C, D, E	5862-06	C, D, E
STRAIGHT, 7¼" (19.5CM)	5860-07	C, D, E	5862-07	C, D, E



MALIS	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 6¼" (14CM)	5870-07		5872-07	
STRAIGHT, 6¼" (17CM)	5870-08		5872-08	
STRAIGHT, 7¼" (19.5CM)	5870-09		5872-09	
STRAIGHT, 9¼" (21.5CM)	5870-10		5872-10	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS		
INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





RHOTON	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 7¼" (19.5CM)	5880-07		5881-07		5882-07		5883-07	
STRAIGHT, 8¾" (22CM)	5880-08		5881-08		5882-08		5883-08	
STRAIGHT, 9¾" (24.7CM)	5880-09		5881-09		5882-09		5883-09	




SEMKIN	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 5½" (13.5CM)	5890-05		5892-06	
STRAIGHT, 6¾" (15CM)	5890-06		5892-06	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS

INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





BIPOLAR BAYONET FORCEPS WITH IRRIGATION — FOOT-SWITCHING, 2-PIN CONNECTION



MADE IN THE USA



CUSHING	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 6¾" (17CM)	5531-06		5532-06		5533-06		5534-06	
BAYONET, 7¾" (19.5CM)	5531-07		5532-07		5533-07		5534-07	
BAYONET, 8¾" (22CM)	5531-08		5532-08		5533-08		5534-08	
BAYONET, 10¾" (25.3CM)	5531-09		5532-09		5533-09		5534-09	



SCOVILLE—GREENWOOD	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 6¾" (17CM)	5540-06		5542-06	
BAYONET, 7¾" (19.5CM)	5540-07		5542-07	
BAYONET, 8¾" (22CM)	5540-08		5542-08	




YASARGIL	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 7" (17.8CM)	5550-06		5552-06	
BAYONET, 8" (20.3CM)	5550-07		5552-07	
BAYONET, 9" (22.8CM)	5550-08		5552-08	
BAYONET, 10¾" (25.3CM)	5550-10		5552-10	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS

INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





MALIS	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANG. TIP		NON-INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ ANG. TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 6¾" (17CM)	5570-07		CURVED UP 5571-07U	D	5572-07		CURVED UP 5573-07U	D
			CURVED DN 5571-07D	D			CURVED DN 5573-07D	D
BAYONET, 7¾" (19.5CM)	5570-08		CURVED UP 5571-08U	D	5572-08		CURVED UP 5573-08U	D
			CURVED DN 5571-08D	D			CURVED DN 5573-08D	D
BAYONET, 8¾" (22CM)	5570-09		CURVED UP 5571-09U	D	5572-09		CURVED UP 5573-09U	D
			CURVED DN 5571-09D	D			CURVED DN 5573-09D	D
BAYONET, 10¾" (25.3CM)	5570-10		CURVED UP 5571-10U	D	5572-10		CURVED UP 5573-09U	D
			CURVED DN 5571-10D	D			CURVED DN 5573-09D	D




SPECIALTY MALIS— 1		INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
		LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 8¾" (22CM)	90° BLUNT CURVE UP, .75 MM	5575-08		5576-08	
BAYONET, 8¾" (22CM)	90° BLUNT CURVE DN, .75 MM	5575-09		5576-09	
BAYONET, 10¾" (25.3CM)	90° BLUNT CURVE UP, .75 MM	5575-10		5576-10	
BAYONET, 10¾" (25.3CM)	90° BLUNT CURVE DN, .75 MM	5575-11		5576-11	



SPECIALTY MALIS— 2		INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
		LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 8¾" (22CM)	CIRCULAR LOOP, 5 MM Ø	5577-05		5578-05	
BAYONET, 8¾" (22CM)	CIRCULAR LOOP, 7 MM Ø	5577-07		5578-07	
BAYONET, 10¾" (25.3CM)	CIRCULAR LOOP, 10 MM Ø	5577-10		5578-10	
BAYONET, 10¾" (25.3CM)	DUCK-BILLED, 1 MM Ø	5577-11		5578-11	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS		
INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





RHOTON	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 9" (22.8CM)	5580-08		5582-08	
BAYONET, 9¾" (24.7CM)	5580-09		5582-09	



HARDY	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
BAYONET, 9" (22.8CM)	5590-08		5592-08	

BIPOLAR STRAIGHT FORCEPS WITH IRRIGATION FOOT-SWITCHING, 2-PIN CONNECTION



CUSHING STANDARD	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 7¾" (19.5CM)	5432-07		5434-07	
STRAIGHT, 8¾" (22CM)	5432-08		5434-08	
STRAIGHT, 10¾" (25.3CM)	5432-10		5434-10	
STRAIGHT, 12¾" (32.3CM)	5432-12		5434-12	




MALIS		INSULATED W/ STRAIGHT TIP	INSULATED W/ STRAIGHT TIP	NON-INSULATED W/ STRAIGHT TIP	NON-INSULATED W/ STRAIGHT TIP
		LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 5¾" (14.5CM)	.75MM Tip	5470-05		5472-05	
STRAIGHT, 5¾" (14.5CM)	.25MM tip	5470-06		5472-06	
STRAIGHT, 6¾" (17CM)	.25MM tip	5470-07		5472-07	
STRAIGHT, 6¾" (17CM)	.45MM Tip	5470-08		5472-08	
STRAIGHT 7¾" (19.5CM)	.25MM tip	5470-09		5472-09	

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS

INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





RHOTON	INSULATED W/ STRAIGHT TIP		NON-INSULATED W/ STRAIGHT TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 7¼" (19.5CM)	5480-07	.7MM Tip	5482-07	.7MM Tip
STRAIGHT, 9" (22.8CM)	5480-08	.7MM Tip	5482-08	.7MM Tip
STRAIGHT, 9¾" (24.7CM)	5480-09	.7MM Tip	5482-09	.7MM Tip

SPECIAL MICRO & OPHTHALMIC BIPOLAR FORCEPS



ADSON	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANGLED TIP		NON-INSULATED W/ STRAIGHT TIAP		NON-INSULATED W/ ANGLED TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 4¾" (10CM)	5710-04	C	5711-04	C	5712-04	C	5713-04	C
STARIGHT 5½" (14CM)	5710-05	C	5711-05	C	5712-05	C	5713-05	C
STRAIGHT, 5½" (14CM)	5710-06	1X2 TEETH			5712-06		1X2 TEETH	



JEWELERS	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANGLED TIP		NON-INSULATED W/ STRAIGHT TIAP		NON-INSULATED W/ ANGLED TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 4¼" (10.7CM)	5720-04	E	5721-04	E	5722-04	E	5723-04	E
STRAIGHT, 5½" 14CM)	5720-05	E	5721-05	E	5722-05	E	5723-05	E




IRIS	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANGLED TIP		NON-INSULATED W/ STRAIGHT TIAP		NON-INSULATED W/ ANGLED TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 3¾" (8CM)	5730-04	E	5731-04	E	5732-04	E	5733-04	E
STRAIGHT, 5¾" (14.5CM)	5730-05	E	5731-05	E	5732-05	E	5733-05	E

SAMPLE LIST NUMBER ORDERING INSTRUCTIONS


INDICATE	LIST #	TIP SIZE
	5000-00	+ A, B, C, D, E

*TIP SIZE SPECIFICATIONS - SELECT AS NEEDED





MCPHERSON	INSULATED W/ STRAIGHT TIP		INSULATED W/ ANGLED TIP		NON-INSULATED W/ STRAIGHT TIAP		NON-INSULATED W/ ANGLED TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 4¼" (10.7CM)	5740-04	D, E	5741-04	D, E	5742-04	D, E	5743-04	D, E
STRAIGHT, 5¾" (13.3CM)	5740-05	D, E	5741-05	D, E	5742-05	D, E	5743-05	D, E



COAPTATION	INSULATED W/ ANGLED TIP		NON-INSULATED W/ ANGLED TIP	
	LIST NO.	*TIP SIZE	LIST NO.	*TIP SIZE
STRAIGHT, 4¾" (22.8CM)	5751-04	D, E	5753-04	D,E
STRAIGHT, 5¼" (13.3CM)	5751-05	D,E	5753-05	D, E



ELMED BIPOLAR OUTPUT TESTER

A unique, simple and economical test box to verify operation of the bipolar coagulation system in the surgical suite.

List # 52-5297

BIPOLAR CABLES

LIST #	DESCRIPTION
5990	BIPOLAR CABLE w/ 2-PIN INST. CONNECTOR & DUAL BANANA PLUG GENERATOR CONNECTOR
5990	BIPOLAR CABLE w/ 2-PIN INST. CONNECTOR & SINGLE ELMED/WOLF GENERATOR CONNECTOR
5300	ELMED BC 50 M/M - BIPOLAR RF MICRO/MACRO GENERATOR, 50 WATTS



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE