



LIPOSUCTION



WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS,
IL 60139

WE SUBSCRIBE TO
COST CONTAINMENT
AND PROTECTION
OF THE
ENVIRONMENT



THEREFORE, WE
MANUFACTURE
REUSABLE PRODUCTS
FOR A CLEANER
WORLD

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM

TABLE OF CONTENTS

	PAGE
CANNULA HUBS/CONNECTIONS	3
INFILTRATION CANNULAS	5
GENERAL SUCTION CANNULA	9
HARVESTER CANNULAS	26
INJECTOR CANNULAS	29
SPECIALTY CANNULA	33
SINGLE-USE INSTRUMENTS	36
LIPOSUCTION ACCESSORIES	37




LIPOSUCTION



MADE IN THE USA

CANNULA HUBS/CONNECTIONS



LUER LOCK



LUER LOCK
(STAINLESS STEEL)



ULTRA LUER LOCK



THREADED



60CC TOMMEY
HUB

MACHINE CONNECTIONS

381000-01

CANNULAS AND HARVESTERS AS A COMPLETE INSTRUMENT READY TO USE WITH AN ASPIRATION MACHINE FOR LIPOSUCTION AND LARGE VOLUME FAT HARVESTING PROCEDURES



FIXED HANDLE



LIPOSUCTION



MADE IN THE USA

LUER LOCK HUB/CONNECTOR:

THE UNIVERSAL LUER LOCK HUB ENSURES DIRECT COMPATIBILITY AND SECURE ATTACHMENT WITH STANDARD LUER LOCK SYRINGES.

ULTRA LUER LOCK HUB/CONNECTOR:

AN UPGRADED, HIGH-PERFORMANCE LUER LOCK HUB MODIFICATION DESIGNED FOR ULTIMATE STABILITY AND POWER. ITS ADVANCED ENGINEERING ENSURES A COMPLETELY LEAK-PROOF, WOBBLE-FREE CONNECTION THAT ACTIVELY GUARDS AGAINST FLUID CONTAMINATION.

THREADED HUB/CONNECTOR:

WHILE THREADED CONNECTOR HUBS ARE PREDOMINANTLY ENGINEERED FOR MACHINE-ASSISTED LIPOSUCTION, LIBRA'S LIGHTWEIGHT THREADED CANNULAS OFFER THE ADDED VERSATILITY OF MANUAL SYRINGE COMPATIBILITY. BY UTILIZING THE PROPRIETARY LIBRA THREADED-TO-LUER ADAPTER, CLINICIANS CAN SEAMLESSLY TRANSITION THESE ADVANCED CANNULAS FOR USE WITH STANDARD LUER LOCK SYRINGES.

60CC TOMMEY HUB/CONNECTOR:

THE 60CC TOOMEY HUB IS THE PREMIER CONFIGURATION FOR SURGEONS WHO PRIORITIZE HIGH-VOLUME FAT HARVESTING AND GRAFTING. THIS SPECIALIZED HUB STYLE PROVIDES A DIRECT, STABLE CONNECTION BETWEEN THE CANNULA AND A STANDARD 60CC TOOMEY TIP SYRINGE.

FIXED HANDLE CANNULAS:

THE HIGH-PERFORMANCE ERGONOMIC HANDLE INTEGRATES SEAMLESSLY WITH THE CANNULA OF YOUR CHOICE. ENGINEERED SPECIFICALLY FOR HIGH-VOLUME TISSUE HARVESTING, THE UNIQUE HANDLE GEOMETRY ENSURES A SECURE, POWERFUL GRIP. DESIGNED TO MINIMIZE HAND FATIGUE, IT PROVIDES EXCEPTIONAL COMFORT THROUGHOUT PROLONGED SURGICAL PROCEDURES. ADDITIONALLY, THE INTEGRATED THUMB REST OFFERS PRECISE TACTILE CONTROL OVER INSERTION AND THRUST DYNAMICS.




INFILTRATION CANNULAS



MADE IN THE USA

ELMED'S INFILTRATION CANNULAS ARE DESIGNED FOR SAFE AND EFFICIENT DELIVERY OF TUMESCENT SOLUTION. THE SPIRAL MICRO HOLES ON CANNULA ENSURE AN EVEN DISTRIBUTION OF LOCAL ANESTHETIC WHICH RESULTS IN A LONGER ANESTHETIC EFFECTS THIS LESS PAINFUL AND MORE PLEASANT LIPOSUCTION EXPERIENCE FOR PATIENT.



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA

TUMESCENT INFILTRATOR

THIS CANNULA IS USED SPECIFICALLY FOR THE INFUSION OF TUMESCENT SOLUTION. THE CONFIGURATION OF THE PORTS ALLOWS EVEN DISTRIBUTION OF THE FLUID.



60CC TOMMEY HUB

38101-01	Ø 2.5 MM X 15 CM
38101-02	Ø 2.5 MM X 20 CM
38101-03	Ø 3.0 MM X 20 CM
38101-04	Ø 3.0 MM X 25 CM

THREADED HUB

38102-01	Ø 2.0 MM X 15 CM
38102-02	Ø 3.0 MM X 15 CM
38102-03	Ø 3.0 MM X 20 CM
38102-04	Ø 3.0 MM X 30 CM
38102-05	Ø 4.0 MM X 20 CM
38102-06	Ø 4.0 MM X 30 CM

ULTRA LUER LOCK

38103-01	Ø 1.0 MM X 5 CM
38103-02	Ø 1.0 MM X 7 CM
38103-03	Ø 1.5 MM X 7 CM
38103-04	Ø 1.5 MM X 9 CM

LUER LOCK

38104-01	Ø 1.0 MM X 4 CM
38104-02	Ø 1.0 MM X 5 CM
38104-03	Ø 1.5 MM X 7 CM
38104-04	Ø 1.5 MM X 9 CM
38104-05	Ø 1.5 MM X 10 CM
38104-06	Ø 1.5 MM X 12 CM
38104-07	Ø 2.0 MM X 10 CM
38104-08	Ø 2.0 MM X 12 CM
38104-09	Ø 2.0 MM X 15 CM
38104-10	Ø 2.0 MM X 20 CM
38104-11	Ø 2.5 MM X 15 CM
38104-12	Ø 2.5 MM X 20 CM
38104-13	Ø 3.0 MM X 15 CM
38104-14	Ø 3.0 MM X 20 CM

LUER LOCK (STAINLESS STEEL)

38105-01	Ø 1.0 MM X 4 CM
38105-02	Ø 1.0 MM X 5 CM
38105-03	Ø 1.0 MM X 7 CM
38105-04	Ø 1.5 MM X 7 CM
38105-05	Ø 1.5 MM X 9 CM
38105-06	Ø 1.5 MM X 10 CM
38105-07	Ø 1.5 MM X 12 CM
38105-08	Ø 2.0 MM X 10 CM
38105-09	Ø 2.0 MM X 12 CM
38105-10	Ø 2.0 MM X 15 CM
38105-11	Ø 2.0 MM X 20 CM
38105-12	Ø 2.5 MM X 15 CM
38105-13	Ø 2.5 MM X 20 CM
38105-14	Ø 3.0 MM X 15 CM
38105-15	Ø 3.0 MM X 20 CM

POWER CANNULA

38106-01	Ø 2.0 MM X 15 CM
38106-02	Ø 2.0 MM X 20 CM
38106-03	Ø 2.0 MM X 25 CM
38106-04	Ø 2.0 MM X 30 CM
38106-05	Ø 2.0 MM X 35 CM
38106-06	Ø 3.0 MM X 15 CM
38106-07	Ø 3.0 MM X 20 CM
38106-08	Ø 3.0 MM X 25 CM
38106-09	Ø 3.0 MM X 30 CM
38106-10	Ø 3.0 MM X 35 CM
38106-11	Ø 4.0 MM X 15 CM
38106-12	Ø 4.0 MM X 20 CM
38106-13	Ø 4.0 MM X 25 CM
38106-14	Ø 4.0 MM X 30 CM
38106-15	Ø 4.0 MM X 35 CM
38106-16	Ø 5.0 MM X 15 CM
38106-17	Ø 5.0 MM X 20 CM
38106-18	Ø 5.0 MM X 25 CM
38106-19	Ø 5.0 MM X 30 CM
38106-20	Ø 5.0 MM X 35 CM
38106-21	Ø 6.0 MM X 15 CM
38106-22	Ø 6.0 MM X 20 CM
38106-23	Ø 6.0 MM X 25 CM
38106-24	Ø 6.0 MM X 30 CM
38106-25	Ø 6.0 MM X 35 CM




INFILTRATION CANNULAS



MADE IN THE USA

TUMESCENT BAYONET INFILTRATOR
THE BAYONET SHAPED HANDLES
ALLOW AN EASIER MANIPULATION



60CC TOMMEY HUB

38101-11	Ø 2.0 MM X 15 CM
38101-12	Ø 2.5 MM X 20 CM
38101-13	Ø 3.0 MM X 20 CM
38101-14	Ø 3.0 MM X 25 CM

THREADED HUB

38102-11	Ø 2.0MM X 15 CM
38102-12	Ø 3.0MM X 15 CM
38102-13	Ø 3.0MM X 20 CM
38102-14	Ø 3.0MM X 30 CM
38102-15	Ø 4.0MM X 20 CM
38102-16	Ø 4.0MM X 30 CM

ULTRA LUER LOCK

38103-11	Ø 1.0 MM X 5 CM
38103-12	Ø 1.0 MM X 7 CM
38103-13	Ø 1.5 MM X 7 CM
38103-14	Ø 1.5 MM X 9 CM

LUER LOCK

38104-11	Ø 1.0 MM X 4 CM
38104-12	Ø 1.0 MM X 5 CM
38104-13	Ø 1.0 MM X 7 CM
38104-14	Ø 1.5 MM X 7 CM
38104-15	Ø 1.5 MM X 9 CM
38104-16	Ø 1.5 MM X 10 CM
38104-17	Ø 1.5 MM X 12 CM
38104-18	Ø 2.0 MM X 10 CM
38104-19	Ø 2.0 MM X 12 CM
38104-20	Ø 2.0 MM X 15 CM
38104-21	Ø 2.0 MM X 20 CM
38104-22	Ø 2.5 MM X 15 CM
38104-23	Ø 2.5 MM X 20 CM
38104-24	Ø 3.0 MM X 15 CM
38104-25	Ø 3.0 MM X 20 CM

LUER LOCK (STAINLESS STEEL)

38105-21	Ø 1.0 MM X 4 CM
38105-22	Ø 1.0 MM X 5 CM
38105-23	Ø 1.0 MM X 7 CM
38105-24	Ø 1.5 MM X 7 CM
38105-25	Ø 1.5 MM X 9 CM
38105-26	Ø 1.5 MM X 10 CM
38105-27	Ø 1.5 MM X 12 CM
38105-28	Ø 2.0 MM X 10 CM
38105-29	Ø 2.0 MM X 12 CM
38105-30	Ø 2.0 MM X 15 CM
38105-31	Ø 2.0 MM X 20 CM
38105-32	Ø 2.5 MM X 15 CM
38105-33	Ø 2.5 MM X 20 CM
38105-34	Ø 3.0 MM X 15 CM
38105-35	Ø 3.0 MM X 20 CM

POWER CANNULA

38106-31	Ø 2 MM X 15 CM
38106-32	Ø 2 MM X 20 CM
38106-33	Ø 2 MM X 25 CM
38106-34	Ø 2 MM X 30 CM
38106-35	Ø 2 MM X 35 CM
38106-36	Ø 3 MM X 15 CM
38106-37	Ø 3 MM X 20 CM
38106-38	Ø 3 MM X 25 CM
38106-39	Ø 3 MM X 30 CM
38106-40	Ø 3 MM X 35 CM
38106-41	Ø 4 MM X 15 CM
38106-42	Ø 4 MM X 20 CM
38106-43	Ø 4 MM X 25 CM
38106-44	Ø 4 MM X 30 CM
38106-45	Ø 4 MM X 35 CM
38106-46	Ø 5 MM X 15 CM
38106-47	Ø 5 MM X 20 CM
38106-48	Ø 5 MM X 25 CM
38106-49	Ø 5 MM X 30 CM
38106-50	Ø 5 MM X 35 CM
38106-51	Ø 6 MM X 15 CM
38106-52	Ø 6 MM X 20 CM
38106-53	Ø 6 MM X 25 CM
38106-54	Ø 6 MM X 30 CM
38106-55	Ø 6 MM X 35 CM



INFILTRATION CANNULAS



MADE IN THE USA

TUMESCENT INFILTRATION ANGLED

THIS CANNULA IS USED SPECIFICALLY FOR THE INFUSION OF TUMESCENT SOLUTION. THE CONFIGURATION OF THE PORTS ALLOWS EVEN DISTRIBUTION OF THE FLUID.



60CC TOMMEY HUB

38101-21	Ø 2.0 MM X 15 CM
38101-22	Ø 2.5 MM X 20 CM
38101-23	Ø 3.0 MM X 20 CM
38101-24	Ø 3.0 MM X 25 CM

THREADED

38102-21	Ø 2.0MM X 15 CM
38102-22	Ø 3.0MM X 15 CM
38102-23	Ø 3.0MM X 20 CM
38102-24	Ø 3.0MM X 30 CM
38102-25	Ø 4.0MM X 20 CM
38102-26	Ø 4.0MM X 30 CM

ULTRA LUER LOCK

38103-21	Ø 1.0 MM X 5 CM
38103-22	Ø 1.0 MM X 7 CM
38103-23	Ø 1.5 MM X 7 CM
38103-24	Ø 1.5 MM X 9 CM

LUER LOCK

38104-31	Ø 1.0 MM X 4 CM
38104-32	Ø 1.0 MM X 5 CM
38104-33	Ø 1.0 MM X 7 CM
38104-34	Ø 1.5 MM X 7 CM
38104-35	Ø 1.5 MM X 9 CM
38104-36	Ø 1.5 MM X 10 CM
38104-37	Ø 1.5 MM X 12 CM
38104-38	Ø 2.0 MM X 10 CM
38104-39	Ø 2.0 MM X 12 CM
38104-40	Ø 2.0 MM X 15 CM
38104-41	Ø 2.0 MM X 20 CM
38104-42	Ø 2.5 MM X 15 CM
38104-43	Ø 2.5 MM X 20 CM
38104-44	Ø 3.0 MM X 15 CM
38104-45	Ø 3.0 MM X 20 CM

LUER LOCK (STAINLESS STEEL)

38105-41	Ø 1.0 MM X 4 CM
38105-42	Ø 1.0 MM X 5 CM
38105-43	Ø 1.0 MM X 7 CM
38105-44	Ø 1.5 MM X 7 CM
38105-45	Ø 1.5 MM X 9 CM
38105-46	Ø 1.5 MM X 10 CM
38105-47	Ø 1.5 MM X 12 CM
38105-48	Ø 2.0 MM X 10 CM
38105-49	Ø 2.0 MM X 12 CM
38105-50	Ø 2.0 MM X 15 CM
38105-51	Ø 2.0 MM X 20 CM
38105-52	Ø 2.5 MM X 15 CM
38105-53	Ø 2.5 MM X 20 CM
38105-54	Ø 3.0 MM X 15 CM
38105-55	Ø 3.0 MM X 20 CM

POWER CANNULA

38106-61	Ø 2 MM X 15 CM
38106-62	Ø 2 MM X 20 CM
38106-63	Ø 2 MM X 25 CM
38106-64	Ø 2 MM X 30 CM
38106-65	Ø 2 MM X 35 CM
38106-66	Ø 3 MM X 15 CM
38106-67	Ø 3 MM X 20 CM
38106-68	Ø 3 MM X 25 CM
38106-69	Ø 3 MM X 30 CM
38106-70	Ø 3 MM X 35 CM
38106-71	Ø 4 MM X 15 CM
38106-72	Ø 4 MM X 20 CM
38106-73	Ø 4 MM X 25 CM
38106-74	Ø 4 MM X 30 CM
38106-75	Ø 4 MM X 35 CM
38106-76	Ø 5 MM X 15 CM
38106-77	Ø 5 MM X 20 CM
38106-78	Ø 5 MM X 25 CM
38106-79	Ø 5 MM X 30 CM
38106-80	Ø 5 MM X 35 CM
38106-81	Ø 6 MM X 15 CM
38106-82	Ø 6 MM X 20 CM
38106-83	Ø 6 MM X 25 CM
38106-84	Ø 6 MM X 30 CM
38106-85	Ø 6 MM X 35 CM



INFILTRATION CANNULAS



MADE IN THE USA

FACELIFT INFILTRATION CANNULA

THIS CANNULA IS SPECIFICALLY FOR THE INFUSION OF TUMESCENT SOLUTION. THE CONFIGURATION OF THE PORTS ALLOWS EVEN DISTRIBUTION OF THE FLUID.



60CC TOMMEY HUB

38101-31	Ø 2.0 MM X 15 CM
38101-32	Ø 2.5 MM X 20 CM
38101-33	Ø 3.0 MM X 20 CM
38101-34	Ø 3.0 MM X 25 CM

THREADED

38102-31	Ø 2.0MM X 15 CM
38102-32	Ø 3.0MM X 15 CM
38102-33	Ø 3.0MM X 20 CM
38102-34	Ø 3.0MM X 30 CM
38102-35	Ø 4.0MM X 20 CM
38102-36	Ø 4.0MM X 30 CM

ULTRA LUER LOCK

38103-31	Ø 1.0 MM X 5 CM
38103-32	Ø 1.0 MM X 7 CM
38103-33	Ø 1.5 MM X 7 CM
38103-34	Ø 1.5 MM X 9 CM

LUER LOCK

38104-51	Ø 1.0 MM X 4 CM
38104-52	Ø 1.0 MM X 5 CM
38104-53	Ø 1.0 MM X 7 CM
38104-54	Ø 1.5 MM X 7 CM
38104-55	Ø 1.5 MM X 9 CM
38104-56	Ø 1.5 MM X 10 CM
38104-57	Ø 1.5 MM X 12 CM
38104-58	Ø 2.0 MM X 10 CM
38104-59	Ø 2.0 MM X 12 CM
38104-50	Ø 2.0 MM X 15 CM
38104-51	Ø 2.0 MM X 20 CM
38104-52	Ø 2.5 MM X 15 CM
38104-53	Ø 2.5 MM X 20 CM
38104-54	Ø 3.0 MM X 15 CM
38104-55	Ø 3.0 MM X 20 CM

LUER LOCK (STAINLESS STEEL)

38105-61	Ø 1.0 MM X 4 CM
38105-62	Ø 1.0 MM X 5 CM
38105-63	Ø 1.0 MM X 7 CM
38105-64	Ø 1.5 MM X 7 CM
38105-65	Ø 1.5 MM X 9 CM
38105-66	Ø 1.5 MM X 10 CM
38105-67	Ø 1.5 MM X 12 CM
38105-68	Ø 2.0 MM X 10 CM
38105-69	Ø 2.0 MM X 12 CM
38105-70	Ø 2.0 MM X 15 CM
38105-71	Ø 2.0 MM X 20 CM
38105-72	Ø 2.5 MM X 15 CM
38105-73	Ø 2.5 MM X 20 CM
38105-74	Ø 3.0 MM X 15 CM
38105-75	Ø 3.0 MM X 20 CM

POWER CANNULA

38106-91	Ø 2 MM X 15 CM
38106-92	Ø 2 MM X 20 CM
38106-93	Ø 2 MM X 25 CM
38106-94	Ø 2 MM X 30 CM
38106-95	Ø 2 MM X 35 CM
38106-96	Ø 3 MM X 15 CM
38106-97	Ø 3 MM X 20 CM
38106-98	Ø 3 MM X 25 CM
38106-99	Ø 3 MM X 30 CM
38106-100	Ø 3 MM X 35 CM
38106-101	Ø 4 MM X 15 CM
38106-102	Ø 4 MM X 20 CM
38106-103	Ø 4 MM X 25 CM
38106-104	Ø 4 MM X 30 CM
38106-105	Ø 4 MM X 35 CM
38106-106	Ø 5 MM X 15 CM
38106-107	Ø 5 MM X 20 CM
38106-108	Ø 5 MM X 25 CM
38106-109	Ø 5 MM X 30 CM
38106-110	Ø 5 MM X 35 CM
38106-111	Ø 6 MM X 15 CM
38106-112	Ø 6 MM X 20 CM
38106-113	Ø 6 MM X 25 CM
38106-114	Ø 6 MM X 30 CM
38106-115	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

ELMED'S GENERAL SUCTION CANNULA ARE IDEAL TOOL FOR LARGE VOLUME HARVESTING WITH EASE PERFECTION. THESE CANNULAS ARE DESIGNED TO FACILITATE YOU TO GET PERFECTLY DESIGNED GOALS IN MINIMUM TIME.



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA

STEVENS STYLE CANNULA



60CC TOMMEY HUB

38201-01	Ø 2.0 MM X 15 CM
38201-02	Ø 2.5 MM X 20 CM
38201-03	Ø 3.0 MM X 20 CM
38201-04	Ø 3.0 MM X 25 CM

THREADED

38202-01	Ø 3.0 MM X 20 CM
38202-02	Ø 3.0 MM X 30 CM
38202-03	Ø 4.0 MM X 30 CM
38202-04	Ø 5.0 MM X 30 CM
38202-05	Ø 6.0 MM X 30 CM

ULTRA LUER LOCK

38203-01	Ø 1.0 MM X 5 CM
38203-02	Ø 1.0 MM X 7 CM
38203-03	Ø 1.5 MM X 7 CM
38203-04	Ø 1.5 MM X 9 CM
38203-05	Ø 1.5 MM X 7 CM

LUER LOCK

38204-01	Ø 1.0 MM X 4 CM
38204-02	Ø 1.0 MM X 5 CM
38204-03	Ø 1.0 MM X 7 CM
38204-04	Ø 1.5 MM X 7 CM
38204-05	Ø 1.5 MM X 9 CM
38204-06	Ø 1.5 MM X 10 CM
38204-07	Ø 1.5 MM X 12 CM
38204-08	Ø 2.0 MM X 10 CM
38204-09	Ø 2.0 MM X 12 CM
38204-10	Ø 2.0 MM X 15 CM
38204-11	Ø 2.0 MM X 20 CM
38204-12	Ø 2.5 MM X 15 CM
38204-13	Ø 2.5 MM X 20 CM
38204-14	Ø 3.0 MM X 15 CM
38204-15	Ø 3.0 MM X 20 CM

LUER LOCK (STAINLESS STEEL)

38205-01	Ø 1.0 MM X 4 CM
38205-02	Ø 1.0 MM X 5 CM
38205-03	Ø 1.0 MM X 7 CM
38205-04	Ø 1.5 MM X 7 CM
38205-05	Ø 1.5 MM X 9 CM
38205-06	Ø 1.5 MM X 10 CM
38205-07	Ø 1.5 MM X 12 CM
38205-08	Ø 2.0 MM X 10 CM
38205-09	Ø 2.0 MM X 12 CM
38205-10	Ø 2.0 MM X 15 CM
38205-11	Ø 2.0 MM X 20 CM
38205-12	Ø 2.5 MM X 15 CM
38205-13	Ø 2.5 MM X 20 CM
38205-14	Ø 3.0 MM X 15 CM
38205-15	Ø 3.0 MM X 20 CM

POWER CANNULA

38206-01	Ø 2 MM X 15 CM
38206-02	Ø 2 MM X 20 CM
38206-03	Ø 2 MM X 25 CM
38206-04	Ø 2 MM X 30 CM
38206-05	Ø 2 MM X 35 CM
38206-06	Ø 3 MM X 15 CM
38206-07	Ø 3 MM X 20 CM
38206-08	Ø 3 MM X 25 CM
38206-09	Ø 3 MM X 30 CM
38206-10	Ø 3 MM X 35 CM
38206-11	Ø 4 MM X 15 CM
38206-12	Ø 4 MM X 20 CM
38206-13	Ø 4 MM X 25 CM
38206-14	Ø 4 MM X 30 CM
38206-15	Ø 4 MM X 35 CM
38206-16	Ø 5 MM X 15 CM
38206-17	Ø 5 MM X 20 CM
38206-18	Ø 5 MM X 25 CM
38206-19	Ø 5 MM X 30 CM
38206-20	Ø 5 MM X 35 CM
38206-21	Ø 6 MM X 15 CM
38206-22	Ø 6 MM X 20 CM
38206-23	Ø 6 MM X 25 CM
38206-24	Ø 6 MM X 30 CM
38206-25	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

COBRA BIBEVEL STYLE

MORE AGGRESSIVE THAN THE STANDARD, PYRAMID OR MERCEDES, THIS CANNULA HAS A TAPERED TIP WITH TWP PORTS AND WORKS EFFECTIVELY WITHIN DENSE FIBROUS TISSUE



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA



60CC TOMMEY HUB

38201-11	Ø 2.5 MM X 15 CM
38201-12	Ø 2.5 MM X 20 CM
38201-13	Ø 3.0 MM X 15 CM
38201-14	Ø 3.0 MM X 20 CM
38201-15	Ø 3.0 MM X 25 CM
38201-16	Ø 4.0 MM X 25 CM
38201-17	Ø 4.0 MM X 30 CM
38201-18	Ø 4.0 MM X 35 CM
38201-19	Ø 4.5 MM X 30 CM
38201-20	Ø 4.5 MM X 35 CM

THREADED HUB

38202-11	Ø 3.0 MM X 15 CM
38202-12	Ø 3.0 MM X 25 CM
38202-13	Ø 4.0 MM X 15 CM
38202-14	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-11	Ø 2.0 MM X 15 CM
38203-12	Ø 2.5 MM X 15 CM
38203-13	Ø 2.5 MM X 20 CM
38203-14	Ø 3.0 MM X 15 CM
38203-15	Ø 3.0 MM X 20 CM
38203-16	Ø 3.0 MM X 25 CM

LUER LOCK

38204-21	Ø 2.0 MM X 15 CM
38204-22	Ø 2.5 MM X 15 CM
38204-23	Ø 2.5 MM X 20 CM
38204-24	Ø 3.0 MM X 15 CM
38204-25	Ø 3.0 MM X 20 CM
38204-26	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-21	Ø 2.0 MM X 15 CM
38205-22	Ø 2.5 MM X 15 CM
38205-23	Ø 2.5 MM X 20 CM
38205-24	Ø 3.0 MM X 15 CM
38205-25	Ø 3.0 MM X 20 CM
38205-26	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-31	Ø 2MM X 15 CM
38206-32	Ø 2MM X 20 CM
38206-33	Ø 2MM X 25 CM
38206-34	Ø 2MM X 30 CM
38206-35	Ø 2MM X 35 CM
38206-36	Ø 3MM X 15 CM
38206-37	Ø 3MM X 20 CM
38206-38	Ø 3MM X 25 CM
38206-39	Ø 3MM X 30 CM
38206-40	Ø 3MM X 35 CM
38206-41	Ø 4MM X 15 CM
38206-42	Ø 4MM X 20 CM
38206-43	Ø 4MM X 25 CM
38206-44	Ø 4MM X 30 CM
38206-45	Ø 4MM X 35 CM
38206-46	Ø 5MM X 15 CM
38206-47	Ø 5MM X 20 CM
38206-48	Ø 5MM X 25 CM
38206-49	Ø 5MM X 30 CM
38206-50	Ø 5MM X 35 CM
38206-51	Ø 6MM X 15 CM
38206-52	Ø 6MM X 20 CM
38206-53	Ø 6MM X 25 CM
38206-54	Ø 6MM X 30 CM
38206-55	Ø 6MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

COBRA STYLE ROUND TIP

LESS AGGRESSIVE THAN THE COBRA BIBEVEL THIS CELL FRIENDLY CANNULA HAS A ROUNDED TIP WITH TWO PORTS AND WORKS EFFECTIVELY WITHIN DENSE FIBROUS TISSUE DURING LIPOSUCTION PROCEDURES



60CC TOMMEY HUB

38201-21	Ø 2.5 MM X 15 CM
38201-22	Ø 2.5 MM X 20 CM
38201-23	Ø 3.0 MM X 15 CM
38201-24	Ø 3.0 MM X 20 CM
38201-25	Ø 3.0 MM X 25 CM
38201-26	Ø 4.0 MM X 25 CM
38201-27	Ø 4.0 MM X 30 CM
38201-28	Ø 4.0 MM X 35 CM
38201-29	Ø 4.5 MM X 30 CM
38201-30	Ø 4.5 MM X 35 CM

THREADED HUB

38202-21	Ø 3.0 MM X 15 CM
38202-22	Ø 3.0 MM X 25 CM
38202-23	Ø 4.0 MM X 15 CM
38202-24	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-21	Ø 2.0 MM X 15 CM
3820322	Ø 2.5 MM X 15 CM
38203-23	Ø 2.5 MM X 20 CM
38203-24	Ø 3.0 MM X 15 CM
38203-25	Ø 3.0 MM X 20 CM
38203-26	Ø 3.0 MM X 25 CM

LUER LOCK

38204-31	Ø 2.0 MM X 15 CM
38204-32	Ø 2.5 MM X 15 CM
38204-33	Ø 2.5 MM X 20 CM
38204-34	Ø 3.0 MM X 15 CM
38204-35	Ø 3.0 MM X 20 CM
38204-36	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-31	Ø 2.0 MM X 15 CM
38205-32	Ø 2.5 MM X 15 CM
38205-33	Ø 2.5 MM X 20 CM
38205-34	Ø 3.0 MM X 15 CM
38205-.35	Ø 3.0 MM X 20 CM
38205-36	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-61	Ø 2 MM X 15 CM
38206-62	Ø 2 MM X 20 CM
38206-63	Ø 2 MM X 25 CM
38206-64	Ø 2 MM X 30 CM
38206-65	Ø 2 MM X 35 CM
38206-66	Ø 3 MM X 15 CM
38206-67	Ø 3 MM X 20 CM
38206-68	Ø 3 MM X 25 CM
38206-69	Ø 3 MM X 30 CM
38206-70	Ø 3 MM X 35 CM
38206-71	Ø 4 MM X 15 CM
38206-72	Ø 4 MM X 20 CM
38206-73	Ø 4 MM X 25 CM
38206-74	Ø 4 MM X 30 CM
38206-75	Ø 4 MM X 35 CM
38206-76	Ø 5 MM X 15 CM
38206-77	Ø 5 MM X 20 CM
38206-78	Ø 5 MM X 25 CM
38206-79	Ø 5 MM X 30 CM
38206-80	Ø 5 MM X 35 CM
38206-81	Ø 6 MM X 15 CM
38206-82	Ø 6 MM X 20 CM
38206-83	Ø 6 MM X 25 CM
38206-84	Ø 6 MM X 30 CM
38206-85	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

STANDARD, SINGLE HOLE

THIS IS THE MOST COMMON CELL FRIENDLY CANNULA DESIGN FOR SOFT TISSUE REMOVAL. THE SINGLE PORT IS LOCATED 180° FROM THE THUMB REST ON THE HUB.



60CC TOMMEY HUB

38201-31	Ø 2.5 MM X 15 CM
38201-32	Ø 2.5 MM X 20 CM
38201-33	Ø 3.0 MM X 15 CM
38201-34	Ø 3.0 MM X 20 CM
38201-35	Ø 3.0 MM X 25 CM
38201-36	Ø 4.0 MM X 25 CM
38201-37	Ø 4.0 MM X 30 CM
38201-38	Ø 4.0 MM X 35 CM
38201-39	Ø 4.5 MM X 30 CM
38201-40	Ø 4.5 MM X 35 CM

THREADED HUB

38202-31	Ø 3.0 MM X 15 CM
38202-32	Ø 3.0 MM X 25 CM
38202-33	Ø 4.0 MM X 15 CM
38202-34	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-31	Ø 1.5 MM X 12 CM
38203-32	Ø 2.0 MM X 15 CM
38203-33	Ø 2.0 MM X 20 CM
38203-34	Ø 2.5 MM X 15 CM
38203-35	Ø 2.5 MM X 20 CM
38203-36	Ø 3.0 MM X 15 CM
38203-37	Ø 3.0 MM X 20 CM
38203-38	Ø 3.0 MM X 25 CM

LUER LOCK

38204-41	Ø 1.5 MM X 12 CM
38204-42	Ø 2.0 MM X 15 CM
38204-43	Ø 2.0 MM X 20 CM
38204-44	Ø 2.5 MM X 15 CM
38204-45	Ø 2.5 MM X 20 CM
38204-46	Ø 3.0 MM X 15 CM
38204-47	Ø 3.0 MM X 20 CM
38204-48	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-41	Ø 1.5 MM X 12 CM
38205-42	Ø 2.0 MM X 15 CM
38205-43	Ø 2.0 MM X 20 CM
38205-44	Ø 2.5 MM X 15 CM
38205-45	Ø 2.5 MM X 20 CM
38205-46	Ø 3.0 MM X 15 CM
38205-47	Ø 3.0 MM X 20 CM
38205-48	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-91	Ø 2 MM X 15 CM
38206-92	Ø 2 MM X 20 CM
38206-93	Ø 2 MM X 25 CM
38206-94	Ø 2 MM X 30 CM
38206-95	Ø 2 MM X 35 CM
38206-96	Ø 3 MM X 15 CM
38206-97	Ø 3 MM X 20 CM
38206-98	Ø 3 MM X 25 CM
38206-99	Ø 3 MM X 30 CM
38206-100	Ø 3 MM X 35 CM
38206-101	Ø 4 MM X 15 CM
38206-102	Ø 4 MM X 20 CM
38206-103	Ø 4 MM X 25 CM
38206-104	Ø 4 MM X 30 CM
38206-105	Ø 4 MM X 35 CM
38206-106	Ø 5 MM X 15 CM
38206-107	Ø 5 MM X 20 CM
38206-108	Ø 5 MM X 25 CM
38206-109	Ø 5 MM X 30 CM
38206-110	Ø 5 MM X 35 CM
38206-111	Ø 6 MM X 15 CM
38206-112	Ø 6 MM X 20 CM
38206-113	Ø 6 MM X 25 CM
38206-114	Ø 6 MM X 30 CM
38206-115	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

PYRAMID HOLES CANNULA

THIS COMMONLY USED DESIGN HAS THREE PORTS LOCATED 180° FROM THE THUMB PRINT. IT'S SLIGHTLY MORE AGGRESSIVE CANNULA THAN THE STANDARD DURING LIPOSUCTION PROCEDURES



60CC TOMMEY HUB

38201-41	Ø 2.5 MM X 15 CM
38201-42	Ø 2.5 MM X 20 CM
38201-43	Ø 3.0 MM X 15 CM
38201-44	Ø 3.0 MM X 20 CM
38201-45	Ø 3.0 MM X 25 CM
38201-46	Ø 4.0 MM X 25 CM
38201-47	Ø 4.0 MM X 30 CM
38201-48	Ø 4.0 MM X 35 CM
38201-49	Ø 4.5 MM X 30 CM
38201-50	Ø 4.5 MM X 35 CM

THREADED HUB

38202-41	Ø 3.0 MM X 15 CM
38202-42	Ø 3.0 MM X 25 CM
38202-43	Ø 4.0 MM X 15 CM
38202-44	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-41	Ø 1.5 MM X 12 CM
38203-42	Ø 2.0 MM X 15 CM
38203-43	Ø 2.5 MM X 15 CM
38203-44	Ø 2.5 MM X 20 CM
38203-45	Ø 3.0 MM X 15 CM
38203-46	Ø 3.0 MM X 20 CM
38203-47	Ø 3.0 MM X 25 CM

LUER LOCK

38204-51	Ø 1.5 MM X 12 CM
38204-52	Ø 2.0 MM X 15 CM
38204-53	Ø 2.5 MM X 15 CM
38204-54	Ø 2.5 MM X 20 CM
38204-55	Ø 3.0 MM X 15 CM
38204-56	Ø 3.0 MM X 20 CM
38204-57	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-51	Ø 1.5 MM X 12 CM
38205-52	Ø 2.0 MM X 15 CM
38205-53	Ø 2.0 MM X 20 CM
38205-54	Ø 2.5 MM X 15 CM
38205-55	Ø 2.5 MM X 20 CM
38205-56	Ø 3.0 MM X 15 CM
38205-57	Ø 3.0 MM X 20 CM
38205-58	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-121	Ø 2 MM X 15 CM
38206-122	Ø 2 MM X 20 CM
38206-123	Ø 2 MM X 25 CM
38206-124	Ø 2 MM X 30 CM
38206-125	Ø 2 MM X 35 CM
38206-126	Ø 3 MM X 15 CM
38206-127	Ø 3 MM X 20 CM
38206-128	Ø 3 MM X 25 CM
38206-129	Ø 3 MM X 30 CM
38206-130	Ø 3 MM X 35 CM
38206-131	Ø 4 MM X 15 CM
38206-132	Ø 4 MM X 20 CM
38206-133	Ø 4 MM X 25 CM
38206-134	Ø 4 MM X 30 CM
38206-135	Ø 4 MM X 35 CM
38206-136	Ø 5 MM X 15 CM
38206-137	Ø 5 MM X 20 CM
38206-138	Ø 5 MM X 25 CM
38206-139	Ø 5 MM X 30 CM
38206-140	Ø 5 MM X 35 CM
38206-141	Ø 6 MM X 15 CM
38206-142	Ø 6 MM X 20 CM
38206-143	Ø 6 MM X 25 CM
38206-144	Ø 6 MM X 30 CM
38206-145	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

MERCEDES

THIS CIRCUMFERENTIAL THREE PORT CONFIGURATION CELL FRIENDLY LIPOSUCTION CANNULA IS EFFECTIVE FOR MULTI-PLANE SOFT TISSUE HARVESTING.



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA



60CC TOMMEY HUB

38201-51	Ø 2.5 MM X 15 CM
38201-52	Ø 2.5 MM X 20 CM
38201-53	Ø 3.0 MM X 15 CM
38201-54	Ø 3.0 MM X 20 CM
38201-55	Ø 3.0 MM X 25 CM
38201-56	Ø 4.0 MM X 25 CM
38201-57	Ø 4.0 MM X 30 CM
38201-58	Ø 4.0 MM X 35 CM
38201-59	Ø 4.5 MM X 30 CM
38201-60	Ø 4.5 MM X 35 CM

THREADED HUB

38202-51	Ø 3.0 MM X 15 CM
38202-52	Ø 3.0 MM X 25 CM
38202-53	Ø 4.0 MM X 15 CM
38202-54	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-51	Ø 1.5 MM X 12 CM
38203-52	Ø 2.0 MM X 15 CM
38203-53	Ø 2.5 MM X 15 CM
38203-54	Ø 2.5 MM X 20 CM
38203-55	Ø 3.0 MM X 15 CM
38203-56	Ø 3.0 MM X 20 CM
38203-57	Ø 3.0 MM X 25 CM

LUER LOCK

38204-61	Ø 1.5 MM X 12 CM
38204-62	Ø 2.0 MM X 15 CM
38204-63	Ø 2.5 MM X 15 CM
38204-64	Ø 2.5 MM X 20 CM
38204-65	Ø 3.0 MM X 15 CM
38204-66	Ø 3.0 MM X 20 CM
38204-67	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-61	Ø 1.5 MM X 12 CM
38205-62	Ø 2.0 MM X 15 CM
38205-63	Ø 2.0 MM X 20 CM
38205-64	Ø 2.5 MM X 15 CM
38205-65	Ø 2.5 MM X 20 CM
38205-66	Ø 3.0 MM X 15 CM
38205-67	Ø 3.0 MM X 20 CM
38205-68	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-151	Ø 2 MM X 15 CM
38206-152	Ø 2 MM X 20 CM
38206-153	Ø 2 MM X 25 CM
38206-154	Ø 2 MM X 30 CM
38206-155	Ø 2 MM X 35 CM
38206-156	Ø 3 MM X 15 CM
38206-157	Ø 3 MM X 20 CM
38206-158	Ø 3 MM X 25 CM
38206-159	Ø 3 MM X 30 CM
38206-160	Ø 3 MM X 35 CM
38206-161	Ø 4 MM X 15 CM
38206-162	Ø 4 MM X 20 CM
38206-163	Ø 4 MM X 25 CM
38206-164	Ø 4 MM X 30 CM
38206-165	Ø 4 MM X 35 CM
38206-166	Ø 5 MM X 15 CM
38206-167	Ø 5 MM X 20 CM
38206-168	Ø 5 MM X 25 CM
38206-169	Ø 5 MM X 30 CM
38206-170	Ø 5 MM X 35 CM
38206-171	Ø 6 MM X 15 CM
38206-172	Ø 6 MM X 20 CM
38206-173	Ø 6 MM X 25 CM
38206-174	Ø 6 MM X 30 CM
38206-175	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

3 SPIRAL PORT CANNULA (CARRAWAY STYLE)

SPIRAL THREE PORT DESIGN ON THIS CELL FRIENDLY LIPOSUCTION CANNULA IS EFFECTIVE FOR MULTI-PLANE ADIPOSE TISSUE HARVESTING.



LUER LOCK

LUER LOCK
(STAINLESS
STEEL)

ULTRA LUER
LOCK

THREADED

60CC
TOMMEY
HUB

POWER CANNULA



60CC TOMMEY HUB

38201-61	Ø 2.5 MM X 15 CM
38201-62	Ø 2.5 MM X 20 CM
38201-63	Ø 3.0 MM X 15 CM
38201-64	Ø 3.0 MM X 20 CM
38201-65	Ø 3.0 MM X 25 CM
38201-66	Ø 4.0 MM X 25 CM
38201-67	Ø 4.0 MM X 30 CM
38201-68	Ø 4.0 MM X 35 CM
38201-69	Ø 4.5 MM X 30 CM
38201-70	Ø 4.5 MM X 35 CM

THREADED HUB

38202-61	Ø 3.0 MM X 15 CM
38202-62	Ø 3.0 MM X 25 CM
38202-63	Ø 4.0 MM X 15 CM
38202-64	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-61	Ø 1.5 MM X 12 CM
38203-62	Ø 2.0 MM X 15 CM
38203-63	Ø 2.5 MM X 15 CM
38203-64	Ø 2.5 MM X 20 CM
38203-65	Ø 3.0 MM X 15 CM
38203-66	Ø 3.0 MM X 20 CM
38203-67	Ø 3.0 MM X 25 CM

LUER LOCK

38204-71	Ø 1.5 MM X 12 CM
38204-72	Ø 2.0 MM X 15 CM
38204-73	Ø 2.5 MM X 15 CM
38204-74	Ø 2.5 MM X 20 CM
38204-75	Ø 3.0 MM X 15 CM
38204-76	Ø 3.0 MM X 20 CM
38204-77	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-71	Ø 1.5 MM X 12 CM
38205-72	Ø 2.0 MM X 15 CM
38205-73	Ø 2.0 MM X 20 CM
38205-74	Ø 2.5 MM X 15 CM
38205-75	Ø 2.5 MM X 20 CM
38205-76	Ø 3.0 MM X 15 CM
38205-77	Ø 3.0 MM X 20 CM
3820578	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-181	Ø 2 MM X 15 CM
38206-182	Ø 2 MM X 20 CM
38206-183	Ø 2 MM X 25 CM
38206-184	Ø 2 MM X 30 CM
38206-185	Ø 2 MM X 35 CM
38206-186	Ø 3 MM X 15 CM
38206-187	Ø 3 MM X 20 CM
38206-188	Ø 3 MM X 25 CM
38206-189	Ø 3 MM X 30 CM
38206-190	Ø 3 MM X 35 CM
38206-191	Ø 4 MM X 15 CM
38206-192	Ø 4 MM X 20 CM
38206-193	Ø 4 MM X 25 CM
38206-194	Ø 4 MM X 30 CM
38206-195	Ø 4 MM X 35 CM
38206-196	Ø 5 MM X 15 CM
38206-197	Ø 5 MM X 20 CM
38206-198	Ø 5 MM X 25 CM
38206-199	Ø 5 MM X 30 CM
38206-200	Ø 5 MM X 35 CM
38206-201	Ø 6 MM X 15 CM
38206-202	Ø 6 MM X 20 CM
38206-203	Ø 6 MM X 25 CM
38206-204	Ø 6 MM X 30 CM
38206-205	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

SIX SPIRAL PORT (ADJUST BELOW 3 SPIRAL PORTS)

SIX PORT SPIRAL CONFIGURATION ON THIS CELL FRIENDLY LIPOSUCTION CANNULA HARVESTS MAXIMUM TISSUE VOLUME IN MINIMAL TIME.



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA



60CC SMART HUB

38201-91	Ø 2.5 MM X 15 CM
38201-92	Ø 2.5 MM X 20 CM
38201-93	Ø 3.0 MM X 15 CM
38201-94	Ø 3.0 MM X 25 CM
38201-95	Ø 4.0 MM X 25 CM

THREADED HUB

38202-81	Ø 3.0 MM X 15 CM
38202-82	Ø 3.0 MM X 25 CM
38202-83	Ø 4.0 MM X 15 CM
38202-84	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-81	Ø 2.5 MM X 15 CM
38203-82	Ø 3.0 MM X 15 CM
38203-83	Ø 3.0 MM X 25 CM
38203-84	Ø 4.0 MM X 15 CM
38203-85	Ø 4.0 MM X 25 CM

LUER LOCK

38204-91	Ø 2.5 MM X 15 CM
38204-92	Ø 3.0 MM X 15 CM
38204-93	Ø 3.0 MM X 25 CM
38204-94	Ø 4.0 MM X 15 CM
38204-95	Ø 4.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-231	Ø 2.5 MM X 15 CM
38205-232	Ø 3.0 MM X 15 CM
38205-233	Ø 3.0 MM X 25 CM
38205-234	Ø 4.0 MM X 15 CM
38205-235	Ø 4.0 MM X 25 CM

POWER CANNULA

38206-901	2MM OD X 15 CM
38206-902	2MM OD X 20 CM
38206-903	2MM OD X 25 CM
38206-904	2MM OD X 30 CM
38206-905	2MM OD X 35 CM
38206-906	3MM OD X 15 CM
38206-907	3MM OD X 20 CM
38206-908	3MM OD X 25 CM
38206-909	3MM OD X 30 CM
38206-910	3MM OD X 35 CM
38206-911	4MM OD X 15 CM
38206-912	4MM OD X 20 CM
38206-913	4MM OD X 25 CM
38206-914	4MM OD X 30 CM
38206-915	4MM OD X 35 CM
38206-916	5MM OD X 15 CM
38206-917	5MM OD X 20 CM
38206-918	5MM OD X 25 CM
38206-919	5MM OD X 30 CM
38206-920	5MM OD X 35 CM
38206-921	6MM OD X 15 CM
38206-922	6MM OD X 20 CM
38206-923	6MM OD X 25 CM
38206-924	6MM OD X 30 CM
38206-925	6MM OD X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

THREE PORT CANNULA

THIS BODY SCULPTING INSTRUMENT IS EXCELLENT FOR HARVESTING SMALL ALIQUOTS OF TISSUE DURING LIPOSUCTION PROCEDURES.



60CC TOMMEY HUB

38201-71	Ø 2.5 MM X 15 CM
38201-72	Ø 2.5 MM X 20 CM
38201-73	Ø 3.0 MM X 15 CM
38201-74	Ø 3.0 MM X 20 CM
38201-75	Ø 3.0 MM X 25 CM
38201-76	Ø 4.0 MM X 25 CM
38201-77	Ø 4.0 MM X 30 CM
38201-78	Ø 4.0 MM X 35 CM
38201-79	Ø 4.5 MM X 30 CM
38201-80	Ø 4.5 MM X 35 CM

THREADED HUB

38202-71	Ø 3.0 MM X 15 CM
38202-72	Ø 3.0 MM X 25 CM
38202-73	Ø 4.0 MM X 15 CM
38202-74	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38203-71	Ø 1.5 MM X 12 CM
38203-72	Ø 2.0 MM X 15 CM
38203-73	Ø 2.5 MM X 15 CM
38203-74	Ø 2.5 MM X 20 CM
38203-75	Ø 3.0 MM X 15 CM
38203-76	Ø 3.0 MM X 20 CM
38203-77	Ø 3.0 MM X 25 CM

LUER LOCK

38204-81	Ø 1.5 MM X 12 CM
38204-82	Ø 2.0 MM X 15 CM
38204-83	Ø 2.5 MM X 15 CM
38204-84	Ø 2.5 MM X 20 CM
38204-85	Ø 3.0 MM X 15 CM
38204-86	Ø 3.0 MM X 20 CM
38204-87	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38205-81	Ø 1.5 MM X 12 CM
38205-82	Ø 2.0 MM X 15 CM
38205-83	Ø 2.0 MM X 20 CM
38205-84	Ø 2.5 MM X 15 CM
38205-85	Ø 2.5 MM X 20 CM
38205-86	Ø 3.0 MM X 15 CM
38205-87	Ø 3.0 MM X 20 CM
38205-88	Ø 3.0 MM X 25 CM

POWER CANNULA

38206-211	Ø 2 MM X 15 CM
38206-212	Ø 2 MM X 20 CM
38206-213	Ø 2 MM X 25 CM
38206-214	Ø 2 MM X 30 CM
38206-215	Ø 2 MM X 35 CM
38206-216	Ø 3 MM X 15 CM
38206-217	Ø 3 MM X 20 CM
38206-218	Ø 3 MM X 25 CM
38206-219	Ø 3 MM X 30 CM
38206-220	Ø 3 MM X 35 CM
38206-221	Ø 4 MM X 15 CM
38206-222	Ø 4 MM X 20 CM
38206-223	Ø 4 MM X 25 CM
38206-224	Ø 4 MM X 30 CM
38206-225	Ø 4 MM X 35 CM
38206-226	Ø 5 MM X 15 CM
38206-227	Ø 5 MM X 20 CM
38206-228	Ø 5 MM X 25 CM
38206-229	Ø 5 MM X 30 CM
38206-230	Ø 5 MM X 35 CM
38206-231	Ø 6 MM X 15 CM
38206-232	Ø 6 MM X 20 CM
38206-233	Ø 6 MM X 25 CM
38206-234	Ø 6 MM X 30 CM
38206-235	Ø 6 MM X 35 CM



GENERAL SUCTION CANNULAS



MADE IN THE USA

BASKET SHAPED CANNULA

THE CELL FRIENDLY BASKET CANNULA IS A VERSATILE LIPOSCULPTING TOOL. IT CAN OPERATE AS AN EQUALIZER WITHOUT VACUUM TO TREAT ADIPOSE TISSUE AND SKIN IRREGULARITIES IN REVISION LIPOSUCTION AND FAT TRANSFER PROCEDURES. THE BASKET CANNULA IS ALSO ADEPT UNDER VACUUM AS A FAT HARVESTING CANNULA. FOR FINISHING, THIS CANNULA CAN BE USED AS AN INJECTOR FOR EQUALIZING AND ADDING ADIPOSE TISSUE TO SCULPT SMOOTHER CONTOURS .



60CC TOMMEY HUB

38201-81	Ø 3.0 MM X 15 CM
38201-82	Ø 3.0 MM X 20 CM
38201-83	Ø 3.0 MM X 25 CM
38201-84	Ø 4.0 MM X 25 CM
38201-85	Ø 4.0 MM X 30 CM
38201-86	Ø 4.0 MM X 35 CM



HOLE ON TIP-STYLE I

POWER CANNULA

38206-251	Ø 2 MM X 15 CM	38206-260	Ø 3 MM X 35 CM	38206-268	Ø 5 MM X 25 CM
38206-252	Ø 2 MM X 20 CM	38206-261	Ø 4 MM X 15 CM	38206-269	Ø 5 MM X 30 CM
38206-253	Ø 2 MM X 25 CM	38206-262	Ø 4 MM X 20 CM	38206-270	Ø 5 MM X 35 CM
38206-254	Ø 2 MM X 30 CM	38206-263	Ø 4 MM X 25 CM	38206-271	Ø 6 MM X 15 CM
38206-255	Ø 2 MM X 35 CM	38206-264	Ø 4 MM X 30 CM	38206-272	Ø 6 MM X 20 CM
38206-256	Ø 3 MM X 15 CM	38206-265	Ø 4 MM X 35 CM	38206-273	Ø 6 MM X 25 CM
38206-257	Ø 3 MM X 20 CM	38206-266	Ø 5 MM X 15 CM	38206-274	Ø 6 MM X 30 CM
38206-259	Ø 3 MM X 30 CM				



HOLE ON TIP-STYLE II

POWER CANNULA

38206-281	Ø 2 MM X 15 CM	38206-290	Ø 3 MM X 35 CM	38206-298	Ø 5 MM X 25 CM
38206-282	Ø 2 MM X 20 CM	38206-291	Ø 4 MM X 15 CM	38206-299	Ø 5 MM X 30 CM
38206-283	Ø 2 MM X 25 CM	38206-292	Ø 4 MM X 20 CM	38206-300	Ø 5 MM X 35 CM
38206-284	Ø 2 MM X 30 CM	38206-293	Ø 4 MM X 25 CM	38206-301	Ø 6 MM X 15 CM
38206-285	Ø 2 MM X 35 CM	38206-294	Ø 4 MM X 30 CM	38206-302	Ø 6 MM X 20 CM
38206-286	Ø 3 MM X 15 CM	38206-295	Ø 4 MM X 35 CM	38206-303	Ø 6 MM X 25 CM
38206-287	Ø 3 MM X 20 CM	38206-296	Ø 5 MM X 15 CM	38206-304	Ø 6 MM X 30 CM
38206-289	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA



DEEP SHARP CUT DISSECTOR

POWER CANNULA

38206-311	Ø 2 MM X 15 CM	38206-320	Ø 3 MM X 35 CM	38206-328	Ø 5 MM X 25 CM
38206-312	Ø 2 MM X 20 CM	38206-321	Ø 4 MM X 15 CM	38206-329	Ø 5 MM X 30 CM
38206-313	Ø 2 MM X 25 CM	38206-322	Ø 4 MM X 20 CM	38206-330	Ø 5 MM X 35 CM
38206-314	Ø 2 MM X 30 CM	38206-323	Ø 4 MM X 25 CM	38206-331	Ø 6 MM X 15 CM
38206-315	Ø 2 MM X 35 CM	38206-324	Ø 4 MM X 30 CM	38206-332	Ø 6 MM X 20 CM
38206-316	Ø 3 MM X 15 CM	38206-325	Ø 4 MM X 35 CM	38206-333	Ø 6 MM X 25 CM
38206-317	Ø 3 MM X 20 CM	38206-326	Ø 5 MM X 15 CM	38206-334	Ø 6 MM X 30 CM
38206-318	Ø 3 MM X 25 CM	38206-327	Ø 5 MM X 20 CM	38206-335	Ø 6 MM X 35 CM
38206-319	Ø 3 MM X 30 CM				



BI-PORT CANNULA

LUER LOCK (STAINLESS STEEL)

	10CM	15CM	20CM	25CM
2.0 Ø	38205-91	38205-101	38205-111	38205-121
2.5 Ø	38205-92	38205-102	38205-112	38205-122
3.0 Ø	38205-93	38205-103	38205-113	38205-123
3.5 Ø	38205-94	38205-104	38205-114	38205-124
4.0 Ø	38205-95	38205-105	38205-115	38205-125

POWER CANNULA

38206-351	Ø 2 MM X 15 CM	38206-360	Ø 3 MM X 35 CM	38206-368	Ø 5 MM X 25 CM
38206-352	Ø 2 MM X 20 CM	38206-361	Ø 4 MM X 15 CM	38206-369	Ø 5 MM X 30 CM
38206-353	Ø 2 MM X 25 CM	38206-362	Ø 4 MM X 20 CM	38206-370	Ø 5 MM X 35 CM
38206-354	Ø 2 MM X 30 CM	38206-363	Ø 4 MM X 25 CM	38206-371	Ø 6 MM X 15 CM
38206-355	Ø 2 MM X 35 CM	38206-364	Ø 4 MM X 30 CM	38206-372	Ø 6 MM X 20 CM
38206-356	Ø 3 MM X 15 CM	38206-365	Ø 4 MM X 35 CM	38206-373	Ø 6 MM X 25 CM
38206-357	Ø 3 MM X 20 CM	38206-366	Ø 5 MM X 15 CM	38206-374	Ø 6 MM X 30 CM
38206-358	Ø 3 MM X 25 CM	38206-367	Ø 5 MM X 20 CM	38206-375	Ø 6 MM X 35 CM
38206-359	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



TAPERED TIP SINGLE CENTRAL HOLE

POWER CANNULA

38206-391	Ø 2 MM X 15 CM	38206-400	Ø 3 MM X 35 CM	38206-408	Ø 5 MM X 25 CM
38206-392	Ø 2 MM X 20 CM	38206-401	Ø 4 MM X 15 CM	38206-409	Ø 5 MM X 30 CM
38206-393	Ø 2 MM X 25 CM	38206-402	Ø 4 MM X 20 CM	38206-410	Ø 5 MM X 35 CM
38206-394	Ø 2 MM X 30 CM	38206-403	Ø 4 MM X 25 CM	38206-411	Ø 6 MM X 15 CM
38206-395	Ø 2 MM X 35 CM	38206-404	Ø 4 MM X 30 CM	38206-412	Ø 6 MM X 20 CM
38206-396	Ø 3 MM X 15 CM	38206-405	Ø 4 MM X 35 CM	38206-413	Ø 6 MM X 25 CM
38206-397	Ø 3 MM X 20 CM	38206-406	Ø 5 MM X 15 CM	38206-414	Ø 6 MM X 30 CM
38206-398	Ø 3 MM X 25 CM	38206-407	Ø 5 MM X 20 CM	38206-415	Ø 6 MM X 35 CM
38206-399	Ø 3 MM X 30 CM				



TAPERED TIP TWO SPIRAL HOLE

POWER CANNULA

38206-431	Ø 2 MM X 15 CM	38206-440	Ø 3 MM X 35 CM	38206-448	Ø 5 MM X 25 CM
38206-432	Ø 2 MM X 20 CM	38206-441	Ø 4 MM X 15 CM	38206-449	Ø 5 MM X 30 CM
38206-433	Ø 2 MM X 25 CM	38206-442	Ø 4 MM X 20 CM	38206-450	Ø 5 MM X 35 CM
38206-434	Ø 2 MM X 30 CM	38206-443	Ø 4 MM X 25 CM	38206-451	Ø 6 MM X 15 CM
38206-435	Ø 2 MM X 35 CM	38206-444	Ø 4 MM X 30 CM	38206-452	Ø 6 MM X 20 CM
38206-436	Ø 3 MM X 15 CM	38206-445	Ø 4 MM X 35 CM	38206-453	Ø 6 MM X 25 CM
38206-437	Ø 3 MM X 20 CM	38206-446	Ø 5 MM X 15 CM	38206-454	Ø 6 MM X 30 CM
38206-438	Ø 3 MM X 25 CM	38206-447	Ø 5 MM X 20 CM	38206-455	Ø 6 MM X 35 CM
38206-439	Ø 3 MM X 30 CM				



FIG. 3F

POWER CANNULA

38206-471	Ø 2 MM X 15 CM	38206-480	Ø 3 MM X 35 CM	38206-488	Ø 5 MM X 25 CM
38206-472	Ø 2 MM X 20 CM	38206-481	Ø 4 MM X 15 CM	38206-489	Ø 5 MM X 30 CM
38206-473	Ø 2 MM X 25 CM	38206-482	Ø 4 MM X 20 CM	38206-490	Ø 5 MM X 35 CM
38206-474	Ø 2 MM X 30 CM	38206-483	Ø 4 MM X 25 CM	38206-491	Ø 6 MM X 15 CM
38206-475	Ø 2 MM X 35 CM	38206-484	Ø 4 MM X 30 CM	38206-492	Ø 6 MM X 20 CM
38206-476	Ø 3 MM X 15 CM	38206-485	Ø 4 MM X 35 CM	38206-493	Ø 6 MM X 25 CM
38206-477	Ø 3 MM X 20 CM	38206-486	Ø 5 MM X 15 CM	38206-494	Ø 6 MM X 30 CM
38206-478	Ø 3 MM X 25 CM	38206-487	Ø 5 MM X 20 CM	38206-495	Ø 6 MM X 35 CM
38206-479	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK

LUER LOCK (STAINLESS STEEL)

ULTRA LUER LOCK

THREADED

60CC TOMMEY HUB

POWER CANNULA



FIG. 2ECH

POWER CANNULA

38206-501	Ø 2 MM X 15 CM	38206-510	Ø 3 MM X 35 CM	38206-518	Ø 5 MM X 25 CM
38206-502	Ø 2 MM X 20 CM	38206-511	Ø 4 MM X 15 CM	38206-519	Ø 5 MM X 30 CM
38206-503	Ø 2 MM X 25 CM	38206-512	Ø 4 MM X 20 CM	38206-520	Ø 5 MM X 35 CM
38206-504	Ø 2 MM X 30 CM	38206-513	Ø 4 MM X 25 CM	38206-521	Ø 6 MM X 15 CM
38206-505	Ø 2 MM X 35 CM	38206-514	Ø 4 MM X 30 CM	38206-522	Ø 6 MM X 20 CM
38206-506	Ø 3 MM X 15 CM	38206-515	Ø 4 MM X 35 CM	38206-523	Ø 6 MM X 25 CM
38206-507	Ø 3 MM X 20 CM	38206-516	Ø 5 MM X 15 CM	38206-524	Ø 6 MM X 30 CM
38206-508	Ø 3 MM X 25 CM	38206-517	Ø 5 MM X 20 CM	38206-525	Ø 6 MM X 35 CM
38206-509	Ø 3 MM X 30 CM				



TWO SQUARE-IN LINE

LUER LOCK (STAINLESS STEEL)

7CM	10CM	15CM	20CM	25CM
2.0 Ø 38205-131	38205-141	38205-151	38205-161	38205-171
2.5 Ø 38205-132	38205-142	38205-152	38205-162	38205-172
3.0 Ø 38205-133	38205-143	38205-153	38205-163	38205-173
3.5 Ø 38205-134	38205-144	38205-154	38205-164	38205-174
4.0 Ø 38205-135	38205-145	38205-155	38205-165	38205-175

POWER CANNULA

38206-551	Ø 2 MM X 15 CM	38206-560	Ø 3 MM X 35 CM	38206-568	Ø 5 MM X 25 CM
38206-552	Ø 2 MM X 20 CM	38206-561	Ø 4 MM X 15 CM	38206-569	Ø 5 MM X 30 CM
38206-553	Ø 2 MM X 25 CM	38206-562	Ø 4 MM X 20 CM	38206-570	Ø 5 MM X 35 CM
38206-554	Ø 2 MM X 30 CM	38206-563	Ø 4 MM X 25 CM	38206-571	Ø 6 MM X 15 CM
38206-555	Ø 2 MM X 35 CM	38206-564	Ø 4 MM X 30 CM	38206-572	Ø 6 MM X 20 CM
38206-556	Ø 3 MM X 15 CM	38206-565	Ø 4 MM X 35 CM	38206-573	Ø 6 MM X 25 CM
38206-557	Ø 3 MM X 20 CM	38206-566	Ø 5 MM X 15 CM	38206-574	Ø 6 MM X 30 CM
38206-558	Ø 3 MM X 25 CM	38206-567	Ø 5 MM X 20 CM	38206-575	Ø 6 MM X 35 CM
38206-559	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



ONE CENTRAL PORT—TWO PARALLEL PORTS

POWER CANNULA

38206-591	Ø 2 MM X 15 CM	38206-600	Ø 3 MM X 35 CM	38206-608	Ø 5 MM X 25 CM
38206-592	Ø 2 MM X 20 CM	38206-601	Ø 4 MM X 15 CM	38206-609	Ø 5 MM X 30 CM
38206-593	Ø 2 MM X 25 CM	38206-602	Ø 4 MM X 20 CM	38206-610	Ø 5 MM X 35 CM
38206-594	Ø 2 MM X 30 CM	38206-603	Ø 4 MM X 25 CM	38206-611	Ø 6 MM X 15 CM
38206-595	Ø 2 MM X 35 CM	38206-604	Ø 4 MM X 30 CM	38206-612	Ø 6 MM X 20 CM
38206-596	Ø 3 MM X 15 CM	38206-605	Ø 4 MM X 35 CM	38206-613	Ø 6 MM X 25 CM
38206-597	Ø 3 MM X 20 CM	38206-606	Ø 5 MM X 15 CM	38206-614	Ø 6 MM X 30 CM
38206-598	Ø 3 MM X 25 CM	38206-607	Ø 5 MM X 20 CM	38206-615	Ø 6 MM X 35 CM
38206-599	Ø 3 MM X 30 CM				



TWO PARALLEL PORTS—ONE CENTRAL PORT

POWER CANNULA


38206-651	Ø 2 MM X 15 CM	38206-660	Ø 3 MM X 35 CM	38206-668	Ø 5 MM X 25 CM
38206-652	Ø 2 MM X 20 CM	38206-661	Ø 4 MM X 15 CM	38206-669	Ø 5 MM X 30 CM
38206-653	Ø 2 MM X 25 CM	38206-662	Ø 4 MM X 20 CM	38206-670	Ø 5 MM X 35 CM
38206-654	Ø 2 MM X 30 CM	38206-663	Ø 4 MM X 25 CM	38206-671	Ø 6 MM X 15 CM
38206-655	Ø 2 MM X 35 CM	38206-664	Ø 4 MM X 30 CM	38206-672	Ø 6 MM X 20 CM
38206-656	Ø 3 MM X 15 CM	38206-665	Ø 4 MM X 35 CM	38206-673	Ø 6 MM X 25 CM
38206-657	Ø 3 MM X 20 CM	38206-666	Ø 5 MM X 15 CM	38206-674	Ø 6 MM X 30 CM
38206-658	Ø 3 MM X 25 CM	38206-667	Ø 5 MM X 20 CM	38206-675	Ø 6 MM X 35 CM
38206-659	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA

FIG. 3ME



LUER LOCK (STAINLESS STEEL)

	7CM	10CM	15CM	20CM	25CM
2.0 Ø	38205-181	38205-191	38205-201	38205-211	38205-221
2.5 Ø	38205-182	38205-192	38205-202	38205-212	38205-222
3.0 Ø	38205-183	38205-193	38205-203	38205-213	38205-223
3.5 Ø	38205-184	38205-194	38205-204	38205-214	38205-224
4.0 Ø	38205-185	38205-195	38205-205	38205-215	38205-225

POWER CANNULA

38206-681	Ø 2 MM X 15 CM	38206-690	Ø 3 MM X 35 CM	38206-698	Ø 5 MM X 25 CM
38206-682	Ø 2 MM X 20 CM	38206-691	Ø 4 MM X 15 CM	38206-699	Ø 5 MM X 30 CM
38206-683	Ø 2 MM X 25 CM	38206-692	Ø 4 MM X 20 CM	38206-700	Ø 5 MM X 35 CM
38206-684	Ø 2 MM X 30 CM	38206-693	Ø 4 MM X 25 CM	38206-701	Ø 6 MM X 15 CM
38206-685	Ø 2 MM X 35 CM	38206-694	Ø 4 MM X 30 CM	38206-702	Ø 6 MM X 20 CM
38206-686	Ø 3 MM X 15 CM	38206-695	Ø 4 MM X 35 CM	38206-703	Ø 6 MM X 25 CM
38206-687	Ø 3 MM X 20 CM	38206-696	Ø 5 MM X 15 CM	38206-704	Ø 6 MM X 30 CM
38206-688	Ø 3 MM X 25 CM	38206-697	Ø 5 MM X 20 CM	38206-705	Ø 6 MM X 35 CM
38206-689	Ø 3 MM X 30 CM				



FIG. 3HR

POWER CANNULA

38206-721	Ø 2 MM X 15 CM	38206-730	Ø 3 MM X 35 CM	38206-738	Ø 5 MM X 25 CM
38206-722	Ø 2 MM X 20 CM	38206-731	Ø 4 MM X 15 CM	38206-739	Ø 5 MM X 30 CM
38206-723	Ø 2 MM X 25 CM	38206-732	Ø 4 MM X 20 CM	38206-740	Ø 5 MM X 35 CM
38206-724	Ø 2 MM X 30 CM	38206-733	Ø 4 MM X 25 CM	38206-741	Ø 6 MM X 15 CM
38206-725	Ø 2 MM X 35 CM	38206-734	Ø 4 MM X 30 CM	38206-742	Ø 6 MM X 20 CM
38206-726	Ø 3 MM X 15 CM	38206-735	Ø 4 MM X 35 CM	38206-743	Ø 6 MM X 25 CM
38206-727	Ø 3 MM X 20 CM	38206-736	Ø 5 MM X 15 CM	38206-744	Ø 6 MM X 30 CM
38206-728	Ø 3 MM X 25 CM	38206-737	Ø 5 MM X 20 CM	38206-745	Ø 6 MM X 35 CM
38206-729	Ø 3 MM X 30 CM				



FIG. RIZ

POWER CANNULA


38206-751	Ø 2 MM X 15 CM	38206-760	Ø 3 MM X 35 CM	38206-768	Ø 5 MM X 25 CM
38206-752	Ø 2 MM X 20 CM	38206-761	Ø 4 MM X 15 CM	38206-769	Ø 5 MM X 30 CM
38206-753	Ø 2 MM X 25 CM	38206-762	Ø 4 MM X 20 CM	38206-770	Ø 5 MM X 35 CM
38206-754	Ø 2 MM X 30 CM	38206-763	Ø 4 MM X 25 CM	38206-771	Ø 6 MM X 15 CM
38206-755	Ø 2 MM X 35 CM	38206-764	Ø 4 MM X 30 CM	38206-772	Ø 6 MM X 20 CM
38206-756	Ø 3 MM X 15 CM	38206-765	Ø 4 MM X 35 CM	38206-773	Ø 6 MM X 25 CM
38206-757	Ø 3 MM X 20 CM	38206-766	Ø 5 MM X 15 CM	38206-774	Ø 6 MM X 30 CM
38206-758	Ø 3 MM X 25 CM	38206-767	Ø 5 MM X 20 CM	38206-775	Ø 6 MM X 35 CM
38206-759	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



FIG. 5AI

POWER CANNULA

38206-721	Ø 2 MM X 15 CM	38206-730	Ø 3 MM X 35 CM	38206-738	Ø 5 MM X 25 CM
38206-722	Ø 2 MM X 20 CM	38206-731	Ø 4 MM X 15 CM	38206-739	Ø 5 MM X 30 CM
38206-723	Ø 2 MM X 25 CM	38206-732	Ø 4 MM X 20 CM	38206-740	Ø 5 MM X 35 CM
38206-724	Ø 2 MM X 30 CM	38206-733	Ø 4 MM X 25 CM	38206-741	Ø 6 MM X 15 CM
38206-725	Ø 2 MM X 35 CM	38206-734	Ø 4 MM X 30 CM	38206-742	Ø 6 MM X 20 CM
38206-726	Ø 3 MM X 15 CM	38206-735	Ø 4 MM X 35 CM	38206-743	Ø 6 MM X 25 CM
38206-727	Ø 3 MM X 20 CM	38206-736	Ø 5 MM X 15 CM	38206-744	Ø 6 MM X 30 CM
38206-728	Ø 3 MM X 25 CM	38206-737	Ø 5 MM X 20 CM	38206-745	Ø 6 MM X 35 CM
38206-729	Ø 3 MM X 30 CM				



FIG. MR

POWER CANNULA

38206-751	Ø 2 MM X 15 CM	38206-760	Ø 3 MM X 35 CM	38206-768	Ø 5 MM X 25 CM
38206-752	Ø 2 MM X 20 CM	38206-761	Ø 4 MM X 15 CM	38206-769	Ø 5 MM X 30 CM
38206-753	Ø 2 MM X 25 CM	38206-762	Ø 4 MM X 20 CM	38206-770	Ø 5 MM X 35 CM
38206-754	Ø 2 MM X 30 CM	38206-763	Ø 4 MM X 25 CM	37206-771	Ø 6 MM X 15 CM
38206-755	Ø 2 MM X 35 CM	38206-764	Ø 4 MM X 30 CM	38206-772	Ø 6 MM X 20 CM
38206-756	Ø 3 MM X 15 CM	38206-765	Ø 4 MM X 35 CM	38206-773	Ø 6 MM X 25 CM
38206-757	Ø 3 MM X 20 CM	38206-766	Ø 5 MM X 15 CM	38206-774	Ø 6 MM X 30 CM
38206-758	Ø 3 MM X 25 CM	38206-767	Ø 5 MM X 20 CM	38206-775	Ø 6 MM X 35 CM
38206-759	Ø 3 MM X 30 CM				



FIG. P3G

POWER CANNULA

38206-791	Ø 2 MM X 15 CM	38206-800	Ø 3 MM X 35 CM	38206-808	Ø 5 MM X 25 CM
38206-792	Ø 2 MM X 20 CM	38206-801	Ø 4 MM X 15 CM	38206-809	Ø 5 MM X 30 CM
38206-793	Ø 2 MM X 25 CM	38206-802	Ø 4 MM X 20 CM	38206-810	Ø 5 MM X 35 CM
38206-794	Ø 2 MM X 30 CM	38206-803	Ø 4 MM X 25 CM	38206-811	Ø 6 MM X 15 CM
38206-795	Ø 2 MM X 35 CM	38206-804	Ø 4 MM X 30 CM	38206-812	Ø 6 MM X 20 CM
38206-796	Ø 3 MM X 15 CM	38206-805	Ø 4 MM X 35 CM	38206-813	Ø 6 MM X 25 CM
38206-797	Ø 3 MM X 20 CM	38206-806	Ø 5 MM X 15 CM	38206-814	Ø 6 MM X 30 CM
38206-798	Ø 3 MM X 25 CM	38206-807	Ø 5 MM X 20 CM	38206-815	Ø 6 MM X 35 CM
38206-799	Ø 3 MM X 30 CM				



GENERAL SUCTION CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



FIG. Pv

POWER CANNULA

38206-831	Ø 2 MM X 15 CM	38206-840	Ø 3 MM X 35 CM	38206-848	Ø 5 MM X 25 CM
38206-832	Ø 2 MM X 20 CM	38206-841	Ø 4 MM X 15 CM	38206-849	Ø 5 MM X 30 CM
38206-833	Ø 2 MM X 25 CM	38206-842	Ø 4 MM X 20 CM	38206-850	Ø 5 MM X 35 CM
38206-834	Ø 2 MM X 30 CM	38206-843	Ø 4 MM X 25 CM	38206-851	Ø 6 MM X 15 CM
38206-835	Ø 2 MM X 35 CM	38206-844	Ø 4 MM X 30 CM	38206-852	Ø 6 MM X 20 CM
38206-836	Ø 3 MM X 15 CM	38206-845	Ø 4 MM X 35 CM	38206-853	Ø 6 MM X 25 CM
38206-837	Ø 3 MM X 20 CM	38206-846	Ø 5 MM X 15 CM	38206-854	Ø 6 MM X 30 CM
38206-838	Ø 3 MM X 25 CM	38206-847	Ø 5 MM X 20 CM	38206-855	Ø 6 MM X 35 CM
38206-839	Ø 3 MM X 30 CM				



FIG. LE

POWER CANNULA

38206-861	Ø 2 MM X 15 CM	38206-870	Ø 3 MM X 35 CM	38206-878	Ø 5 MM X 25 CM
38206-862	Ø 2 MM X 20 CM	38206-871	Ø 4 MM X 15 CM	38206-879	Ø 5 MM X 30 CM
38206-863	Ø 2 MM X 25 CM	38206-872	Ø 4 MM X 20 CM	38206-880	Ø 5 MM X 35 CM
38206-864	Ø 2 MM X 30 CM	38206-873	Ø 4 MM X 25 CM	38206-881	Ø 6 MM X 15 CM
38206-865	Ø 2 MM X 35 CM	38206-874	Ø 4 MM X 30 CM	38206-882	Ø 6 MM X 20 CM
38206-866	Ø 3 MM X 15 CM	38206-875	Ø 4 MM X 35 CM	38206-883	Ø 6 MM X 25 CM
38206-867	Ø 3 MM X 20 CM	38206-876	Ø 5 MM X 15 CM	38206-884	Ø 6 MM X 30 CM
38206-868	Ø 3 MM X 25 CM	38206-877	Ø 5 MM X 20 CM	38206-885	Ø 6 MM X 35 CM
38206-868	Ø 3 MM X 30 CM				



HARVESTER CANNULAS



MADE IN THE USA

FAT HARVESTING FOR RE-GRAFTING IS A DELICATE, DIFFICULT, AND TIME CONSUMING PROCESS. BUT WITH ELMED'S SPECIFICALLY DESIGNED HARVESTING CANNULAS WILL HELP YOU TO ACHIEVE THIS GOAL IN SUPER EFFICIENT WAY. MICRO FAT HARVESTING - NO PROBLEM WITH ELMED CANNULAS.

MULTIPOINT SPEED HARVESTER- BOTH END SHARP

A MULTIPOINT HARVESTER WITH ELEVATED 1MM MICROPORTS WHICH YIELD SMALL ENOUGH ALIQUOTS OF FAT TISSUE TO ALLOW FOR EASY AND SMOOTH INJECTION WITH THE TULIP 0.7, 0.9 AND 1.2MM INJECTORS.



60CC TOMMEY HUB

38301-01	Ø 2.5 MM X 15 CM
38301-02	Ø 2.5 MM X 20 CM
38301-03	Ø 3.0 MM X 15 CM
38301-04	Ø 3.0 MM X 20 CM
38301-05	Ø 3.0 MM X 25 CM
38301-06	Ø 4.0 MM X 25 CM
38301-07	Ø 4.0 MM X 30 CM
38301-08	Ø 4.0 MM X 35 CM
38301-09	Ø 4.5 MM X 30 CM
38301-10	Ø 4.5 MM X 35 CM

THREADED HUB

38301-01	Ø 2.5 MM X 15 CM
38301-02	Ø 2.5 MM X 20 CM
38301-03	Ø 3.0 MM X 15 CM
38301-04	Ø 3.0 MM X 20 CM
38301-05	Ø 3.0 MM X 25 CM
38301-06	Ø 4.0 MM X 25 CM
38301-07	Ø 4.0 MM X 30 CM
38301-08	Ø 4.0 MM X 35 CM
38301-09	Ø 4.5 MM X 30 CM
38301-10	Ø 4.5 MM X 35 CM

ULTRA LUER LOCK

38303-01	Ø 2.0 MM X 15 CM
38303-02	Ø 2.5 MM X 15 CM
38303-03	Ø 2.5 MM X 20 CM
38303-04	Ø 3.0 MM X 15 CM
38303-05	Ø 3.0 MM X 20 CM
38303-06	Ø 3.0 MM X 25 CM

LUER LOCK

38304-01	Ø 2.0 MM X 15 CM
38304-02	Ø 2.5 MM X 15 CM
38304-03	Ø 2.5 MM X 20 CM
38304-04	Ø 3.0 MM X 15 CM
38304-05	Ø 3.0 MM X 20 CM
38304-06	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38305-01	Ø 2.0 MM X 15 CM
38305-02	Ø 2.5 MM X 15 CM
38305-03	Ø 2.5 MM X 20 CM
38305-04	Ø 3.0 MM X 15 CM
38305-05	Ø 3.0 MM X 20 CM
38305-06	Ø 3.0 MM X 25 CM

POWER CANNULA


38306-01	Ø 2 MM X 15 CM
38306-02	Ø 2 MM X 20 CM
38306-03	Ø 2 MM X 25 CM
38306-04	Ø 2 MM X 30 CM
38306-05	Ø 2 MM X 35 CM
38306-06	Ø 3 MM X 15 CM
38306-07	Ø 3 MM X 20 CM
38306-08	Ø 3 MM X 25 CM
38306-09	Ø 3 MM X 30 CM
38306-10	Ø 3 MM X 35 CM
38306-11	Ø 4 MM X 15 CM
38306-12	Ø 4 MM X 20 CM
38306-13	Ø 4 MM X 25 CM
38306-14	Ø 4 MM X 30 CM
38306-15	Ø 4 MM X 35 CM
38306-16	Ø 5 MM X 15 CM
38306-17	Ø 5 MM X 20 CM
38306-18	Ø 5 MM X 25 CM
38306-19	Ø 5 MM X 30 CM
38306-20	Ø 5 MM X 35 CM
38306-21	Ø 6 MM X 15 CM
38306-22	Ø 6 MM X 20 CM
38306-23	Ø 6 MM X 25 CM
38306-24	Ø 6 MM X 30 CM
38306-25	Ø 6 MM X 35 CM



HARVESTER CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA

MULTIPOINT SPEED HARVESTER-ONE END SHARP

A MULTIPOINT HARVESTER WITH ELEVATED 1MM MICROPORTS WHICH YIELD SMALL ENOUGH ALIQUOTS OF FAT TISSUE TO ALLOW FOR EASY AND SMOOTH INJECTION WITH THE TULIP 0.7, 0.9 AND 1.2MM INJECTORS.



60CC TOMMEY HUB

38301-11	Ø 2.5 MM X 15 CM
38301-12	Ø 2.5 MM X 20 CM
38301-13	Ø 3.0 MM X 15 CM
38301-14	Ø 3.0 MM X 20 CM
38301-15	Ø 3.0 MM X 25 CM
38301-16	Ø 4.0 MM X 25 CM
38301-17	Ø 4.0 MM X 30 CM
38301-18	Ø 4.0 MM X 35 CM
38301-19	Ø 4.5 MM X 30 CM
38301-20	Ø 4.5 MM X 35 CM

THREADED HUB

38302-11	Ø 3.0 MM X 15 CM
38302-12	Ø 3.0 MM X 25 CM
38302-13	Ø 4.0 MM X 15 CM
38302-14	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38303-11	Ø 2.0 MM X 15 CM
38303-12	Ø 2.5 MM X 15 CM
38303-13	Ø 2.5 MM X 20 CM
38303-14	Ø 3.0 MM X 15 CM
38303-15	Ø 3.0 MM X 20 CM
38303-16	Ø 3.0 MM X 25 CM

LUER LOCK

38304-11	Ø 2.0 MM X 15 CM
38304-12	Ø 2.5 MM X 15 CM
38304-13	Ø 2.5 MM X 20 CM
38304-14	Ø 3.0 MM X 15 CM
38304-15	Ø 3.0 MM X 20 CM
38304-16	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38305-11	Ø 2.0 MM X 15 CM
38305-12	Ø 2.5 MM X 15 CM
38305-13	Ø 2.5 MM X 20 CM
38305-14	Ø 3.0 MM X 15 CM
38305-15	Ø 3.0 MM X 20 CM
38305-16	Ø 3.0 MM X 25 CM

POWER CANNULA


38306-11	Ø 2 MM X 15 CM
38306-12	Ø 2 MM X 20 CM
38306-13	Ø 2 MM X 25 CM
38306-14	Ø 2 MM X 30 CM
38306-15	Ø 2 MM X 35 CM
39306-16	Ø 3 MM X 15 CM
38306-17	Ø 3 MM X 20 CM
38306-18	Ø 3 MM X 25 CM
38306-19	Ø 3 MM X 30 CM
38306-20	Ø 3 MM X 35 CM
38306-21	Ø 4 MM X 15 CM
38306-22	Ø 4 MM X 20 CM
38306-23	Ø 4 MM X 25 CM
38306-24	Ø 4 MM X 30 CM
38306-25	Ø 4 MM X 35 CM
38306-26	Ø 5 MM X 15 CM
38306-27	Ø 5 MM X 20 CM
38306-28	Ø 5 MM X 25 CM
38306-29	Ø 5 MM X 30 CM
38306-30	Ø 5 MM X 35 CM
38306-31	Ø 6 MM X 15 CM
38306-32	Ø 6 MM X 20 CM
38306-33	Ø 6 MM X 25 CM
38306-34	Ø 6 MM X 30 CM
38306-35	Ø 6 MM X 35 CM



HARVESTER CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



MULTI PORT HARVESTER - BOTH END BEVEL

A SOPHISTICATED, PRECISION-DESIGNED HARVESTER THAT IS TRULY PATIENT FRIENDLY. THE PORTS IN THE SORESENSEN HARVESTER ARE MACHINED WITH ALTERNATING 60 DEGREE CUTTING EDGES. THIS UNIQUE DESIGN IS LESS TRAUMATIC THAN OTHER MULTI-PORT HARVESTERS AND DELIVERS A FINE-GRAINED, HIGH FAT CONTENT/PURITY YIELD THAT CAN EASILY PASS THROUGH 1.2MM, 0.9MM AND 0.7MM BLUNT INJECTORS

60CC TOMMEY HUB

38301-21	Ø 2.5 MM X 15 CM
38301-22	Ø 2.5 MM X 20 CM
38301-23	Ø 3.0 MM X 15 CM
38301-24	Ø 3.0 MM X 20 CM
38301-25	Ø 3.0 MM X 25 CM
38301-26	Ø 4.0 MM X 25 CM
38301-27	Ø 4.0 MM X 30 CM
38301-28	Ø 4.0 MM X 35 CM
38301-29	Ø 4.5 MM X 30 CM
38301-30	Ø 4.5 MM X 35 CM

THREADED HUB

38302-21	Ø 3.0 MM X 15 CM
38302-22	Ø 3.0 MM X 25 CM
38302-23	Ø 4.0 MM X 15 CM
38302-24	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38303-21	Ø 2.0 MM X 15 CM
38303-22	Ø 2.5 MM X 15 CM
38303-23	Ø 2.5 MM X 20 CM
38303-24	Ø 3.0 MM X 15 CM
38303-25	Ø 3.0 MM X 20 CM
38303-26	Ø 3.0 MM X 25 CM

LUER LOCK

38304-21	Ø 2.0 MM X 15 CM
38304-22	Ø 2.5 MM X 15 CM
38304-23	Ø 2.5 MM X 20 CM
38304-24	Ø 3.0 MM X 15 CM
38304-25	Ø 3.0 MM X 20 CM
38304-26	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38305-21	Ø 2.0 MM X 15 CM
38305-22	Ø 2.5 MM X 15 CM
38305-23	Ø 2.5 MM X 20 CM
38305-24	Ø 3.0 MM X 15 CM
38305-25	Ø 3.0 MM X 20 CM
38305-26	Ø 3.0 MM X 25 CM

POWER CANNULA

38306-21	Ø 2 MM X 15 CM
38306-22	Ø 2 MM X 20 CM
38306-23	Ø 2 MM X 25 CM
38306-24	Ø 2 MM X 30 CM
38306-25	Ø 2 MM X 35 CM
38306-26	Ø 3 MM X 15 CM
38306-27	Ø 3 MM X 20 CM
38306-28	Ø 3 MM X 25 CM
38306-29	Ø 3 MM X 30 CM
38306-30	Ø 3 MM X 35 CM
38306-31	Ø 4 MM X 15 CM
38306-32	Ø 4 MM X 20 CM
38306-33	Ø 4 MM X 25 CM
38306-34	Ø 4 MM X 30 CM
38306-35	Ø 4 MM X 35 CM
38306-36	Ø 5 MM X 15 CM
38306-37	Ø 5 MM X 20 CM
38306-38	Ø 5 MM X 25 CM
38306-39	Ø 5 MM X 30 CM
38306-40	Ø 5 MM X 35 CM
38306-41	Ø 6 MM X 15 CM
38306-42	Ø 6 MM X 20 CM
38306-43	Ø 6 MM X 25 CM
38306-44	Ø 6 MM X 30 CM
38306-45	Ø 6 MM X 35 CM




INJECTOR CANNULAS



MADE IN THE USA

MICRO INJECTORS



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



SHARP EXTRACTOR/INJECTOR—THIS TIP WORKS WELL FOR NEEDLE EXTRACTION OF SOFT TISSUE AND AS A SOFT TISSUE REINJECTION NEEDLE

ULTRA LUER LOCK

38403-01	Ø 1.0 MM X 5 CM
38403-02	Ø 1.0 MM X 7 CM
38403-03	Ø 2.0 MM X 15 CM
38403-04	Ø 2.5 MM X 15 CM
38403-05	Ø 3.0 MM X 15 CM

ULTRA LUER LOCK

38404-01	Ø 1.0 MM X 5 CM
38404-02	Ø 1.0 MM X 7 CM
38404-03	Ø 2.0 MM X 15 CM
38404-04	Ø 2.5 MM X 15 CM
38404-05	Ø 3.0 MM X 15 CM

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1 Ø	38405-01	38405-11	38405-21	38405-31	38405-41
1.4 Ø	38405-02	38405-12	38405-22	38405-32	38405-42
1.8 Ø	38405-03	38405-13	38405-23	38405-33	38405-43
2.0 Ø	38405-04	38405-14	38405-24	38405-34	38405-44



BLUNT EXTRACTOR/INJECTOR —THIS TIP WORKS WELL FOR NEEDLE EXTRACTION OF SOFT TISSUE AND AS A SOFT TISSUE REINJECTION NEEDLE.

ULTRA LUER LOCK

38403-01	Ø 1.0 MM X 5 CM
38403-02	Ø 1.0 MM X 7 CM
38403-03	Ø 2.0 MM X 15 CM
38403-04	Ø 2.5 MM X 15 CM
38403-05	Ø 3.0 MM X 15 CM



BLUNT TIP INJECTOR —ROUNDED BLUNT TIP DESIGN PASSES THROUGH SOFT TISSUE EASILY WITHOUT THE TRAUMA OR POTENTIAL PROBLEMS OF SHARP INSTRUMENTS.

ULTRA LUER LOCK

38403-21	Ø 1.0 MM X 4 CM
38403-22	Ø 1.0 MM X 5 CM
38403-23	Ø 1.0 MM X 7 CM
38403-24	CURVED 12 CM
38403-25	CURVED 15 CM
38403-26	CURVED 20 CM

LUER LOCK

38404-21	Ø 1.0 MM X 4 CM
38404-22	Ø 1.0 MM X 5 CM
38404-23	Ø 1.0 MM X 7 CM
38404-24	CURVED 12 CM
38404-25	CURVED 15 CM
38404-26	CURVED 20 CM

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1 Ø	38405-101	38405-111	38405-121	38405-131	38405-141
1.4 Ø	38405-102	38405-112	38405-122	38405-132	38405-142
1.8 Ø	38405-103	38405-113	38405-123	38405-133	38405-143
2.0 Ø	38405-104	38405-114	38405-124	38405-134	38405-144



INJECTOR CANNULAS



MADE IN THE USA

MICRO INJECTORS



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



FIG. 2B

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1Ø	38405-201	38405-211	38405-221	38405-231	38405-241
1.4Ø	38405-202	38405-212	38405-222	38405-232	38405-242
1.8Ø	38405-203	38405-213	38405-223	38405-233	38405-243
2.0Ø	38405-204	38405-214	38405-224	38405-234	38405-244



FIG. 1C

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1Ø	38405-251	38405-261	38405-271	38405-281	38405-291
1.4Ø	38405-252	38405-262	38405-272	38405-282	38405-292
1.8Ø	38405-253	38405-263	38405-273	38405-283	38405-293
2.0Ø	38405-254	38405-264	38405-274	38405-284	38405-294



FIG. 1P

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1Ø	38405-301	38405-311	38405-321	38405-331	38405-341
1.4Ø	38405-302	38405-312	38405-322	38405-332	38405-342
1.8Ø	38405-303	38405-313	38405-323	38405-333	38405-343
2.0Ø	38405-304	38405-314	38405-324	38405-334	38405-344



FIG. 2E

LUER LOCK (STAINLESS STEEL)


	3CM	5CM	7CM	10CM	12CM
1.1Ø	38405-351	38405-361	38405-371	38405-381	38405-391
1.4Ø	38405-352	38405-362	38405-372	38405-382	38405-392
1.8Ø	38405-353	38405-363	38405-373	38405-383	38405-393
2.0Ø	38405-354	38405-364	38405-374	38405-384	38405-394



INJECTOR CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA

MICRO INJECTORS



COLEMAN LIPOSUCTION CANNULA FIG. MIC I

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM	15CM	20CM	3.5CM
1.1Ø	38405-401	38405-411	38405-421	38405-431	38405-441	38405-451	38405-461	38405-471
1.4Ø	38405-402	38405-412	38405-422	38405-432	38405-442	38405-452	38405-462	38405-472
1.8Ø	38405-403	38405-413	38405-423	38405-433	38405-443	38405-453	38405-463	38405-473
2.0Ø	38405-404	38405-414	38405-424	38405-434	38405-444	38405-454	38405-464	38405-474
2.5Ø	38405-405	38405-415	38405-425	38405-435	38405-445	38405-455	38405-465	38405-475



COLEMAN LIPOSUCTION CANNULA FIG. MIC II

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM	15CM	20CM
1.1Ø	38405-501	38405-511	38405-521	38405-531	38405-541	38405-551	38405-561
1.4Ø	38405-502	38405-512	38405-522	38405-532	38405-542	38405-552	38405-562
1.8Ø	38405-503	38405-513	38405-523	38405-533	38405-543	38405-553	38405-563
2.0Ø	38405-504	38405-514	38405-524	38405-534	38405-544	38405-554	38405-564
2.5Ø	38405-505	38405-515	38405-525	38405-535	38405-545	38405-555	38405-565



COLEMAN LIPOSUCTION CANNULA SPPON FIG. MIC III

THIS TIP WORKS WELL FOR NEEDLE EXTRACTION OF SOFT TISSUE AND AS A SOFT TISSUE RE INJECTION NEEDLE.

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM	15CM	20CM
1.1Ø	38405-601	38405-611	38405-621	38405-631	38405-641	38405-651	38405-661
1.4Ø	38405-602	38405-612	38405-622	38405-632	38405-642	38405-652	38405-662
1.8Ø	38405-603	38405-613	38405-623	38405-633	38405-643	38405-653	38405-663
2.0Ø	38405-604	38405-614	38405-624	38405-634	38405-644	38405-654	38405-664
2.5Ø	38405-605	38405-615	38405-625	38405-635	38405-645	38405-655	38405-665



LITTLE INJECTOR

OPEN AND MORE POINTED BLUNT TIP IS VERY EFFECTIVE IN DENSE AND FIBROUS TISSUE.

Ultra Luer Lock

38403-51	Ø 1.0 mm X 4 cm
38403-52	Ø 1.0 mm X 5 cm
38403-53	Ø 1.0 mm X 7 cm

Luer Lock

38404-51	Ø 1.0 mm X 4 cm
38404-52	Ø 1.0 mm X 5 cm
38404-53	Ø 1.0 mm X 7 cm

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM
1.1Ø	38405-701	38405-711	38405-721	38405-731	38405-741
1.4Ø	38405-702	38405-712	38405-722	38405-732	38405-742
1.8Ø	38405-703	38405-713	38405-723	38405-733	38405-743
2.0Ø	38405-704	38405-714	38405-724	38405-734	38405-744



INJECTOR CANNULAS



MADE IN THE USA



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA

MICRO INJECTORS



FIG. 1B

LUER LOCK (STAINLESS STEEL)

	3CM	5CM	7CM	10CM	12CM	15CM	20CM
1.1Ø	38405-801	38405-811	38405-821	38405-831	38405-841	38405-851	38405-861
1.4Ø	38405-802	38405-812	38405-822	38405-832	38405-842	38405-852	38405-862
1.8Ø	38405-803	38405-813	38405-823	38405-833	38405-843	38405-853	38405-863
2.0Ø	38405-804	38405-814	38405-824	38405-834	38405-844	38405-854	38405-864
2.5Ø	38405-805	38405-815	38405-825	38405-835	38405-845	38405-855	38405-865

CURVED INJECTORS



NEW **CURVED INJECTORS** ARE DESIGNED ESPECIALLY FOR FAT TRANSFER PROCEDURES TO THE BODY. THEIR UNIQUE SHAPE MIMICS THE ANATOMICAL CURVATURE OF THE FEMALE TORSO AND, THUS, EASILY PRODUCES THE NATURAL-LOOKING RESULTS THAT PATIENTS DESIRE. CURVED INJECTION CANNULAS AMELIORATE THE FAT TRANSFER PROCESS BY PROVIDING A DIRECT AND PRECISE ROUTE TO THE INJECTION SITE AND BY LESSENING THE TIME REQUIRED TO CORRECTLY SCULPT CURVES.




SPECIALTY CANNULAS



MADE IN THE USA

ELMED'S SPECIALTY CANNULAS ARE USEFUL TOOLS FOR DIFFERENT TO OPERATION ARE SUCH AS FACE OR NECK.



LUER LOCK LUER LOCK (STAINLESS STEEL) ULTRA LUER LOCK THREADED 60CC TOMMEY HUB POWER CANNULA



V- DISSECTOR

THIS CANNULA IS SPECIFICALLY DESIGNED FOR USE IN THE SUPERFICIAL COMPARTMENT OF THE TRUNCAL AREA. THE CUTTING NOTCH IS VERY USEFUL FOR UNDERMINING IN FACIAL PROCEDURES AND IN DENSE, FIBROUS TISSUE.

60CC TOMMEY HUB

38501-01	Ø 2.5 MM X 15 CM
38501-02	Ø 2.5 MM X 20 CM
38501-03	Ø 3.0 MM X 15 CM
38501-04	Ø 3.0 MM X 20 CM
38501-05	Ø 3.0 MM X 25 CM
38501-06	Ø 4.0 MM X 25 CM
38501-07	Ø 4.0 MM X 30 CM
38501-08	Ø 4.0 MM X 35 CM
38501-09	Ø 4.5 MM X 30 CM
38501-10	Ø 4.5 MM X 35 CM

THREADED HUB

38502-01	Ø 3.0 MM X 15 CM
38502-02	Ø 3.0 MM X 25 CM
38502-03	Ø 4.0 MM X 15 CM
38502-04	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38503-01	Ø 2.0 MM X 15 CM
38503-02	Ø 2.5 MM X 15 CM
38503-03	Ø 2.5 MM X 20 CM
38503-04	Ø 3.0 MM X 15 CM
38503-05	Ø 3.0 MM X 20 CM
38503-06	Ø 3.0 MM X 25 CM

LUER LOCK

38504-01	Ø 2.0 MM X 15 CM
38504-02	Ø 2.5 MM X 15 CM
38504-03	Ø 2.5 MM X 20 CM
38504-04	Ø 3.0 MM X 15 CM
38504-05	Ø 3.0 MM X 20 CM
38504-06	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38505-01	Ø 2.0 MM X 15 CM
38505-02	Ø 2.5 MM X 15 CM
38505-03	Ø 2.5 MM X 20 CM
38505-04	Ø 3.0 MM X 15 CM
38505-05	Ø 3.0 MM X 20 CM
38505-06	Ø 3.0 MM X 25 CM

POWER CANNULA

38506-01	Ø 2 MM X 15 CM
38506-02	Ø 2 MM X 20 CM
38506-03	Ø 2 MM X 25 CM
38506-04	Ø 2 MM X 30 CM
38506-05	Ø 2 MM X 35 CM
38506-06	Ø 3 MM X 15 CM
38506-07	Ø 3 MM X 20 CM
38506-08	Ø 3 MM X 25 CM
38506-09	Ø 3 MM X 30 CM
38506-10	Ø 3 MM X 35 CM
38506-11	Ø 4 MM X 15 CM
38506-12	Ø 4 MM X 20 CM
38506-13	Ø 4 MM X 25 CM
38506-14	Ø 4 MM X 30 CM
38506-15	Ø 4 MM X 35 CM
38506-16	Ø 5 MM X 15 CM
38506-17	Ø 5 MM X 20 CM
38506-18	Ø 5 MM X 25 CM
38506-19	Ø 5 MM X 30 CM
38506-20	Ø 5 MM X 35 CM
38506-21	Ø 6 MM X 15 CM
38506-22	Ø 6 MM X 20 CM
38506-23	Ø 6 MM X 25 CM
38506-24	Ø 6 MM X 30 CM
38506-25	Ø 6 MM X 35 CM



SPECIALTY CANNULAS



MADE IN THE USA



FLAP

THE FLAP IS PRIMARILY USED AS AN UNDERMINING DISSECTOR, AND IN SMALLER DIAMETERS, FOR FACIAL TISSUE REMOVAL. THE FLAP CANNULA IN A 2.1MM DIAMETER IS VERY USEFUL FOR UNDERMINING THE NASOLABIAL FOLD PRIOR TO AUTOLOGOUS TISSUE INJECTIONS, AS WELL AS FOR TISSUE REMOVAL IN THIS ANATOMIC AREA.

60CC TOMMEY HUB

38501-11	Ø 2.5 MM X 15 CM
38501-12	Ø 2.5 MM X 20 CM
38501-13	Ø 3.0 MM X 15 CM
38501-14	Ø 3.0 MM X 20 CM
38501-15	Ø 3.0 MM X 25 CM
38501-16	Ø 4.0 MM X 25 CM
38501-17	Ø 4.0 MM X 30 CM
38501-18	Ø 4.0 MM X 35 CM
38501-19	Ø 4.5 MM X 30 CM
38501-20	Ø 4.5 MM X 35 CM

THREADED HUB

38502-11	Ø 3.0 MM X 15 CM
38502-12	Ø 3.0 MM X 25 CM
38502-13	Ø 4.0 MM X 15 CM
38502-14	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38503-11	Ø 2.0 MM X 15 CM
38503-12	Ø 2.5 MM X 15 CM
38503-13	Ø 2.5 MM X 20 CM
38503-14	Ø 3.0 MM X 15 CM
38503-15	Ø 3.0 MM X 20 CM
38503-16	Ø 3.0 MM X 25 CM

LUER LOCK

38504-11	Ø 2.0 MM X 15 CM
38504-12	Ø 2.5 MM X 15 CM
38504-13	Ø 2.5 MM X 20 CM
38504-14	Ø 3.0 MM X 15 CM
38504-15	Ø 3.0 MM X 20 CM
38504-16	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38505-11	Ø 2.0 MM X 15 CM
38505-12	Ø 2.5 MM X 15 CM
38505-13	Ø 2.5 MM X 20 CM
38505-14	Ø 3.0 MM X 15 CM
38505-15	Ø 3.0 MM X 20 CM
38505-16	Ø 3.0 MM X 25 CM

POWER CANNULA


38506-51	Ø 2 MM X 15 CM
38506-52	Ø 2 MM X 20 CM
38506-53	Ø 2 MM X 25 CM
38506-54	Ø 2 MM X 30 CM
38506-55	Ø 2 MM X 35 CM
38506-56	Ø 3 MM X 15 CM
38506-57	Ø 3 MM X 20 CM
38506-58	Ø 3 MM X 25 CM
38506-59	Ø 3 MM X 30 CM
38506-60	Ø 3 MM X 35 CM
38506-61	Ø 4 MM X 15 CM
38506-62	Ø 4 MM X 20 CM
38506-63	Ø 4 MM X 25 CM
38506-64	Ø 4 MM X 30 CM
38506-65	Ø 4 MM X 35 CM
38506-66	Ø 5 MM X 15 CM
38506-67	Ø 5 MM X 20 CM
38506-68	Ø 5 MM X 25 CM
38506-69	Ø 5 MM X 30 CM
38506-70	Ø 5 MM X 35 CM
38506-71	Ø 6 MM X 15 CM
38506-72	Ø 6 MM X 20 CM
38506-73	Ø 6 MM X 25 CM
38506-74	Ø 6 MM X 30 CM
38506-75	Ø 6 MM X 35 CM



SPECIALTY CANNULAS



MADE IN THE USA



FIELD FLAP

WITH ITS ELEVATED LEADING EDGE, THIS CANNULA IS DESIGNED TO OPEN AND EXPOSE SOFT TISSUE FOR MORE PRECISE SUCTIONING ALONG THE MANDIBLE AND SUB-MENTAL AREAS.

60CC TOMMEY HUB

38501-21	Ø 2.5 MM X 15 CM
38501-22	Ø 2.5 MM X 20 CM
38501-23	Ø 3.0 MM X 15 CM
38501-24	Ø 3.0 MM X 20 CM
38501-25	Ø 3.0 MM X 25 CM
38501-26	Ø 4.0 MM X 25 CM
38501-27	Ø 4.0 MM X 30 CM
38501-28	Ø 4.0 MM X 35 CM
38501-29	Ø 4.5 MM X 30 CM
38501-30	Ø 4.5 MM X 35 CM

THREADED HUB

38502-21	Ø 3.0 MM X 15 CM
38502-22	Ø 3.0 MM X 25 CM
38502-23	Ø 4.0 MM X 15 CM
38502-24	Ø 4.0 MM X 25 CM

ULTRA LUER LOCK

38503-21	Ø 2.0 MM X 15 CM
38503-22	Ø 2.5 MM X 15 CM
38503-23	Ø 2.5 MM X 20 CM
38503-24	Ø 3.0 MM X 15 CM
38503-25	Ø 3.0 MM X 20 CM
38503-26	Ø 3.0 MM X 25 CM

LUER LOCK

38504-21	Ø 2.0 MM X 15 CM
38504-22	Ø 2.5 MM X 15 CM
38504-23	Ø 2.5 MM X 20 CM
38504-24	Ø 3.0 MM X 15 CM
38504-25	Ø 3.0 MM X 20 CM
38504-26	Ø 3.0 MM X 25 CM

LUER LOCK (STAINLESS STEEL)

38505-21	Ø 2.0 MM X 15 CM
38505-22	Ø 2.5 MM X 15 CM
38505-23	Ø 2.5 MM X 20 CM
38505-24	Ø 3.0 MM X 15 CM
38505-25	Ø 3.0 MM X 20 CM
38505-26	Ø 3.0 MM X 25 CM

POWER CANNULA

38506-101	Ø 2 MM X 15 CM
38506-102	Ø 2 MM X 20 CM
38506-103	Ø 2 MM X 25 CM
38506-104	Ø 2 MM X 30 CM
38506-105	Ø 2 MM X 35 CM
38506-106	Ø 3 MM X 15 CM
38506-107	Ø 3 MM X 20 CM
38506-108	Ø 3 MM X 25 CM
38506-109	Ø 3 MM X 30 CM
38506-110	Ø 3 MM X 35 CM
38506-111	Ø 4 MM X 15 CM
38506-112	Ø 4 MM X 20 CM
38506-113	Ø 4 MM X 25 CM
38506-114	Ø 4 MM X 30 CM
38506-115	Ø 4 MM X 35 CM
38506-116	Ø 5 MM X 15 CM
38506-117	Ø 5 MM X 20 CM
38506-118	Ø 5 MM X 25 CM
38506-119	Ø 5 MM X 30 CM
38506-120	Ø 5 MM X 35 CM
38506-121	Ø 6 MM X 15 CM
38506-122	Ø 6 MM X 20 CM
38506-123	Ø 6 MM X 25 CM
38506-124	Ø 6 MM X 30 CM
38506-125	Ø 6 MM X 35 CM



SINGLE USE INSTRUMENTS



MADE IN THE USA



ULTRA LUER LOCK



TUMESCENT INFILTRATOR

THIS CANNULA IS USED SPECIFICALLY FOR THE INFUSION OF TUMESCENT SOLUTION AND FOR HYDRAULIC DISSECTION. THE CONFIGURATION OF THE PORTS ALLOWS EVEN DISTRIBUTION OF THE FLUID.

38601-01	Ø 2.0 MM X 12CM
38601-02	Ø 2.0 MM X 20CM



TONNARD HARVESTER

A MULTIPORT HARVESTER WITH 1MM ELEVATED SHARP MICROPORTS. THE 20 PORT CONFIGURATION ALLOWS FOR HARVESTING BOTH ON THE FORWARD AND BACKWARD STROKES.

38602-01	Ø 3.0 MM X 20CM
38602-02	Ø 2.4 MM X 20CM



TONNARD HARVESTER

EXCELLENT FOR HARVESTING SMALL ALIQUOTS OF TISSUES.

38603-01	Ø 2.5 MM X 15CM
38603-02	Ø 2.0 MM X 15CM



CARRAWAY HARVESTER

SPIRAL THREE PORT DESIGN EFFECTIVE FOR MULTI-PLANE TISSUE HARVESTING.

38604-01	Ø 2.5 MM X 15CM
38604-02	Ø 2.0 MM X 10CM
38604-03	2.0 MM X 15CM
38604-04	Ø 1.5 MM X 8CM
38604-05	Ø 1.5 MM X 10CM



INJECTOR

ROUND BLUNT TIP DESIGN PASSES THROUGH SOFT TISSUE EASILY WITHOUT THE TRAUMA OR POTENTIAL PROBLEMS OF SHARP INJECTORS.

38605-01	Ø 1.5 MM X 9CM
38605-02	Ø 1.5 MM X 7CM
38605-03	Ø 1.0 MM X 5CM
38605-04	Ø 1.0 MM X 4CM



SPOON TIP INJECTOR

BLUNT INJECTOR DESIGNED WITH A MORE OPEN PORT FOR RAPID INJECTION OF LARGER VOLUMES - A TORSO AND BUTTOCK WORKHOUSE.

38606-01	Ø 3.0 MM X 20CM
38606-02	Ø 2.40MM X 15CM



MILLER SPEED HARVESTER

38607-01	Ø 2.0 MM X 15CM
-----------------	-----------------



LIPOSUCTION ACCESSORIES




MADE IN THE USA

HANDLES FOR CANNULAS

ANAEROBIC TRANSFER



LUER LOCK HANDLE | **38701-01**



ADAPTOR HANDLE
STYLE I

38702-01

20CC



ADAPTOR HANDLE
STYLE II

38703-01

20CC



LIPOSUCTION ACCESSORIES



MADE IN THE USA

HANDLES FOR CANNULAS



38704-01

WITH CONNECTION THREAD FOR
CANNULAS HOSE CONNECTION



38705-01

WITH CONNECTION THREAD FOR
CANNULAS HOSE CONNECTION



38706-01

WITH CONNECTION THREAD FOR
CANNULAS HOSE CONNECTION



38707-01

WITH CONNECTION THREAD FOR
CANNULAS HOSE CONNECTION



38750-01

THE POCAR IS A PERCUTANEOUS PUNCH THAT IS USED TO MAKE A CIRCUMFERENTIAL INCISION PRIOR TO INFILTRATION AND TISSUE EXTRACTION. IT IS DESIGNED TO MAINTAIN A TIGHT SEAL WHILE MINIMIZING SCARRING.



LIPOSUCTION ACCESSORIES



MADE IN THE USA

NANO TRANSFER

ANAEROBIC TRANSFER



38751-01 | NANO TRANSFER



38752-01 | MESH FOR NANO TRANSFER

ANAEROBICS TRANSFERS

38760-01	LUER 1.2MM
38760-02	LUER 1.4MM
38760-03	LUER 2.4MM





LIPOSUCTION ACCESSORIES



MADE IN THE USA

38771-01 | LUER TO LUER



DECANT STANDS

38772-02	60CC TO LUER
38773-03	LUER TO LUER



SNAPLOKS

THESE SYRINGE SNAP LOCKS FIT INTO THE LUNGER FLANGES OF THE SYRINGE. WHEN THE PLUNGER IS PULLED OUT THE LOCK SNAPS OPEN ON THE BARREL LIP THUS HOLDING THE SYRINGE UNDER VACUUM.

38901-01	60CC
38901-02	50CC
38901-03	30CC
38901-04	20CC
38901-05	10CC





LIPOSUCTION ACCESSORIES



MADE IN THE USA

ANAEROBIC TRANSFER



UNIVERSAL POWER INJECTOR

38801-01

THE POWER INJECTOR ALLOWS THE SURGEON TO RE-INJECT SOFT TISSUE THROUGH LARGE OR SMALL DIAMETER INJECTORS WITH VERY LITTLE EFFORT AND MAXIMUM CONTROL. THE UNIVERSAL POWER INJECTOR WILL HOLD A 10CC OR A 20CC BD SYRINGE OR A BD SYRINGE OR A 60CC MONOJECT SYRINGE.



UNIVERSAL POWER INJECTOR

38950-01


THE POWER INJECTOR ALLOWS THE SURGEON TO RE-INJECT SOFT TISSUE THROUGH LARGE OR SMALL DIAMETER INJECTORS WITH VERY LITTLE EFFORT AND MAXIMUM CONTROL. THE UNIVERSAL POWER INJECTOR WILL HOLD A 10CC OR A 20CC BD SYRINGE OR A BD SYRINGE OR A 60CC MONOJECT SYRINGE.



LIPOSUCTION ACCESSORIES



MADE IN THE USA



381000-01 | STERILIZATION BOX



SYRINGE DECANTING AND INJECTION

38707-01	60cc
38707-02	50cc
38707-03	30cc
38707-04	20cc
38707-05	10cc