

ELECTROSURGERY

ELMED REUSABLE MONOPOLAR SURGICAL FORCEPS

SCISSORS, HEMOSTATS, & SCALPEL HANDLES

Quality, value priced, reusable, American made surgical instruments by Elmed.

Clinical precision and efficacy in monopolar electrosurgery depend fundamentally on advanced instrument design and manufacturing standards. For nearly half a century, ELMED has led the industry in developing, designing, and manufacturing premium, reusable electrosurgical equipment and accessories. Handcrafted at our facility in Glendal Heights, Illinois, USA, our proprietary forceps are engineered from superior medical-grade stainless steel. To eliminate radiofrequency (RF) leakage, each instrument features a dual-layer, electrostatically applied Nylon/Kynar insulation coating. Furthermore, micro-polished tips minimize eschar accumulation, ensuring long-term reusability and seamless autoclave sterilization.

HANDSWITCHING

With attached cable and plug



FOOTSWITCHING

Multiple cable connector selection







ELMED HAND ACTIVATED MONOPOLAR FORCEPS



MADE IN THE USA

- ATTACHED 10-FOOT, HIGH GRADE SILICONE CABLE, MACHINED PLUG AND DOUBLE-INSULATED NYLON / KYNAR SHANKS
- THE INSTRUMENTS ARE REUSABLE AND FULLY AUTOCLAVABLE
- TIPS ARE MICRO-POLISHED FOR MINIMUM ESCHER BUILDUP

HAND SWITCHING WITH ATTACHED CABLE AND PLUG



TWO PRONG PLUG FOR VALLEY LAB, CONMED, BOVIE, ELMED ETC. GENERATOR




DUAL CONTACT PLUG FOR ELMED, WITH ROUND SOCKET

DESCRIPTION	ITEM NO.	TIP	ITEM NO.	TIP
BAYONET, SERRATED, STRAIGHT TIP, 5" (13 CM)	5501-05		5551-05	
BAYONET, SERRATED, CURVED TIP, 5" (13 CM)	5502-05		5552-05	
BAYONET, SERRATED, STRAIGHT TIP, 6" (15 CM)	5501-06		5551-06	
BAYONET, SERRATED, CURVED TIP, 6" (15 CM)	5502-06		5552-06	
BAYONET, SERRATED, STRAIGHT TIP, 8" (21 CM)	5501-08		5551-08	
BAYONET, SERRATED, CURVED TIP, 8" (21 CM)	5502-08		5552-08	
BAYONET, SERRATED, STRAIGHT TIP, 10" (25 CM)	5501-10		5551-10	
BAYONET, SERRATED, CURVED TIP, 10" (25 CM)	5502-10		5552-10	
POTTS-SMITH, 1X2 TEETH, STRAIGHT TIP, 7" (18 CM)	5510-07	N/A	5560-07	N/A
POTTS-SMITH, 1X2 TEETH, STRAIGHT TIP, 8" (21 CM)	5510-08	N/A	5560-08	N/A
POTTS-SMITH, 1X2 TEETH, STRAIGHT TIP, 10" (25 CM)	5510-10	N/A	5560-10	N/A
POTTS-SMITH, SERRATED, STRAIGHT TIP, 7" (18 CM)	5511-07		5561-07	
POTTS-SMITH, SERRATED, CURVED TIP, 7" (18 CM)	5512-07		5562-07	
POTTS-SMITH, SERRATED, STRAIGHT TIP, 8" (21 CM)	5511-08		5561-08	
POTTS-SMITH, SERRATED, CURVED TIP, 8" (21 CM)	5512-08		5562-08	
POTTS-SMITH, SERRATED, STRAIGHT TIP, 10" (25 CM)	5511-10		5561-10	
POTTS-SMITH, SERRATED, CURVED TIP, 10" (25 CM)	5512-10		5562-10	
CUSHING, 1X2 TEETH, STRAIGHT TIP, 7" (18 CM)	5513-07	N/A	5563-07	N/A
CUSHING, 1X2 TEETH, STRAIGHT TIP, 8" (21 CM)	5513-08	N/A	5563-08	N/A
CUSHING, 1X2 TEETH, STRAIGHT TIP, 10" (25 CM)	5513-10	N/A	5563-10	N/A
CUSHING, SERRATED, STRAIGHT TIP, 7" (18 CM)	5514-07		5564-07	
CUSHING, SERRATED, CURVED TIP, 7" (18 CM)	5515-07		5565-07	
CUSHING, SERRATED, STRAIGHT TIP, 8" (21 CM)	5514-08		5564-08	
CUSHING, SERRATED, CURVED TIP, 8" (21 CM)	5515-08		5565-08	
CUSHING, SERRATED, STRAIGHT TIP, 10" (25 CM)	5514-10		5564-10	
CUSHING, SERRATED, CURVED TIP, 10" (25 CM)	5515-10		5565-10	
GERALD, SERRATED, STRAIGHT TIP, 7" (18 CM)	5521-07		5571-07	
GERALD, 1X2 TEETH, STRAIGHT TIP, 7" (25 CM)	5522-07	N/A	5572-07	N/A
ADSON, SERRATED, STRAIGHT TIP, 4 1/2" (18 CM)	5523-04		5583-04	
ADSON, 1X2 TEETH, STRAIGHT TIP, 4 1/2" (25 CM)	5524-04	N/A	5584-04	N/A
JEWELERS, SMOOTH – STRAIGHT TIP, 3 1/2" (9 CM)	5225-03		5295-03	
JEWELERS, SMOOTH – CURVED TIP, 3 1/2" (9 CM)	5226-03		5296-03	

SAMPLE: LIST NUMBER ORDERING INSTRUCTIONS

TIP SIZE SPECIFICATIONS

INDICATE:	ITEM NO.	TIP SIZE LETTER
SAMPLE:	5000-00	+ A, B, C, D, E



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE




ELMED FOOT-SWITCHING MONOPOLAR FORCEPS



MADE IN THE USA

- DOUBLE-INSULATED NYLON/KYNAR SHANKS AVAILABLE w/ 3 DIFFERENT CABLE CONNECTOR TYPES.
 - THE INSTRUMENTS ARE REUSABLE AND FULLY AUTOCLAVABLE.
 - TIPS ARE MICRO-POLISHED FOR MINIMUM ESCHER BUILDUP.



	ITEM NO.	TIP	ITEM NO.	TIP	ITEM NO.	TIP
BAYONET, SERRATED, STRAIGHT TIP , 5" (13 CM)	5201-05		5301-05		5401-05	
BAYONET, SERRATED, CURVED TIP , 5" (13 CM)	5202-05		5302-05		5402-05	
BAYONET, SERRATED, STRAIGHT TIP , 6" (15 CM)	5201-06		5301-06		5401-06	
BAYONET, SERRATED, CURVED TIP , 6" (15 CM)	5202-06		5302-06		5402-06	
BAYONET, SERRATED, STRAIGHT TIP , 8" (21 CM)	5201-08		5301-08		5401-08	
BAYONET, SERRATED, CURVED TIP , 8" (21 CM)	5202-08		5302-08		5402-08	
BAYONET, SERRATED, STRAIGHT TIP , 10" (25 CM)	5201-10		5301-10		5401-10	
BAYONET, SERRATED, CURVED TIP , 10" (25 CM)	5202-10		5302-10		5402-10	
POTTS-SMITH, 1x2 TEETH, STRAIGHT TIP, 7" (18 CM)	5210-07	N/A	5310-07	N/A	5410-07	N/A
POTTS-SMITH, 1x2 TEETH, STRAIGHT TIP, 8" (21 CM)	5210-08	N/A	5310-08	N/A	5410-08	N/A
POTTS-SMITH, 1x2 TEETH, STRAIGHT TIP, 10" (25 CM)	5210-10	N/A	5310-10	N/A	5410-10	N/A
POTTS-SMITH, SERRATED, STRAIGHT TIP, 7" (18 CM)	5211-07		5311-07		5411-07	
POTTS-SMITH, SERRATED, CURVED TIP, 7" (18 CM)	5212-07		5312-07		5412-07	
POTTS-SMITH, SERRATED, STRAIGHT TIP, 8" (21 CM)	5211-08		5311-08		5411-08	
POTTS-SMITH, SERRATED, CURVED TIP, 8" (21 CM)	5212-08		5312-08		5412-08	
POTTS-SMITH, SERRATED, STRAIGHT TIP, 10" (25 CM)	5211-10		5311-10		5411-10	
POTTS-SMITH, SERRATED, CURVED TIP, 10" (25 CM)	5212-10		5312-10		5412-10	
CUSHING, 1x2 TEETH, STRAIGHT TIP, 7" (18 CM)	5213-07	N/A	5313-07	N/A	5413-07	N/A
CUSHING, 1x2 TEETH, STRAIGHT TIP, 8" (21 CM)	5213-08	N/A	5313-08	N/A	5413-08	N/A
CUSHING, 1x2 TEETH, STRAIGHT TIP, 10" (25 CM)	5213-10	N/A	5313-10	N/A	5413-10	N/A
CUSHING, SERRATED, STRAIGHT TIP, 7" (18 CM)	5214-07		5314-07		5414-07	
CUSHING, SERRATED, CURVED TIP, 7" (18 CM)	5215-07		5315-07		5415-07	
CUSHING, SERRATED, STRAIGHT TIP, 8" (21 CM)	5214-08		5314-08		5414-08	
CUSHING, SERRATED, CURVED TIP, 8" (21 CM)	5215-08		5315-08		5415-08	
CUSHING, SERRATED, STRAIGHT TIP, 10" (25 CM)	5214-10		5314-10		5414-10	
CUSHING, SERRATED, CURVED TIP, 10" (25 CM)	5215-10		5315-10		5415-10	
GERALD, SERRATED, STRAIGHT TIP, 7" (18 CM)	5221-07		5321-07		5421-07	
GERALD, 1x2 TEETH, STRAIGHT TIP, 7" (25 CM)	5222-07	N/A	5322-07	N/A	5422-07	N/A
ADSON, SERRATED, STRAIGHT TIP, 4 1/2" (18 CM)	5223-04		5323-04		5423-04	
ADSON, 1x2 TEETH, STRAIGHT TIP, 4 1/2" (25 CM)	5224-04	N/A	5324-04	N/A	5424-04	N/A
JEWELERS, SMOOTH – STRAIGHT TIP, 3 1/2" (9 CM)	5225-03		5325-03		5425-03	
JEWELERS, SMOOTH – CURVED TIP, 3 1/2" (9 CM)	5226-03		5326-03		5426-03	

SAMPLE: LIST NUMBER ORDERING INSTRUCTIONS

TIP SIZE SPECIFICATIONS

INDICATE:	ITEM NO.	TIP SIZE LETTER	2.0 mm	1.5 mm	1.0 mm	0.5 mm	NEEDLE
SAMPLE:	5000-00	+ A, B, C, D, E	A	B	C	D	E



TO ORDER CALL: (224) 353-6446





MONOPOLAR HEMOSTATIC FORCEPS, SCALPEL HANDLES & SCISSORS



MADE IN THE USA

THE HEMOSTATS, SCALPEL HANDLES AND SCISSORS ARE FITTED WITH A BANANA RF-SOCKET FOR
FOOTSWITCH OPERATION.



ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
5574	HEMOSTATIC FORCEPS, STRAIGHT, INSULATED, 7" (18 CM)	5584	SCISSORS, CURVED, INSULATED HANDLE, 7" (18 CM)
5575	HEMOSTATIC FORCEPS, CURVED, INSULATED, 7" (18 CM)	5585	SCISSORS, STRAIGHT, INSULATED HANDLE, 8" (20 CM)
5576	HEMOSTATIC FORCEPS, STRAIGHT, INSULATED, 8" (20 CM)	5586	SCISSORS, CURVED, INSULATED HANDLE, 8" (20 CM)
5577	HEMOSTATIC FORCEPS, CURVED, INSULATED, 8" (20 CM)	5587	SCISSORS, STRAIGHT, INSULATED HANDLE & OUTSIDE BLADES, 7" (18 CM)
5581	SCALPEL HANDLE NO. 3, INSULATED	5588	SCISSORS, CURVED, INSULATED HANDLE & OUTSIDE BLADES, 7" (18 CM)
5582	SCALPEL HANDLE NO. 4, INSULATED	5589	SCISSORS, STRAIGHT, INSULATED HANDLE & OUTSIDE BLADES, 8" (20 CM)
5583	SCISSORS, STRAIGHT, INSULATED HANDLE, 7" (18 CM)	5590	SCISSORS, CURVED, INSULATED HANDLE & OUTSIDE BLADES, 8" (20 CM)

MONOPOLAR RF-CABLES FOR USE W/ FOOT SWITCHING MONOPOLAR INSTRUMENTS:

ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
5021	MONOPOLAR RF-CABLE FOR INSTRUMENTS WITH BANANA PLUG RF-SOCKET AND UNITS WITH BANANA PLUG RF-SOCKET	5022	MONOPOLAR RF-CABLE FOR INSTRUMENTS WITH OLSEN TYPE SCREW-ON RF-CONNECTOR AND UNITS WITH BANANA PLUG RF-SOCKET
5009	MONOPOLAR RF-CABLE FOR INSTRUMENTS WITH BANANA PLUG RF-SOCKET AND UNITS WITH BOVIE TYPE RF-SOCKET	5023	MONOPOLAR RF-CABLE FOR INSTRUMENTS WITH OLSEN TYPE SCREW-ON RF-CONNECTOR AND UNITS WITH BOVIE TYPE RF-SOCKET



TO ORDER CALL: (224) 353-6446