

SELF RETAINING RETRACTOR SYSTEMS



WE DESIGN, MANUFACTURE, & SELL THE TOOLS THE SURGEONS USE

ELMED INC.
35 N. BRANDON DR.
GLENDALE HEIGHTS,
IL 60139
USA

WE SUBSCRIBE TO
COST
CONTAINMENT AND
PROTECTION OF
THE ENVIRONMENT



THEREFORE, WE
MANUFACTURE
REUSABLE
PRODUCTS FOR A
CLEANER WORLD

PH (224) 353-6446
FAX (224) 653-8178
EMAIL
MEDICAL@ELMED.COM
WWW.ELMED.COM

TABLE OF CONTENTS


	PAGE
RETRACT ROBOT	1-2
NATHANSON LIVER RETRACTORS	3-5
GOOSE NECK RETRACTOR	6
FLEX FRAME RETRACTOR SET	7
MULTI-RING TABLE MOUNT RETRACTOR	8
MULTI-FLEX TABLE MOUNT RETRACTOR	9
FRAME RETRACTOR SYSTEM	10



ELMED RETRACT ROBOT IN LAPAROSCOPIC SURGERY



MADE IN THE USA



VARIOUS INSTRUMENTS SHOWN WITH THE SINGLE-ARM RETRACT- ROBOT:

**1) WITH THE NATHANSON
LIVER RETRACTOR**

**2) WITH THE ARTICULATING
SCOPE HOLDER: QUICK,
STABLE & PRECISE
ADJUSTMENT CAN BE MADE
WITH THE TURN OF A KNOB.**

**3) WITH THE FAN
RETRACTOR; ANY TYPE OF**

THE ELMED STERIL-CLAMP:

PROVIDE SAFE AND STABLE
MOUNTING IN STERILE CONDITION
ON TOP OF THE DRAPE ON ANY
OPERATING TABLE RAIL

ELMED McLUCAS SCOPE/
INSTRUMENT HOLDER W/
ARTICULATION JOINT & QUICK-
RELEASE SOFT-ALLOY SELF-
RETAINING CLAMP

ELMED EXTENSION CLAMP FOR
BARIATRIC SURGERY (FITS INTO
STERILE-CLAMP)

#15090-57

ITEM No.	DESCRIPTION
15088-00	ELMED RETRACT-ROBOT SINGLE ARM INSTRUMENT
15088-01	ELMED RETRACT-ROBOT SINGLE ARM INSTRUMENT W/ EXTRA-LONG MOUNTING ROD (26") (BARIATRIC)
15090-53	ELMED STERIL-CLAMP HEAVY DUTY, 5/8" ROUND HOLE OR 1/2" X 1" SLOT (FOR EXTRA LONG ROBOT ARM)
15090-57	ELMED EXTENSION CLAMP FOR BARIATRIC SURGERY: AUTOCLAVABLE, 5/8" ACCESSORY OPENING
52410-01	ELMED SCOPE/INSTRUMENT HOLDER W/ ARTICULATION JOINT SELF-RETAINING CLAMP (SLEEVE STYLE)
52400-00	ELMED ENDOSCOPIC FIXATION DEVICE CONSISTING OF: ITEM NO's: 15088-00 , 15090-53 , & 52410-01
52400-01	ELMED BARIATRIC ENDOSCOPIC FIXATION DEVICE CONSISTING OF ITEM NO's: 15088-01 , 15090-53 , & 52410-01



WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



ELMED RETRACT ROBOT, IRON INTERN™ STYLE, IN LAPAROSCOPIC SURGERY



MADE IN THE USA



ITEM NO.	DESCRIPTION
15088-08	ELMED RETRACT-ROBOT TABLE RETRACTOR ARM, 22" LONG WITH "SLIP-IN" INSTRUMENT CLAMP, ROTATING LOCK AND SLIP-IN INTERCHANGEABLE UNIVERSAL INSTRUMENT CLAMP. AUTOCLAVABLE
15090-53	ELMED STERIL-CLAMP HEAVY DUTY, 5/8" ROUND HOLE OR 1/2" X 1" SLOT (FOR EXTRA LONG ROBOT ARM)
15090-57	ELMED EXTENSION CLAMP FOR BARIATRIC SURGERY: AUTOCLAVABLE, 5/8" ACCESSORY OPENING
15088-08-01	SCOPE HOLDER CLAMP WITH 8MM POST
15088-08-02	INSTRUMENT HOLDER CLAMP WITH 8MM POST



TO ORDER CALL: (224) 353-6446



NATHANSON LIVER RETRACTORS



MADE IN THE USA



ITEM NO.	DESCRIPTION
52420-00-H	ELMED NATHANSON STYLE LIVER RETRACTORS 6.5MM DIA. SET OF 3 - SIZES: SMALL MEDIUM LARGE (8MM)
52420-01-H	ELMED NATHANSON STYLE LIVER RETRACTORS 6.5MM DIA.- SMALL (8MM)
52420-02-H	ELMED NATHANSON STYLE LIVER RETRACTORS 6.5MM DIA.- MEDIUM (8MM)
52420-03-H	ELMED NATHANSON STYLE LIVER RETRACTORS 6.5MM DIA.- LARGE (8MM)
52430-00-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA. SET OF 3 - SIZES: SMALL MEDIUM LARGE (8MM)
52430-01-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA.- SMALL (8MM)
52430-02-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA.- MEDIUM (8MM)
52430-03-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA.- LARGE (8MM)
52430-04-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA.- X-SMALL (8MM)
52430-05-H	ELMED NATHANSON STYLE LIVER RETRACTORS 5MM DIA.- X-LARGE (8MM)
52460-00-H	ELMED NATHANSON STYLE LIVER RETRACTORS 3.0MM DIA. SET OF 3 - SIZES: SMALL MEDIUM LARGE (8MM)
52460-03-H	ELMED NATHANSON STYLE LIVER RETRACTORS 3.0MM DIA.- SMALL (8MM)
52460-04-H	ELMED NATHANSON STYLE LIVER RETRACTORS 3.0MM DIA.- MEDIUM (8MM)
52460-05-H	ELMED NATHANSON STYLE LIVER RETRACTORS 3.0MM DIA.- LARGE (8MM)
52460-06-H	ELMED NATHANSON STYLE LIVER RETRACTORS 3.0MM DIA.- X-SMALL (8MM)





NATHANSON LIVER RETRACTORS FOR ROBOTIC USE



MADE IN THE USA



ITEM NO.	DESCRIPTION
52409-00	SET OF 6.5MM NATHANSON RETRACTORS, ROBOTIC SHAFT 9CM: SMALL, MEDIUM, LARGE
52409-01	SMALL 6.5MM NATHANSON RETRACTOR, ROBOTIC SHAFT 9CM
52409-02	MEDIUM 6.5MM NATHANSON RETRACTOR, ROBOTIC SHAFT 9CM
52409-03	LARGE, 6.5MM, ROBOTIC SHAFT 9CM
52409-04	SET OF 5.0MM NATHANSON RETRACTORS, ROBOTIC SHAFT 9CM: SMALL, MEDIUM, LARGE
52409-05	SMALL 5.0MM NATHANSON RETRACTOR, ROBOTIC SHAFT 9CM
52409-06	MEDIUM 5.0MM NATHANSON RETRACTOR, ROBOTIC SHAFT 9CM
52409-07	LARGE 5.0MM NATHANSON RETRACTOR, ROBOTIC SHAFT 9CM
52412-00	SET OF 6.5MM NATHANSON RETRACTORS, ROBOTIC SHAFT 12CM: SMALL, MEDIUM, LARGE
52412-01	SMALL 6.5MM NATHANSON RETRACTOR, ROBOTIC SHAFT 12CM
52412-02	MEDIUM 6.5MM NATHANSON RETRACTOR, ROBOTIC SHAFT 12CM
52412-03	LARGE 6.5MM NATHANSON RETRACTOR, ROBOTIC SHAFT 12CM
52450-00	SET OF 5.0MM NATHANSON HOOKS: SMALL, MEDIUM, LARGE
52450-01	SMALL 5.0MM NATHANSON RETRACTOR
52450-02	MEDIUM 5.0MM NATHANSON RETRACTOR
52450-03	LARGE 5.0MM NATHANSON RETRACTOR
52450-04	EXTRA SMALL 5.0MM NATHANSON RETRACTOR



TO ORDER CALL: (224) 353-6446



NATHANSON LIVER RETRACTORS FOR ROBOTIC USE



MADE IN THE USA



52440-00	SET OF 6.5MM NATHANSON HOOKS: SMALL, MEDIUM, LARGE
52440-01	SMALL 6.5MM NATHANSON RETRACTOR
52440-02	MEDIUM 6.5MM NATHANSON RETRACTOR
52440-03	LARGE, 6.5MM NATHANSON RETRACTOR, LONG TIP, ROBOTIC SHAFT 9CM
15088-08-03	INSTRUMENT HOLDER CLAMP WITH HEX FITTING
15088-08-04	UNIVERSAL INSTRUMENT HOLDER CLAMP WITH HEX FITTING



ELMED GOOSE NECK RETRACTOR IN LAPAROSCOPIC SURGERY



MADE IN THE USA




FEATURES OF THE GOOSE-NECK SYSTEM:

- ALL COMPONENTS ARE MADE OF A HIGH GRADE STAINLESS STEEL.
- THE COMPLETE SYSTEM, INCLUDING THE STERIL-CLAMP™, MAY BE FLASH AUTOCLAVED
- THE SELF-RETAINING JOINTS ROTATE & TWIST IN EVERY DESIRABLE POSITION
- THE INSTRUMENT SHAFT FITS INTO TABLE ACCESSORY CLAMPS SUCH AS THOSE USED FOR STIRRUPS & LEG HOLDERS, AS WELL AS ALL OF THE OTHER OVER-THE-DRAPE CLAMPS WITH AN OPENING OF 5/8"
- THE VERTICAL SHAFT (5/8") MAY BE LOWERED OF THE HEIGHTENED WITHIN THE CLAMP FOR FLEXIBLE ADJUSTMENT
- TWO DIFFERENT MODELS WITH DIFFERENT INSTRUMENT HOLDERS AVAILABLE

ITEM No.	DESCRIPTION
15090-60	ELMED GOOSE-NECK RETRACTOR SYSTEM WITH SELF-RETAINING SCOPE CLAMP WITH ARTICULATING JOINT FOR HOLDING TELESCOPE INCLUDES 15090-53: STERIL-CLAMP™ FOR ATTACHING TO TABLE
15090-61	ELMED GOOSE-NECK RETRACTOR SYSTEM WITH UNIVERSAL MOUNTING CLAMP FOR HOLDING ALL TYPES OF INSTRUMENTS INCLUDES 15090-53: STERIL-CLAMP™ FOR ATTACHING TO TABLE



RETRACT ROBOT & GOOSE NECK RETRACTOR STERILIZATION CONTAINERS



ITEM No.	DESCRIPTION
55040-03	ELMED ALUMINUM STORAGE/STERILIZATION CONTAINER W/ PERFORATED LID & BOTTOM, W/ BRACKETS FOR MOUNTING RETRACT-ROBOT #15088-00 & #15090-53. (FOR WRAPPING) DIMENSIONS: 26"L x 10"W x 7"H (63.5CM x 25CM x 13CM).
55040-02	ELMED ALUMINUM STORAGE/STERILIZATION CONTAINER W/ PERFORATED LID & BOTTOM, W/ BRACKETS FOR MOUNTING RETRACT-ROBOT #15088-01 OR 15088-03 & STERIL-CLAMP #15090-53. (FOR WRAPPING.) DIMENSIONS: 30"L x 10"W x 6"H (76CM x 25CM x 15CM).






FLEX FRAME RETRACTOR SET



MADE IN THE USA



18-88502 | COMPLETE “OMNI-STYLE” STANDARD RETRACTOR SET

CONSISTING OF:

ITEM NO.	DESCRIPTION	DIMENSIONS		QTY
18-10244	TABLE CLAMP W/VERTICAL POST			1
18-4040	RETRACTOR CLAMP			6
18-4020	VISION FRAME			1
18-3335	MAYO SWIVEL	2 ¾ X 2 IN	70 X 51 MM	2
18-3784	DEAVER SWIVEL	2 ¾ X 4 ½ IN	70 X 114 MM	1
18-3796	HARRINGTON “SWEETHEART”	2 X 5 IN	51 X 127 MM	1
18-3862	KELLY SWIVEL “MODIFIED”	3 X 2 ½ IN	76 X 64 MM	1
18-3864	KELLY SWIVEL “MODIFIED”	3 X 4 IN	76 X 102 MM	1
18-3788	MALLEABLE SWIVEL, PERFORATED	2 X 6 IN	51 X 152 MM	1
18-3790	MALLEABLE SWIVEL, PERFORATED	3 X 6 ½ IN	76 X 165 MM	1
18-3772	FENCE SWIVEL, PERFORATED	4 X 6 IN	102 X 152 MM	1
18-3774	SPLANCHNIC SWIVEL	2 X 6 IN	51 X 152 MM	1
18-3424	STERILIZATION CASE			1



TO ORDER CALL: (224) 353-6446



MULTI-RING TABLE MOUNT RETRACTOR




MADE IN THE USA



15800-00 GENERAL SURGERY FRAME RETRACTOR SYSTEM

CONSISTING OF:



QTY	ITEM NO.	DESCRIPTION	
1	15800-10	TABLE POST 1"(2CM) x 19"(49CM)	A
1	15800-12	HORIZONTAL FLEX BAR 3/4"(2CM) DIAMETER	B
1	15800-15	S.S. POST COUPLING 1" x 1"(2.5CM x 2.5CM)	C
1	15800-20	OVAL MEDIUM RING 9"x9"x19.5 I.D.	D
1	15800-30	SEGMENTED MEDIUM RING 9.5"x19.5" I.D.	E
2	15800-40	RATCHET MECHANISM	F
6	15800-41	TILT RATCHET MECHANISM	G
2	15800-50	KELLY RETRACTOR 2"x3"	H
1	15800-51	KELLY RETRACTOR 2"x4"	
1	15800-52	KELLY RETRACTOR 2"x5"	
1	15800-53	KELLY RETRACTOR 2"x6"	I
1	15800-60	BALFOUR RETRACTOR 3 3/4" x 2 3/4"	
1	15800-65	HARRINGTON RETRACTOR 2.5"x1"x6"	J
1	15800-70	MALLEABLE RETRACTOR 1.5"x6" (38MM x 152MM)	K
2	15800-71	MALLEABLE RETRACTOR 2"x6" (51MM x 152MM)	




WE DESIGN, MANUFACTURE & SELL THE TOOLS THE SURGEONS USE



MULTI-FLEX TABLE MOUNT RETRACTOR



MADE IN THE USA



ITEM NO.	DESCRIPTION
15900-01	MULTI-FRAME TABLE MOUNT RETRACTOR SYSTEM
CONSISTING OF:	
15099-91	MULTI-FLEX ARM WITH QUICK RELEASE, ROUND ENDO SCOPE AND INSTRUMENT CLAMP, 1 x 6' (2.5CM x 15CM) MOUNTING BAR
1800-11	TABLE POST, 1" (2.5CM) x 19"(49CM)
15800-15	S.S. POST COUPLING CLAMP, 1" x 1" (25.CM x 2.5CM)



TO ORDER CALL: (224) 353-6446



FRAME RETRACTOR SYSTEM



MADE IN THE USA



ITEM No.	DESCRIPTION
15090-100	COMPLETE DUAL FLEX-ARM RETRACTOR SYSTEM, INCLUDES: (1) 15090-98: 19" LONG REACH ARM W/ SWIVEL BALL ASSEMBLY (1) 15090-53: STERICLAMP RAIL CLAMP (2) 15090-63: FLEX ARM W/ FIXED INSTRUMENT (SCOPE) CONNECTOR
15090-101	COMPLETE DUAL FLEX ARM X - LONG RETRACTOR SYSTEM, INCLUDES: (1) 15090-99 X - LONG REACH MOUNTING POST / WITH SWIVEL BALL ASSEMBLY 22" LONG 5/8" DIA., VERTICAL BAR AND 18" LONG 1/2" DIA., HORIZONTAL BAR (1) 15090-53: STERICLAMP RAIL CLAMP (2) 15090-63: FLEX ARM W/ FIXED INSTRUMENT (SCOPE) CONNECTOR





SERVICE AND REPAIR



MADE IN THE USA

IN NEED OF INSTRUMENT REPAIRS? SEND THEM TO THE *ER!*

ELMED REPAIR

60 YEARS OF EXPERIENCE AND KNOWLEDGE IN DESIGNING, ENGINEERING, MANUFACTURING AND REPAIRING FOR THE MEDICAL COMMUNITY HAS GIVEN ELMED, INC. ONE OF THE BEST REPUTATIONS IN THE INDUSTRY. WE REPAIR WHAT OTHERS ARE UNABLE OR UNWILLING TO DO AND WE REPAIR ALL COMPETITIVE BRANDS. OUR WORK IS PERFORMED BY SKILLED TECHNICIANS, ALL OF WHICH HAVE BEEN WITH ELMED FOR THE PAST 10-30 YEARS.

COMPLETE IN-HOUSE MACHINE SHOP, POWDER COATING, LASER WELDING, & MARKING FINISHING, ELECTRONICS AND OPTICAL DEPARTMENTS

AREAS OF SPECIALTY

- LAPAROSCOPIC AND ENDOSCOPIC INSTRUMENTS
- ELECTROSURGICAL UNITS AND ACCESSORIES
- FIBER OPTIC LIGHT CABLES AND RETRACTORS
- RIGID ENDOSCOPES & FLEXIBLE SINUSCOPES
- STAINLESS STEEL HAND INSTRUMENTS
- LIGHT SOURCES AND HEADLIGHTS
- DEBRIDER UNITS, HAND PIECES AND BLADES
- SELF RETAINING RETRACTOR SYSTEMS

HOW TO PACKAGE AND SHIP INSTRUMENTS

- IMPORTANT:** ELMED CANNOT ACCEPT INSTRUMENTS THAT HAVE NOT BEEN PROPERLY CLEANED, DISINFECTED OR STERILIZED.
- INCLUDE CUSTOMER NAME AND ALL CONTACT INFORMATION (NAME, ADDRESS, PHONE NO., FAX NO., EMAIL ADDRESS, ETC.)
- INCLUDE CUSTOMER P.O. (PURCHASE ORDER), IF USED
- INCLUDE REASON FOR RETURN - BRIEF DESCRIPTION OF REPAIR PROBLEM IF AVAILABLE
- WRAP INSTRUMENT FOR PROTECTION DURING TRANSPORT

PLEASE FILL OUT, DETACH, AND INCLUDE THE FORM BELOW WITH YOUR REPAIR:

FROM:	
FACILITY:	CONTACT PERSON:
ADDRESS:	
CITY, STATE, ZIP	
EMAIL:	PHONE NO.
PO. No.	QTY.
INSTRUMENT DESCRIPTION:	

(NOTE: PLEASE INCLUDE MODEL NO. AND SERIAL NO. IN INSTRUMENT DESCRIPTION FOR SCOPES)

FAST SERVICE: INCOMING REPAIRS COMING IN BY NEXT DAY AIR WILL BE COMPLETED WITHIN 24 HRS (WHEN PARTS ARE AVAILABLE)

SHIP TO: **ELMED** 35 N BRANDON DR., GLENDALE HEIGHTS, IL 60139